B-Smart™

Instructions for Use

A disposable in-line manometer for monitoring injection pressure during administration of nerve blocks

Indications for Use: A disposable in-line manometer for monitoring injection pressure during administration of peripheral nerve blocks.

Contraindications: There are no known contraindications.
The disposable manometer is sterile in undamaged package.
For single use only. **Do not reuse, clean, reprocess or re-sterilize.**
Store at Controlled Room Temperature.

Directions:

1. Use aseptic technique.
2. Open package and remove the B-Smart™ manometer.
3. Attach the female end of B-Smart™ manometer to the syringe and the luer lock to the tubing of the block needle.
4. Prime System (Syringe, B-Smart™ and Needle). The B-Smart™ piston must move upwards during priming in order for the device to be ready for use.
5. After a nerve is located using a typical technique, begin the injection after confirming that the needle is not placed into a blood vessel.
6. Throughout the injection, the piston indicator will continuously provide information on the pressure during injection.
7. The pressure ranges specified on the piston indicator are: <15 psi (white), 15-20 psi (yellow); >20 psi (orange).

Technical Characteristics: Priming volume 1ml max.

Caution: Opening injection pressure varies among different tissues, being the highest when the needle tip is lodged in low compliance tissues (e.g., roots of brachial plexus, tendon) and lowest when injected in the soft connective tissues (e.g., adipose tissues, perineural space). In animal and fresh human cadaver models, injection of local anesthetic through a nerve block needle (22 gauge - 25 gauge) resulting in opening injection pressures >20 psi when intraneural (intrafascicular) and <15 psi when perineural (extrafascicular). No data on opening injection pressure during intraneural (intrafascicular) injection is available in living humans, however needle-nerve contact is often characterized by opening pressure >15 psi before an injection can commence. Based on the available data, keeping the injection pressure during nerve block injection at <15 psi throughout the injection procedure is recommended.

Intraneural (intrafascicular) injection may increase the risk of nerve injury during administration of nerve blocks. The B-Smart™ Injection pressure Monitor will not reliably distinguish between harmful injection sites, effective injection sites, or sites where the injection will be ineffective. Once the opening pressure is reached and injectate begins to flow, the pressure information provided by B-Smart™ will depend on the needle length, needle internal diameter, viscosity of injectate, rate of injectate flow, and the specific characteristics of the tissues into which the injection is made.

Precautions:

1. The B-Smart™ has no user serviceable parts. Do not attempt to repair or alter the device.
2. The device is disposable and must be discarded after use in accordance with hospital, administrative and/or local government policy.
3. Do not use the device if the package is open or damaged. Do not use if sterile barrier is damaged.
4. Re-use of the device can increase a risk of infection and can compromise its functionality.

Warranty Statement

MACOSTA MEDITEA S.r.l. declines any responsibility for damages caused to the patient and/or the personnel responsible if the medical device is used in an incorrect or improper manner, or following incorrect handling and, above all, if used by unqualified personnel. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.

Caution: THE DEVICE MUST BE USED EXCLUSIVELY BY SPECIALIZED MEDICAL PERSONNEL. FEDERAL LAW (USA) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

To order contact Concert Medical, LLC at (781) 261-7400.

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Pour commander, contacter Concert Medical, LLC au (781) 261-7400.
MACOSTA MEDITEA Srl - 41037 MIRANDOLA (MO) - ITALIE