

APEX® Compounding System Non-Patch Related Security Updates

ISSUE: cifs-generic-0005: Weak LAN Manager hashing permitted (NTLMv1 allowed)

ACTION:

- Configure the policy value for Computer Configuration -> Windows Settings -> Security Settings -> Local Policies -> Security Options -> "Network security: LAN Manager authentication level" to "Send NTLMv2 response only. Refuse LM & NTLM".
- Configure the policy value for Computer Configuration >> Windows Settings >> Security Settings >> Local Policies >> Security Options >> "Network security: Minimum session security for NTLM SSP based (including secure RPC) clients" to "Require NTLMv2 session security" and "Require 128-bit encryption" (all options selected).

Adobe Reader at 20.012.20043

Adobe Flash removal

kb4577586 - needed to be pulled from the MS catalog for Windows 1607 x64

Last Updated: 3/26/21