

xevonta[®]

CHALLENGE THE THINKING

Choice of a suitable dialyzer


DV|VTeZ _ of a dialyzer is R_ important step e` UVTUZ_X the uremic toxin elimination strategy `WR patient. There are many different properties Rdd TReMU with a dialyzer design. xevonta Z developed by SR|R_TZ_X eYV membrane compositionç e` target an increased removal of uremic toxins + especially for the class of middle molecules. Other properties Z_T]f UV+

- Biocompatible polysulfone membrane
- High membrane permeability
- @ j XV_+ W gamma sterilization


What effect TR_ the dialyzer surface area have on clearance ?

Transport processes, either diffusive or convective, are the underlying mechanisms for elimination of substances during dialysis. : ToVRdv Z_ U|Rj kVc of dRTM dKv TR_ SV aR|R]M e` bfR_EZj dV^` gR] ` WTVceRZ_ ^` JVtf]VdZEYzd Zd UMaV_UV_e` _ ^` JVtf]V TYRdRTeVZdZdCgRZRSZdZ ` WdYV of SdR_TVd Z_ eYV U|WdV_e aRdV_e a`]d|R_U RTTVdd Z_ eYV S] ` UTZtf ZdA clinical study with twelve stable dialysis patients using twelve stable dialysis patients using different surface sizes of xevonta high-flux showed increases\UsPkt/V-values and >-reduction rates when larger U|Rj kVcsurface-sized were used " .

Influence ` Wdialyzer surface-size on spKt/V "




Influence ` Wdialyzer surface-size on B2M reduction rate "



What role does membrane permeability play in hemodialysis ?


During hemodialysis¹ the formation of a protein layer sometimes referred to as secondary membrane formation or membrane fouling² can effect the hydraulic permeability of a dialyzer.² Membrane permeability can therefore be a key determinator of performance.³

Membrane performance based on surface size and Kuf as an index of permeability of multiple membrane brands: Summary from different public sources



What is the in vitro performance of small molecules ?

Conditions: Blood flow = 400 ml/min; Dialysate flow = 500 ml/min; Ultrafiltration flow = 0 ml/min



¹ *Optiflux* + *Yeast extract* *Urea* *Creatinine* *Phosphate*
² *Optiflux* + *Yeast extract* *Urea* *Creatinine* *Phosphate*
³ *Optiflux* + *Yeast extract* *Urea* *Creatinine* *Phosphate*
⁴ *Optiflux* + *Yeast extract* *Urea* *Creatinine* *Phosphate*
⁵ *Optiflux* + *Yeast extract* *Urea* *Creatinine* *Phosphate*

REFERENCES:

¹ Yamamoto K, Hiwatari M, Kohori F, Sakai K, Fukuda M, Hiyoji T. Membrane fouling and dialysate flow pattern in an internal filtration enhancing dialyzer. J Artif Organs. 2005;8:198–205.

² Yamamoto K, Hiwatari M, Kohori F, Sakai K, Fukuda M, Hiyoji T. Membrane fouling and dialysate flow pattern in an internal filtration enhancing dialyzer. J Artif Organs. 2005;8:198–205.

³ Yamamoto K, Hiwatari M, Kohori F, Sakai K, Fukuda M, Hiyoji T. Membrane fouling and dialysate flow pattern in an internal filtration enhancing dialyzer. J Artif Organs. 2005;8:198–205.

xevonta[®]

TECHNICAL & IN-VITRO PERFORMANCE DATA

Measuring conditions and physical data according to ISO 8637

UF-coefficient: , 3° g/V 3° U Hct. 32%, total protein 60X, T = 37°C; Sieving coefficients: Q_B = 300 ml/min, Q_D = 500 ml/min

		Hi 18				Hi 20				Hi 23			
		500	500	500	800	500	500	500	800	500	500	500	800
		200	300	400	500	200	300	400	500	200	300	400	500
dialysate flow (Q _D) ml/min	Urea	198	281	341	414	199	287	349	427	199	290	354	439
blood flow (Q _B) ml/min	Creatinine	19%	263	304	372	196	271	316	390	197	276	324	403
Clearance: ultrafiltration flow (Q _F) = 0 ml/min	Phosphate	194	263	297		196	271	309		198	277	\$#!	
	Vitamin B ₁₂	155	184	210	239	161	195	220	259	166	204	227	272
S. C. (sieving coefficient)													
Q _B = 300 ml/min	β ₂ -Microglobulin												> 0.8
Q _F = 60 ml/min	Albumin												< 0.001
Ultrafiltration coefficient ml/h/mmHg	Q _B = 300 ml/min			99				111					124
KoA Urea (QB = 300 ml/min, B5 . & ! ^ Z^ Z_fi		" %&				" (" %				" *! !			
AZ^ Z_X G]f^ Vç3] ` UbUN (ml)		" "!				" #&				" %"			
Membrane material						amembris polysulfone							
Surface area (m ²)			1.8				2.0				2.3		
Sterilization						Oxygen-free gamma							
Wall thickness / inner diameter (µm)						35 /195							
Units per box						20							
Art. No.		7204403				7204404				7204405			

Data on file



7 c^` d\Z_Wc^ ReZ _cT _eRTe 3Z3dRf_ > VUZR] :_TZR
" !)!!!)%#! ' ' c` cV^ RZ fd Re ceUZ d1 SScRf_ZT ^

h h h Z33dRf_F D2Z ^ Z5Rj kVcd
Cl ^_j Z w#! '' 3Z3dRf_ > VUZR] :_TZR3VeY]VV^ A2c2]] dXVe dVdVgM