Infusomat® Space Pump IV Set

with Universal Spike, Backcheck Valve, CARESITE® Luer Access Device, and Anti Siphon Valve (ASV)

REF **480256**

LENGTH 122 in. (310 cm) (approx.)

PRIMING VOLUME 23 mL (approx.)

15
Drops/mL

CAUTIONS:

- CONTAINER HEIGHT SHOULD NOT EXCEED 40 in. ABOVE ASV.
- DO NOT CONNECT TO PATIENT WHILE PRIMING.
- SET CONTAINS REMOVABLE ASV.
- DO NOT TAMPER WITH LINE LOADING GUIDE ON PUMP SEGMENT.

Instructions for Use: Use Aseptic Technique.

- To access CARESITE Luer Access Device:** Swab top of luer access device vigorously for 15 seconds with 70% Isopropyl Alcohol and allow to air dry prior to use. Use a luer connector without a needle. After medication administration, flush per institutional protocol, attaching completely between flushes.
- Replace set per CDC guidelines.

- Press to prime set. Invert and tap components, while fluid is flowing. Ensure air is expelled.
- Monitor distal end of set. Press STOP to discontinue priming.
CAUTION: AFTER PRIMING IS COMPLETE, ENSURE THERE IS NO FLUID FLOW IN DRIP CHAMBER. IF DROPS ARE SEEN, CLOSE ROLLER CLAMP AND RELOAD OR REPLACE SET.
- Attach to desired connection. Proceed with pump operating instructions.

- CAUTION: VERIFY GREEN ANTI-FREE FLOW CLIP IS FULLY INSERTED AND FLASHING YELLOW TRIANGLE CEASES BEFORE CLOSING PUMP DOOR. IF RESISTANCE IS ENCOUNTERED WHEN CLOSING PUMP DOOR, REINSERT CLIP.**
- Push tubing at each end of pump into tubing channels. Run finger along IV set in pump to push tubing into air-in-line sensor. Close door. Press to proceed with line confirmation and pump operating instructions on pump display.
CAUTION: AFTER LOADING IS COMPLETE, ENSURE THERE IS NO FLUID FLOW IN DRIP CHAMBER. IF DROPS ARE SEEN, CLOSE ROLLER CLAMP AND RELOAD OR REPLACE SET. UNRESTRICTED FLUID FLOW MAY OCCUR IF SET IS NOT PROPERLY INSTALLED IN PUMP.

FOR PUMP ADMINISTRATION WITH ASV

PUMP PRIMING WITH PRIME FUNCTION ENABLED

- Close roller clamp. Fully spike container. Hang fluid container. Gently squeeze and release drip chamber until approximately 2/3 full.
NOTE: Vent cap should remain in the closed position unless infusing from a glass bottle.
- Press to turn on pump. Press to open door. Press when prompted to open door.
- NOTE:** Ensure yellow triangle in pump is flashing prior to set loading. If it is not, press down on opening green lever until a "click" is heard and it will then start to blink. Ensure Line Loading Guide is straight and above pump segment throughout loading procedure (stars on pump segment of tubing are in a straight line).
- Insert IV set from right to left with *white rectangular clip* closest to IV container inserted first. Then insert *white pronged clip*.
- Insert green anti-free flow clip in opened green slot with hook facing down and outward until opening green lever locks and tubing is secure.
CAUTION: VERIFY GREEN ANTI-FREE FLOW CLIP IS FULLY INSERTED AND FLASHING YELLOW TRIANGLE CEASES BEFORE CLOSING PUMP DOOR. IF RESISTANCE IS ENCOUNTERED WHEN CLOSING PUMP DOOR, REINSERT CLIP.
- Push tubing at each end of pump into tubing channels. Run finger along IV set in pump to push tubing into air-in-line sensor. Close door. Press to proceed with line confirmation.
CAUTION: UNRESTRICTED FLUID FLOW MAY OCCUR IF SET IS NOT PROPERLY INSTALLED IN PUMP.
- Tighten ASV and remove end cap.

FOR IV SET ADMINISTRATION WITH GRAVITY PRIMING WITH ASV

PUMP ADMINISTRATION with GRAVITY PRIMING and FLEXIBLE CONTAINER

- Close roller clamp. Fully spike container. Hang fluid container. Gently squeeze and release drip chamber until approximately 2/3 full.
NOTE: Vent cap should remain in the closed position unless infusing from a glass bottle.
- Remove ASV from set and discard.
- CAUTION: DO NOT CONNECT TO PATIENT WHILE PRIMING.**
Prime set by slowly opening roller clamp and allow fluid to fill tubing. Invert and tap components while fluid is flowing. Close roller clamp. Ensure air is expelled.
- Attach new sterile ASV to distal end of set and remove cap.
- Open roller clamp and apply gentle pressure to flexible container until drops form at end of ASV. Close roller clamp.
- Press to turn on pump. Press to open door. Press when prompted to open door.
- NOTE:** Ensure yellow triangle in pump is flashing prior to set loading. If it is not, press down on opening green lever until a "click" is heard and it will then start to blink. Ensure Line Loading Guide is straight and above pump segment throughout loading procedure (stars on pump segment of tubing are in a straight line).
- Insert IV set from right to left with *white rectangular clip* closest to IV container inserted first. Then insert *white pronged clip*.
- Insert green anti-free flow clip in opened green slot with hook facing down and outward until opening green lever locks and tubing is secure.

- Attach to desired connection. Proceed with pump operating instructions.

PUMP ADMINISTRATION with GRAVITY PRIMING and BOTTLE/SEMI-RIGID CONTAINER

- Perform steps 1 - 4 and 6 - 10 of Gravity Priming and Flexible Container section.
- Press for pump to prime ASV. Once ASV is fully primed, press STOP.
- Attach to desired connection. Proceed with pump operating instructions.

GRAVITY ADMINISTRATION AFTER SET REMOVAL FROM PUMP with ASV ATTACHED

- Close roller clamp. Disconnect set from patient. Remove set from pump. Remove ASV from set and discard. Connect set to patient.
- Open green anti-free flow clip and adjust roller clamp for desired flow rate.

WITHOUT ANTI SIPHON VALVE (ASV)

PUMP ADMINISTRATION with ASV REMOVED - Perform steps 1 - 3 and 6 - 11 of Gravity Priming and Flexible Container section.

GRAVITY ADMINISTRATION CAUTION: DO NOT USE ASV FOR GRAVITY INFUSIONS.

- Perform steps 1 - 3 of Gravity Priming and Flexible Container section. Attach to desired connection.
- Slowly open roller clamp and adjust for desired flow rate.

STERILE EO

Sterilized using ethylene oxide

Single sterile barrier system

Non-pyrogenic

Do not resterilize

Single use

Consult Electronic Instructions for Use
www.bbraunusa.com

Not made with natural rubber latex or DEHP.
Patent: www.bbraunusa.com/patents

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B. Braun Medical Inc.