



# EXCEL<sup>®</sup> Plus IV Container



## Helpful Hints (Page 1)

Refer to product labeling for full prescribing information.



### REMOVING THE OVERWRAP

Start a tear at any of the depressions on the notched edge. Notches can be found on both sides of the overwrap. Tear the overwrap by running your index finger through the seam of the overwrap.



### MEDICATION ADMINISTRATION

**STEP 1:** Prior to adding medications, swab the injection port. Allow to air dry (The EXCEL Plus medication port is not made with natural rubber latex).

**STEP 2:** Support the base of the medication port. Using a syringe and 18-22 gauge needle, inject medication through the center of injection site keeping fingers behind the flange. Mix thoroughly by inverting container several times.

### SPIKING THE CONTAINER

**STEP 3:** Remove the foil seal on the set port in two steps: lift the tab from side of port; peel off towards additive port.

**STEP 4:** Hold the flange of the set port between your thumb and forefinger. Push spike straight into port while twisting.

Alternatively, cradle EXCEL Plus in the forearm with the hanger at the elbow and hold the finger guard between the thumb and forefinger and insert the spike straight into the set port diaphragm. If positioning the container in your left arm, turn the container label down so the spike port is next to the left thumb and forefinger. If positioning the container in your right arm, turn the bag label up so the spike port is next to the right thumb and forefinger.



### SCANNING

**OPTION A:** Hang the container on an IV pole and scan at a 45°-90° degree angle, 2"-6" from the container.

**OPTION B:** Alternatively, hold the container by the hanger tab and scan at a 45°-90° degree angle, 2"-6" from the container.

*Do not squeeze the container in the middle while attempting to scan.*

# EXCEL<sup>®</sup> Plus IV Container

## Helpful Hints (Page 2)

Refer to product labeling for full prescribing information.



### SHIPPING AND STORAGE INSTRUCTIONS

- Do not stack individual cases more than 6 containers high
- Avoid freezing and excessive heat
- Maintain upright orientation of cases

### STACKING IN BINS

When stacking loose containers in a bin, do not stack them more than five (5) high.

### REMOVING AIR

- Insert an 18-22 gauge needle of appropriate length (at least 5/8 inch) through the medication port until the bevel of the needle just penetrates the medication port and inner diaphragm.
- Tilt container slightly to elevate fluid and remove as much air as possible by squeezing container continuously.
- Maintain pressure and remove needle immediately to prevent air backflow into the container.



Supporting you, your patients, and environmental safety, with the only full line of IV solutions not made with DEHP, PVC, or natural rubber latex.

[Learn more at bbrauncares.com](https://www.bbrauncares.com)



B. Braun Medical Inc. | Bethlehem, PA  
1-800-227-2862 | [bbraunusa.com](https://www.bbraunusa.com)

Rx Only ©2022 B. Braun Medical Inc. Bethlehem, PA. All rights reserved. 20-0408B\_02/22