

Improve Infusion Safety and Reduce Alarm Fatigue in the ICU

Cathy Sullivan, MSN, RN, FNP, CCRN, Director Patient Care Services

Mount Sinai Beth Israel

Presented at AACN NTI 2014

Objective

Our objective was to improve medication safety and reduce alarm fatigue through wireless integration of smart pumps, real-time monitoring and new data on dosing trends and practices.

Methodology

We implemented B.Braun Outlook® 400ES wireless smart pumps across a 4 hospital system including:

- A single, uniform drug library
- DoseTrac® real time monitoring & retrospective reporting
- DoseTrac® Analysis Service



A 2012 DoseTrac Analysis identified our top drugs associated with 51% of our override alerts at Beth Israel Hospital. We made modifications in dosing limits and practices, resulting in an overall 70% reduction in these alerts (Phase 1).

In 2013 a second DoseTrac Analysis was conducted comparing trends across the system, identifying our top drugs associated with 67% of system alerts. Additional improvements were made across the 4 hospital system (Phase 2).


Results Phase 1

2012 Beth Israel Analysis

- 100% drug library utilization in ICUs
- Only 7% corrections and only 3 insulin alerts
- 51% of overrides associated with key drugs
- Drug library adjustments made
- Practice issues identified and addressed

2013 Results

**70% REDUCTION
IN BETH ISRAEL ALERTS**




Results Phase 2

2013 System Wide Analysis

- Only 7% corrections across system
- 67% of system overrides associated with key drugs
- System wide drug library adjustments made
- System practice deviations (bolus dosing) identified and adjusted to enhance patient safety

2014 Results

**58% REDUCTION
IN SYSTEM ALERTS**



* Through April 2014

Conclusion

We significantly reduced our alerts and improved safety across our health system through the unique insight of the DoseTrac real-time monitoring, reporting and analysis service. We maintain our gold standard of care through scheduled reports, ongoing monitoring & routine data reviews with our vendor.