

## Back Door

Part Number	Name	Description	Qty.
<b>SMPS-MC3 Power Supply</b>			
34562337	Power Supply SMPS_MC3	New Power supply	1
34562338	Retrofit kit	Needed for initial change /Update	1
34562339	PCB Alarm supply	Buzzer/Capacitors	1
34562340	PCB Relay	K3/K4	1
34562343	Fan	On bottom of power supply	1
34562342	Fuse T 1.25 A (TE5)	F403	10
34562341	Fuse T 3.15 A (TE5)	F301 - F303 - F304 - F601 - F602	10
34570179	Fuse M10 A [6.3X32 MM]	F600	10
34570187	Fuse M20 A [6.3X32 MM]	F5-F6 [For Heater]	10
34570470	Fuse T 5 A (TR5)	F302 / F402	10
34570462	Fuse T 6.25 A [6.3X32]	F1-F2 [For Mains Voltage]	10
<b>SMPS-MC2 Power Supply</b>			
3461076	Switch Mode Power Supply MC [110/240V]	Power Supply (Capacitors on Top Left)	1
34570179	Fuse M10 A [6.3X32 MM]	F600	10
34570187	Fuse M20 A [6.3X32 MM]	F5-F6 [For Heater]	10
34771328	Fuse F12 A [6.3X32 MM]	F5A / F5B [For 110/120 V Heater] (Old Power Supply)	10
34570470	Fuse T 5 A (TR5)	F302 / F402	10
34570462	Fuse T 6.25 A [6.3X32]	F1-F2 [For Mains Voltage]	10
34518819	Fuse T 1.25 A (TR5)	F403	10
34570160	Fuse T 3.15 A (TR5)	F301 - F303 - F304 - F601 - F602	10
34514775	Relay [12V]	K3/K4 (12V - SMPS)	1
34562223	RIBBON CABLE 20PIN (SMPS-MC,P24)	Power Supply ----> Basic Board	1
34519912	Capacitor, 1 F	C114/C115	2
34517610	Fan 80 X 80 mm	Bottom of power supply	1
34560580	Buzzer	H104 Add	1
<b>Power cord</b>			
3451059A	Mains Switch	"On-Off" Switch for Machine	1
34770003	Protection Cap for Mains Switch	Cover for Mains Switch	1
3461069	Power Cord Plug to back door plug	Mains Cord for SN # 244242 ↑ Includes all seals and screws	1
34561064	Mains Cord	The long green ground wire to the post And the 2 connectors to SMPS	1
34561161	Mains Cord	Power Cord (USA) Plug to the SMPS Bare wires SN#244241 ↓	1
34518029	Strain Relief	For Mains Cord	1
<b>DF FILTER</b>			
34561633	Coupling Without Magnet	For DF Filter**	1
34990445	O-ring 12.37 x 2.62	O ring for Coupling Without Magnet	5
34562265	Silicone tubing 6X1.75/Blue	For DF Filter option Replaces White tubing	2m
34565256	Silicone tubing 6X1.75/Red	For DF Filter option Replaces White tubing	2m
7102102	DF Filter Kit - EVO	Kit to Add DF Filter to Version 9	1
7107368	DIACAP ULTRA DF FILTER CA & US	Replacement filter	1
<b>MISC.</b>			
3451441A	Tubing Holder / Fan Cover	Holder for Crank for Blood Pump and fan cover	1
3456130A	Fan 92 x 92 mm	Fan on Back Door	1
3457028A	Dust Filter	Filter for Fan	5
3451835A	Crank for blood pump	Handle type Old Style	1
34560157	Crank for blood pump	Peg Type	1
JH14GA	Jug Holder 14 Ga	Disinfect jug holder	1
<b>Base Housing</b>			
3451489A	SWIVEL CASTER+	Rear wheels without brake	1
3451490A	SWIVEL CASTER WITH BRAKE+	Front wheel with brake	1
7102244	OPTION ACCU DIALOG + SW9.XX	Upgrade of Ver 9 Battery Back up	1
3456012A	TUBING CONNECTOR 11XM12 +	Stainless steel connector for inlet hose ( Old Style )	1
34560130	TUBING CONNECTOR 11XM12 +	Plastic for drain hose SN 247830	1
34565140	PVC TUBING 10X3, Red Cold Water feed	Inlet water cold water only	1
34565159	PVC TUBING 10X3, BLACK	Inlet/Drain Cold water only	1
34561749	SILICONE TUBING 8X3,2 / HIGH TEMP,RED	High temp inlet water hose Red	1
3456516A	SILICONE TUBING 8X3,2 / HIGH TEMP, RED	Inlet Hightemp kit	1
P00025	Female Garden hose connector 3/8"	Inlet water hose connection	1
34562278	TUBING CONNECTOR WATER INLET	Connector for inlet W/ "C" Clip SN 247830 ↑	1
34562279	TUBING CONNECTOR WATER OUTLET	Connector for drain hose W/ "C" Clip SN 247830 ↑	1
34773770	Y-TUBING CONNECTOR 6 MM	For DF option ( TSDE or TSDE_S, top DF filter, bottom DF filter )	5
34570365	TUBING CLAMP SNP - 6 [ 11.4-13.0 mm ]	Black Clamps	20
34561781	Base Platform	Machine base for ver 9.xx	1