

## Concentrate Supply

Part Number	Name	Description	Qty.
<b>WANDS - CENTRAL &amp; NON-CENTRAL</b>			
34562260	Male Connector for Suction Rod	Coupling w/ Red-Blue Caps [Non-Central]	1
38245442	SEAL NIPPLE ZKV.	Wall Box Connector [Central]	1
34560920	Suction Coupling	Attaches to Coil Tubing [Central]	1
34570438	O-ring [8.0 x 1.25]	O-Ring for Suction Coupling [Central]	5
34560939	Suction Nipple	Attaches to PE Tube [Central]	1
34570004	O-ring [9.3 x 2.4]	O-Rings for Suction Nipple & Male connector	5
3451728A	PE Tube [4x1]	Hard Plastic Tubing for Jugs	1
34565280	Silicone tubing [3.2 x 1.6]	Attaches PE tube to Suction Nipple	3m
S20110	Delivery Coil - BICARB	Blue Coil Tubing [Central]	1
S20100	Delivery Coil - ACID	Red Coil Tubing [Central]	1
34565213	Silicone tubing [2.9 x 3]	Silicone Tubing for Wands [Non-Central]	10m
34561498	Disposable Filter [2.6 mm]	Filters for Wands	5
3451726A	Feed-Through Connector	Connector Under Front Door for Wands	1
7105196	OPTION CCS DIALOG+ Install kit	Central concentrate block attached to machine	1
<b>RETAINER TUBES</b>			
3456077A	Retainer Tubes (green) for Suction Rods	Green Tubes PE Tubing Rests In	1
34565345	Silicone Tubing [12 x 2.5]	Holds Retainer Tubes on Bracket	1m
34514120	Retainer for Suction Rods	White Housing Wand Couplings Attach To	1
34517138	Reed Sensor [w/ contacts & O-ring 3.68 x 1.78]	Sensor To Make Sure Wands are Attached	1
<b>CONTAINERS</b>			
604010	Concentrate Container 10 L NATURE 365G	Empty jug no Label	1
34514562	Suction Tubing "YELLOW" Complete	Cap. Tubing, and filter for citric jug	1
34565205	Disinfection Tubing PUR tubing 3.6 X 1.1	Tubing only	1
34517774	Canister Lid, Blue	Bicarb w/ hole	1
604032	Stopper For Canister Blue 45 OV	Solid cap	1
S20064	Label BICARB Jug	Blue	1
34517766	Canister Lid, Red	Acid w/ hole	1
S20065	Label Acid Jug	RED	1
34512675	SUCTION PIPE	Metal Filter for Citric Jug	1