

## DF SUB RACK

Part Number	Name	Description	Qty.
<b>DF Block</b>			
34562339	DF Block II	Retro fit kit for update to new DF Block	
34561765	DF Block	Block that Conductivity Sensors Attach to	1
3477103A	Tubing connector [3.5 x 1/8"]	Tubing Connectors on block for KP and BICP	2
34770852	O-Ring [11.1 x 1.6]	O-Ring for tubing connector	5
3451916A	Tubing connector [8 x 1/8"]	Tubing Connectors on block for FPE	5
34770852	O-Ring [11.1 x 1.6]	O-Ring for tubing connector	5
34770925	Plug	Cap for RVB / RVK	5
34770852	O-Ring [11.1 x 1.6]	O-Ring for Plug	5
34771050	Pressure Spring [500 MM HG]	RVB / RVK	10
34772928	Centre Disc	Disc that goes with "Pressure Spring" RVB/RVK	3
34771085	Seal	Goes under Centre Disk	5
3456102A	Bicarbonate Conductivity Sensor Vr.2	BICLF [with temperature sensor]	1
3456103A	END Conductivity Sensor Vr.2	ENDLF [with temperature sensor]	1
34770950	O-Ring [8.0 x 1.6]	O-Ring for BICLF / ENDLF	20
3451614B	Temperature Sensor	TSBIC / TSD / TSD-S / TSDE/TSDE_S/TSE	1
34562060	Temperature Sensor inc. Wire	Temperature sensor with wires attached	1
34570373	Plug (Long) TSBIC/TSD	Plug holds sensor inside housing of the Conductivity cell	10
34570381	PLUG (SHORT) TSD_S/TSDE/TSDE_S/TSE	Plug holds sensor inside housing	10
3477104A	O-Ring [8.3 x 2.4]	O-Ring for TSDS	10
<b>GEAR PUMPS</b>			
34560718	Motor For Gear Pump (SW 1.7.4)	Motor For EP, FPE and FPA pumps	1
34562256	Pump Cover	White plastic cover FPE,FPA,EP	1
7704130	Assembly Tool	Tool for setting gap when replacing gear pump magnets	1
34562310	Retrofit kit Gear pump	To replace the magnets	10
34562330	Gear pump FPA	New FPA Starting SN # 280707 ↑	1
34562331	Gear pump FPE	New FPE Starting SN # 280707 ↑	1
7704134	Distance Gauge Gear pump	Tool for installing New Gear Pumps	1
34561722	Pressure Controller [1.3bar]	RVFPE/RVFPA Throttle Valves	1
34570683	O-Ring [4.47 x 1.78]	O-Ring for RVFPE	10
34514988	Non-Return Valve	RVUFP	1
<b>PISTON PUMPS</b>			
3456055B	Motor for Piston Pump	BICP - KP - UFP	1
34560602	Pump Body for Piston Pump	BICP - KP - UFP	1
7703848	Assembly Jig	Spacing Tool for Assembly of Pumps	1
34519246	Reed Sensor	For Bell Joint With [1 Magnets] Serial Number < 100413/ 150224	1
34560572	Bell Joint for Reed Sensor	BICP - KP - UFP [1 Magnets] Silver	1
34561480	Hall Sensor	For Bell Joint With [2 Magnets] Serial Number > 100414/ 150225	1
34561471	Bell Joint for Hall Sensor	BICP - KP - UFP [2 Magnets] Black	1
34570799	Vibration Damper	For KP, BICP, and UF pump	1
34570861	T-TUBING CONNECTOR 6-4-6 MM	(UFP,FPA,LA)	1
34570705	Y-TUBING CONNECTOR 8 MM	(WA,KP,BICP) (WA,BICP,KP) Bypass for pumps	1
3456222A	Wire Harness Valves and pumps		1