

Front Door

Part Number	Name	Description	Qty.
BLOOD PUMP			
34562319	Blood Pump	With cover and seal But NO ROLLER	1
34562318	Blood Pump Roller	Not backward compatible	1
34562317	Cover for Blood Pump	Not backward compatible	1
3456079A	Blood Pump Without Roller	Comes w/Door Cover	1
3456080A	Roller for Blood Pump [8x12]		1
34570616	Latch	Blood Pump Roller Release	1
3451669F	Cover for Blood Pump	Door for Blood Pump (no hinge)	1
3451981F	Cover for Blood Pump with Pins	Door for Blood Pump (with hinge)	1
34519823	Hinge Pin with Internal Thread	Hinge For Cover	1
3456003B	Metal Disc [10mm]	For Blood Pump Cover	1
34560025	Magnet [D5]	Blood Pump Door Magnet	1
34517138	Reed Sensor [w/ Contact & O-Ring 3.68 x 1.78]	Sensor for When Cover is Closed	1
34561358	Direction of Rotation Board for Blood Pump		1
34516670	Seal for blood pump		1
HEPARIN PUMP			
3451432D	Heparin Pump Compact	Heparin Pump	1
34514732	Syringe Holder	Hook to Hold Syringe With Spring and "c" clip	1
34570667	Clamp	Clip for Heparin Pump (shaped like a "W")	1
34570659	Unlocking Lever Set	Lock for Heparin Pump	1
34562226	RIBBON CABLE 20PIN (HEP,P6)	Ribbon from Heparin pump to Basic Board	1
34570640	Dummy plug	Plugs to cover holes in HP for the syringe holder	1
PRESSURE SENSORS & BOARDS			
34570845	Manometer Connection, BLUE	PV	5
34570853	Manometer Connection, RED	PA / PBE	5
34516409	Filter	Internal Water Transducers PV/PA	1
34565299	Silicon tubing 3 x 2 mm	Manometer connection to the Filter	3m
34562287	FKM Tubing 1.4X1.5	Replaces Silicone tubing between Filter and Board	5m
34561641	Pressure Monitoring Module PA/PV	Evo II (No PBE) "1 red Manometer"	1
34562224	Ribbon Cable 10-Pin	LR- P11 on Basic Board	1
3456081B	Valve Block Level Regulation (SP)	Evo I (Has PBE) "2 red Manometer"	1
34560823	Diaphragm Pump	Evo I	1
34514107	Rotary Magnet (Currentless Closed)+	SAKV Clamp	1
3456005E	SAD with Red Detector, VR. 4	SAD New RVD started SN # 254760 ↑	1
34561803	SAD cover	Cover Behind Front Door	1
34562228	Ribbon Cable for SAD	SAD --> Basic Board	1
HOLDERS & Clamps			
S8020004	Dialyzer Holder	Attaches to IV pole	1
SLDPLATE	Streamline Mounting Bracket	Metal plate for Dialyzer Holder	1
NC10-0527A	Dialyzer Holder Bracket for Streamline	Plastic Holder - Attaches to Mounting Bracket	1
34570543	Double Tubing Holder [6.8 mm]	Above Blood Pump Left side	1
34570519	Single Holder, big	Above Drip Chamber	1
34570527	Single Holder, small	Drip Chamber Holder / Above S.A.D.	1
34570560	Single Tubing holder [6 mm]	Small Tubing Holder Under Blood Pump [vertical]	1
34570535	Single Tubing Holder [6.8 mm]	Small Tubing Holder Above Blood Pump [Horizontal]	1
34514317	Cover with Tubing Holder /+	Under SAD [SAKV-SG]	1
34570896	Cover cap for SAKV	Cap for	20
MISC. Exterior			
3451415A	Cover Complete+	Suction Rod Cover	1
3451489A	Swivel caster	Wheel w/o Break [Push In]	1
3451490A	Swivel caster with brake	Wheel w/ Break [Push in]	1
34561595	Front Pannel, Bottom	2 Holes centered	1
34570772	Countersunk screw [M4x30 (Torx)]	Screws for Screen and Front Panel	6
3451962C	Infusion Pole	IV Pole	1
3451726A	Feed-Through Connector	Connectors for under door for concentrate wands	1
3451439A	Holder, Top	Hold top of IV Pole	1
3451438A	Holder, Bottom	Holds Bottom of IV Pole (machine without Solcart)	1
7704090	Door Latch	For holding the door open	1
S20046BH	Priming Bucket	Bucket on side of machine	1
38540009	Front Pannel, Bottom	4 Holes for the filter connections	1
38540008	Cross member Bracket for filter connection	Goes above the front panel bottom for filters	1