

## Monitor

Part Number	Name	Description	Qty.
<b>BOARDS &amp; MONITOR COMPONENTS</b>			
3461013B	TFT Display 15"	Display Monitor	1
34570772	Countersunk Screw M4X30	Screws for Monitor and Panel Bottom (Torx 15)	6
3451429A	Optical Status Display Board OSD	Green, Red, and Yellow LED's on side of monitor	1
3457069A	OSD Cover	Cover for OSD board, Comes with seal and screws	2
34562184	Backlight Inverter Board NEW	Serial number 105218 & ↑. Anything ↓ needs Retrofit Kit	1
34561757	Retrofit Kit TFT Display G150XGE	For Retrofitting to Chimei TFT Display (1x backlight inverter board 34562184) (cables, brackets Screws and cleaning cloth) Must Also Order: 3461013B	1
34562282	Front Panel Board (FPB)+ Ver 3	FPB SN# 245466 & higher Board with intergrated BIB	1
34562283	RETROFIT KIT FPB	Retrofit FPB Ver 3 ( included ) to replace Ver 1 & 2 FPB	1
34562299	SERV.KIT VI - FPB BOARD		
34514260	Touch Screen for LCD+	Glass Touch Screen	1
34515755	Seal Cord 1100 mm	Seal for touch screen/ Monitor	1
3451460A	Touch Controller Board TCB Version 3	Small Board touch screen ribbon plugs into.	1
3451491B	FRONT PANEL COMPLETE WITH MEMBRANE	Bezel with Membrane	1
<b>RIBBON CABLES &amp; CONNECTORS</b>			
34560971	Ribbon Cable 10 Conductor (FPB-BIB)	Front Panel Board > Backlight Inverter Board	1
34560890	Ribbon Cable 10 Conductors (Motherboard-TCB)	Touch Control Board > Motherboard	1
34560912	Ribbon Cable 12 Conductors (FPB-OSD)	Front Panel Board > OSD	1
34562225	RIBBON CABLE 26PIN (FPB,MB)	Front Panel Board > Motherboard	1
34562222	RIBBON CABLE 16PIN (FPB,P14)	Front Panel Board ----> Basic Board	1
34562297	Monitor Cable BSW MB	LVDS > LCD Video cable	1
34560998	Ribbon Cable 20 Conductors (LVDS-LCD)	LVDS > LCD Old video cable	1
34561285	Membrane	Label under screen that covers Blood pump +,-	1
34561366	Membrane for ADIMEA Option	Label under screen that covers Blood pump +,-	1
34562276	Membrane From SW 9.1A	Membrane with the Green + on the sticker	1
7105238	Multi functional tray	Flat Try behind Monitor	1
39174778	BOX STACKABLE FOR DIALOG PLUS	Bucket Style Tray behind monitor	1