

Rinse Bridge

Part Number	Name	Description	Qty.
Version 4 or 5 compatible Rinse Bridge Components			
34561382	Dialyzer coupling, Red	Red Dialyzer Coupling	1
34990445	O-ring 12.37 X 2.62	O-ring for Dialyzer Coupling	5
34561501	Disposable filter 4.6 mm	Filter on Red Dialyzed Coupling Tubing	5
34561390	Dialyzer coupling, Blue	Blue Dialyzer coupling	1
34990445	O-ring 12.37 X 2.62	O-ring for Dialyzer Coupling	5
7102867	Sample Port Row	Sample Port on Blue Dialyzer Coupling Tubing Push on connection	1
3456173A	Sample Port	Luer Loc connection	1
34773762	O-ring 9.25 x 1.78	O-Ring for Metal Connectors for Dializer couplings	10
34565191	Silicone Tubing 5 X 3	Tubing for Dialyzer Couplings	10 meters
34561692	Blood Leak Detector V4	BLD	1
3451898B	Repair Kit Disinfection Valve	For the Complete Replacement of all Small Internal Parts for the Disinfection Valve (2 x O-ring 2.7 x 1.5) (2 x O-ring 5 x 1.6) (3 x O-ring 6.4 x 1.3) (2 x Cover) (4 x Countersunk screw M3x8) (2 x Valve tappet) (2 x Pressure spring) (1 x Eccentric roller) (1 x Eccentric shaft)	1 kit
34517324	Valve Tappet	For VD	2
34570012	O-ring 2.7 x 1.5	For Valve Tappet	10
34517332	Pressure Spring	For Valve Tappet	2
34560840	Eccentric Shaft	Shaft that VD turns to open	1
34570457	O-ring 5.0 x 1.6	For Eccentric Shaft	10
34560858	Eccentric roller	Roller on eccentric shaft that moves tappets	1
34570870	Flag for Servo Motor	Grey Flag for Ver. 4 Rinse Bridge	1
3451660A	Fork Light Barrier Supervisor sensor	For Servo Motor (Light Barrier Sensor)	1
34514562	SUCTION TUBING "YELLOW" COMPLETE	Yellow Lid, Hose and Metal Filter for disinfection Jug	1
34565272	PVC tubing 3.5 x 1	Tubing for Disinfection Jug	10m
34512675	Suction Pipe	Metal Filter on end of hose going into Citric Jug	1
3477106A	O-ring 6.4 x 1.3	For Cover on Valve Tappet	5
34517138	Reed Sensor with Contact and O-Ring	Sensor for Dialyzer Coupling	1
3456124A	Pressure Sensor	PDA	1
34773762	O-ring 9.25 x 1.78	O-ring for PDA	10
34562311	Digital Servomotor	Works with Digital Servomotor Board ONLY (SN 267600 ↑)	1
34562312	Digital Servomotor Board	Works with Digital Servomotor ONLY (SN 267600 ↑)	1
3451902C	Solenoid Valve [connector/flange]	VDE	1
34770852	O-ring 11.1 x 1.6	O-ring for VDE	5
Version 5 Specific Rinse Bridge [Oval] SN 201787 ↑			
34561846	Rinsing Bridge Housing Version 5 incl. seal	Oval Plastic housing only	1
34561862	Tubing connector (male)	Dialyzer connection " L shape " (Incl. O-ring 12.5 x 1.5 and O-ring 5.6 x 1.8)	1
34562010	Plastic Holder	Mount for VDE & PDA	1
34561994	Valve housing	Grey Housing That Tappets are placed in	1
34561986	Valve housing lid with o-ring 6.4 x 1.3 Ver. 5 only	Lid for Grey Housing Cap for Tappet	1
34562001	Servo Motor Holder	White Plastic Bracket that Holds VD to Rinse Bridge	1
34561854	TUBING CONNECTOR (DISINFECTION PORT)	Port on rinse bridge jug tubing connects to	1
34565205	PUR TUBING 3.6X1.1 MM	Tubing between tappets and Citric jug port internal	1
Version 4 Specific Rinse Bridge [Rectangular] SN 201786 ↓			
3451462A	Rinsing Bridge, Housing V4	Rectangular housing only	1
3477103A	Tubing Connector 3.5 x 1/8" [Plastic]	Rinse Bridge Plastic threaded Feed Through for Vent Tube	2