

Tools & Accessories

Part Number	Name	Description	Qty.
Tools			
77091B	Software 9.1B	SW. 9.1B	1
77091E	Software 9.1E	SW. 9.1E	1
34562274	Retrofit kit 9.1B/9.1E upgrade		1
7704134	Distance Gauge Gear pump	Tool for installing New Gear Pumps	1
7703864	Field Service Utility V 2.00	FSU V 2.00	1
7701756	Extraction Tool Set	Remove Pins from Connectors	1
7704102	TUBEING EXTRACTING TOOL	For removing tubing off connectors	1
7703848	Assembly Jig	Spacing Tool for Piston Pumps	1
7702736	Extraction Tool for Hinge Pin	Remove blood pump cover pins	1
7703627	Touch-Up Pen Light Grey+	Touch-up Paint For Bezel	1
7703619	Touch-Up Pen Pale Green+	Touch-up Paint For machine	1
34561811	<p style="text-align: center;"><u>Wire repair kit, consisting of</u></p> [10 x 30 cm wire 0.25 mm ² with MMNL ST/BU 0.25male/female (tin-plated)] [10 x 30 cm wire 0.25 mm ² with MMNL ST/BU 0.25male/female (gold-plated)] [10 x 30 cm wire 0.5 mm ² with MMNL ST/BU 0.25male/female (tin-plated)] [10 x 30 cm wire 0.25 mm ² with locking clip both female(tin-/gold-plated)] [10 x 30 cm wire 0.5 mm ² with MODU 4 both female(tin-/gold-plated)] [25 x 3M™ Scotchlok™ wire connector UB2A]	Used to Replace any Broken Pins or Wires	1
7701497	SPECIAL SPANNER M18	For tightening the locking ring on DCI Cable	1
7703708	Grease S.K.D. 4002 (100G)	Blood Pump Shaft Lube	1
7702345	Silicone Grease 410 Medium 100G Tube	O-Ring Lube	1
7702353	Silicone Rubber Adhesive RTV 162Q White	RTV Adhesive for Boards	1
PM's			
770085A	FLOWMETER 1,1 LTR/MIN	B Braun flow Meter	1
7704104	CONNECTION SET DIG.MANOMETER/FLOWMETER	Includes B Braun Digital Manometer & all connections	1
7704106	B Braun DIGITAL MANOMETER	Pressure meter Only	1
34562307	Maintenance kit NEW	Pm Kit Including Tubing	1
7703090	Expansion Chamber 500 ML	PM Tool ABPM	1
770203A	Five-Fold Stopcock Incl. Tubing's	PM Tool ABPM	1
7702493	Template 1.4 / 1.5 MM	PM Tool SAKV gap check "Go/No Go"	1
7703805	Male to Male Connector	Adapter to connect both dialyzer couplings for the measurement of the patient leakage current	1
34565191	Silicone Tubing 5 X 3 (10m)	Tubing for test jigs	10m Bag
34570705	Y-TUBING CONNECTOR 8 MM	For test jig	1
CALIBRATION			
34511571	Manometer 0-2,5 BAR	DMV Check/Calibration	1
3451916A	Tubing connector [8 x 1/8"]	Connecting customer meter to DMV to calibrate	5
7701357	Manometer 0-4 BAR	Big gauge with cert.	1
7700911	Reference Solution	Blood Leak Calibration Solution	1
7704090	Door Latch	For holding front door open	1