

TOP & LOW LEVEL

Part Number	Name	Description	Qty.
BOARDS			
34560610	LX 800 Motherboard	TLC	1
34562301	BSW Motherboard	SW 9.1E ↑ (SN# 266394 and ↑)	1
34562295	SERV.KIT II - BSW	Upgrade kit	1
34562302	mSATA mini card	SW 9.1E ↑ With BSW moterboard New	1
34519327	Battery 2032 / 3 Volt	Battery for Mother Board	1
3456088A	Compact Flash Card, 1 GB	CF Card ONLY	1
34562168	Basic Board V3	LLBB	1
34562150	Power Board (Valves)	PBV	1
34560750	Power Board (Motors)	PBM	1
34562303	Digital Board Vr.6 (LLD-III-V5)	SW 9.1E ↑ New (266394 and ↑)	1
34560734	Analog Board	AB	1
3451487E	ABPM Module M2015	ABPM	1
3456127A	Temperature Switch	Turns on Fan When Hot	1
RIBBON CABLES			
34560890	Ribbon cable 10 conductors	Motherboard → TCB	1
34562225	Ribbon cable 26 conductors	Motherboard → FPB	1
34560912	RIBBON CABLE 12 CONDUCTORS (FPB-OSD)	Front Panel Board ----> OSD	1
34560971	RIBBON CABLE (FPB-BIB)	Front Panel Board ----> Backlight inverter Board	1
34562222	RIBBON CABLE 16PIN (FPB,P14)	Front Panel Board ----> Basic Board	1
34562228	Ribbon cable for SAD	Basic Board ----> Saftey Air Detect	1
34562224	Ribbon cable level regulation	Level Regulation ----> Basic Board	1
34562223	RIBBON CABLE 20PIN (SMPS-MC,P24)	Power Supply ---> Basic Board	1
34562297	MONITOR CABLE BSW	LVDS to the LCD Display (266394 and ↑)	1
34562296	SERV.KIT III - HWBSL & DCI OPTION	Adapter for the new BSW Motherboard	1
COMMUNICATIONS			
7102231	WAN - BSL sw. 9.XX	Complete Option for BSW mother board	1
7102230	WAN - BSL Nexadia sw. 9.XX	Complete Option for BSW mother board	1
34562324	Cable BSL for LX800	Adapter cable from BSW to LX800	1
34514554	Bedside Link Board / BSL	BSL [DaVita]	1
7107232	EXTERNAL MFC-LAN CONNECTION CABLE	BSL External Lan Cable [Davita]	1
7107218	DCI Cable - All	DCI Includes adapter cable for BSW	1
770284A	DCI cable, external	External portion of DCI cable	1
34514821	DCI cable, Internal	Internal connected to Motherboard	1
710725A	ReBell Site License	For DCI	1
7107251	Server SW ReBell Tool	For DCI	1
34560548	Card Reader Board	Replacement board (Not for HID Card Reader)	1
34562321	Cable for LX800 Card Reader	Adapter cable for Card Reader from BSW to LX800	1
7105230	HID Card Reader Option for BSW MB	SW 9.1E ↑ with BSW MB	1
34562322	Cable for LX800 Adimea	Adapter cable for Adimea from BSW to LX800	1
34562300	HID Card Reader Board for BSW MB	Replacement board for HID card reader	1
34562294	SERV.KIT I - HID CARD READER		1
7105232	Card Pack	Pack of 5 Patient ID Cards	5
MISC.			
34514929	Loudspeaker+		1