

WATER SUB RACK

Part Number	Name	Description	Qty.
WATER TANK			
34561315	Water Block Without Lid	WAB (Upline Tank)	1
34561323	Lid for Water Block	WAB (Upline Tank)	1
34560459	Level Sensor Complete	NSVB	1
34560246	Float Ball Incl. O-Ring	Sub Component for NSVB	1
3456021B	Heater 120V / 2x900 W	Heater element [incl. seal]	1
34560467	Seal	For Heater	1
EP PUMP			
34560718	Motor for Gear Pump (SW 1.7.4)	Motor For EP, FPE and FPA pumps	1
34562256	Pump Cover	White plastic cover FPE,FPA,EP	1
34562332	Gear Pump EP Higher than SN# 200,000	EP Pump for SN # 200,000 ↑	1
3461097	Gear Pump EP Lower than SN# 200,000	EP Pump for SN # 200,000 ↓	1
7704134	Distance Gauge Gear pump	Tool for installing Gear pump	1
VALVES			
3451906C	Solenoid Valve 2/2 [Connector/Connector]	VVBE-VZ-VBICP-VBKS-VBKO-VVB	1
34561609	Solenoid Valve 2/2 [With Bypass]	VEB (Degassing Valve only)	1
34770852	O-Ring 11.1 x 1.6	O-Rings for Valves	5
34561706	Pressure Controller 0.4 bar	RVDA	1
PRESSURE SENSOR AND REGULATOR			
3456124A	Pressure Sensor	PE Sensor	1
34773762	O-ring [9.25 x 1.78]	O-ring for PE Sensor	10
33193053	HOUSING LOCKING CLIP	Black Electrical connector (Plug)	1
3456168B	Pressure Reducer Valve	DMV	1
34517243	Tubing connector [7 x 1/8"]	For Top and Bottom of DMV	1
3451357A	Stop Screw	Plug For DMV	1
TEMPRATURE SENSORS			
34570730	Housing Complete for Temperature Sensor	Grey housing for TSE (Now white)	1
34570381	PLUG (SHORT) TSD_S/TSDE/TSDE_S/TSE	Plug holds sensor inside housing	10
3451614B	Temperature Sensor	TSBIC / TSD / TSD-S / TSDE/TSDE_S/TSE	1
3477104A	O-ring [8.3 x 2.4]	O-ring for Temperature Sensor	1
34562060	TEMP.SENSORCABLE INC. SEN. SW 9.XX	All Sensors	1
DBK			
34562219	Housing DBK	Block For DBK	1
34771050	Pressure Spring [500 MM HG]	Spring For DBK	10
34771352	Centre Disc (optional)	Disc For DBK	3
34570888	Centre Disc	Disc For DBK (White disk)	3
34770925	Plug	Cap For DBK	5
3451916A	Tubing connector [8 x 1/8"]	Connectors on DBK	5
34770852	O-RING [11.1 x 1.6]	O-Ring for all of DBK	5
MISC			
34561501	Disposable Filter 4.6mm	FEP, FVE	5
34570195	TUBING CLAMP SNP - 1.5 [7.8 - 8.6 mm]	Black Clamps	20
34570209	TUBING CLAMP SNP - 4 [10.3 -11.5 mm]	Black Clamps	20
34570365	TUBING CLAMP SNP - 6 [11.4-13.0 mm]	Black Clamps	20
34570705	Y-TUBING CONNECTOR 8 MM	(VBICP,VVB,DF manifold)	1
34570861	T-TUBING CONNECTOR 6-4-6 MM	(VBKS,VBKO & top of SolCart)	1
7102233	OPTION ADIMEA SW9.XX	With BSW Motherboard connection	1
34562322	Cable for LX800 Adimea	Adapter cable from BSW to LX800 mother board	1
3456063D	KT / V-UV BOARD	Adimea Spare Part	1
34565256	Silicone Tubing [6.0x1.75/red]	Connection to Water Supply (DMV-> VVBE, DMV-> Base)	2m
34570918	Single Ear Clamp [14.0/13.3]	Used to Secure Red Silicone Tubing	10
3456222A	Wire Harness Valves and pumps	For Pumps and Valves	1
38110005	Sensor Harness	PE and NSVB	1