

CARESITE® Luer Access Device

Designed to make access clear and easy for improved patient care

CARESITE® Luer Access Device (LAD)

Clean

Easy to:

Visually Inspect






Clear design helps verify proper flush.

Smooth surfaces permit easy cleaning.

- LAD and extension sets are acceptable for high-pressure use with power injectors.
- For use with power injection.
 LAD 400 psi and 15 mL/sec.
 Extension Sets 300 psi and 10 mL/sec.



Grip

Easy grip barrel designed to minimize slips and touch contamination.

- Available in Y-access site for full IV administration set integration.
- Integrated fluid path designed with minimal dead space.



Connect

Open fluid path designed to minimize areas where bacteria can collect.



CARESITE Performance Specifications	
Flow rate	161-231 mL/min (10-12 L/hr)
Pressure, Power injection	400 psi and 15mL/sec (LAD only)
Retained volume	0.20 mL
Priming volume	0.22 mL
Fluid displacement	0.03 mL (positive)
Durability (# accesses)	216 min
Swabbable	Yes
Microbial barrier (ingress)	Yes
Not made with latex	Yes
Blood compatible	Yes
Lipid compatible	Yes
Luer slip fit	Yes
Materials	Polycarbonate & Silicone
No metal components	Yes
Flow path	Open fluid path around plunger No fluid entrapment areas
CARESITE Ordering Information	

CARESTE Ordering information	
Product order number	415122
System offering	Individual pkgd. unit
Case qty.	200 pkgd. units/case

B. Braun Medical Inc. | Bethlehem, PA 1-800-227-2862 | www.caresite.BBraunUSA.com