

Backed by a Responsive Service and Support Program

When you work with B. Braun, you have access to a comprehensive service and support organization that is staffed with an experienced and responsive team of specialists. The program consists of in-service and customized training for pharmacists and technicians, continuing education opportunities, and immediate 24/7 access to clinical and technical assistance.

For more information, or to schedule a demonstration of the PINNACLE® TPN Management System, call your B. Braun representative at **1-800-227-2862**. For Clinical and Technical Support, call **1-800-854-6851**.

PINNACLE® TPN Management System Compounding Accessories

Streamline Your TPN Process



Bags, Sets, and Accessories designed for safety and efficiency







Accessories for the PINNACLE® TPN Management System



High Flow Y-Adapter

- May be used for any fluid transfer, including lipids
- On/Off clamps

REF No. 2112431

20 per case



Dual-Chamber Bag Fill Set

Attaches to transfer sets for compounding dual-chamber bags

REF No. 2112415

25 per case



Micro-Spike Adapter

Attach to transfer set lead for compounding from glass vials 100 mL or smaller without any leakage.

REF No. 418121 50 per case



Tamper Resistant Clamps

Plastic hand sealer with permanent clamps

For TPN bags 500 per case 1 liter or larger REF No. 2112397

For TPN bags 500 per case

smaller than 1 liter REF No. 2112432



Control Panel Mounts

Provides flexibility within the clean room by allowing the PINNACLE Control Panel to be mounted on the wall or beam next to the laminar flow hood.

Wall Mount One per case REF No. 126909

Swing Arm Mount One per case REF No. 401539



TPN Bag Stand Holder

Assists the Technician or Pharmacist to provide stability to the TPN Bag while adding solutions via a syringe; provides needleless access.

REF No. 435251 One per case



EVA Dual-Chamber Bags

- Easy mixing for homecare patients remove green bar and shake bag
- Needleless access for fluid additions
- Top chamber for lipids; bottom chamber for clear solutions



EVA 3-in-1 Mixing Containers

- Designed for safety and durability
- Needleless access for fluid additions



Transfer Sets

- Available in 6 and 9 lead sets
- Components do not contain DEHP or rubber tree latex

1500 mL Bag Top Chamber: 250 mL Bottom Chamber: 1250	28 per case	250 mL Bag REF No. 2112346	25 per case	Universal Spike 9-Lead REF No. 2112342	10 per case
REF No. 2112406		500 mL Bag	25 per case		
		REF No. 2112347		Universal Spike 6-Lead	10 per case
3000 mL Bag Top Chamber: 500 mL	28 per case	1000 mL Bag REF No. 2112348	40 per case	REF No. 2112341	
Bottom Chamber: 2500 mL REF No. 2112407				Non-Vented	10 per case
		2000 mL Bag REF No. 2112349	40 per case	9-Lead REF No. 2112344	
4000 mL Bag Top Chamber: 500 mL Bottom Chamber: 3500 REF No. 2112408	28 per case	3000 mL Bag REF No. 2112350	40 per case	Non-Vented 6-Lead REF No. 2112343	10 per case
		4000 mL Bag	20 per case		

REF No. 2112351