PINNACLE® TPN Management System

Specifications	
PINNACLE Catalog Number Dimensions (W x H x D)	PC1000
Control Panel	14.76 in. x 11.77 in. x 1.93 in.
Pump Module	28 in. x 19.4 in. x 11.8 in.
Weight	
Control Panel	10.78 lbs.
Pump Module	30 lbs.
Delivery Range	5 mL-5000 mL

Transfer Sets		
REF No.	Description	Units per Case
2112342	9-Lead Vented Set	10
2112341	6-Lead Vented Set	10
2112344	9-Lead Non-Vented Set	10
2112343	6-Lead Non-Vented Set	10

EVA* 3-in-1	Mixing Conta	niners — TPN Bags
REF No.	Size (mL)	Units per Case
2112346	250	25
2112347	500	25
2112348	1000	40
2112349	2000	40
2112350	3000	40
2112351	4000	20

^{*}EVA (Ethyl Vinyl Acetate) does not contain PVC or DEHP

EVA* Dual Chamber Mixing Containers RFF No. Size (ml.) Units

REF No.	Size (mL)	Units per Case
2112406	1500 mL	28
2112407	3000 mL	28
2112408	4000 mL	28

Accessories		
REF No.	Description	Units per Case
418121	Vented Micro Spike Adapte	r 50
2112397	Tamper-Resistant Clamps	500
	for TPN Bags 1 Liter or larger	
2112432	Tamper-Resistant Clamps	500
	for TPN Bags smaller than 1 Lite	er
2112415	Dual Chamber Bag Filling Set	25
2112431	High Flow Y Adaptor Set	20
	,	
435251	TPN Bag Stand Holder	1
126909	LCD Wall Mount	1
401539	LCD Swing Arm Mount	1

Schedule a demo today.

Call your B. Braun Representative at 1–800–227–2862, or visit us at www.tpnmanager.com.





www.bbraunusa.com

PINNACLE® TPN Management System

Control Panel and Compounder Features




Multi-media training and real-time help tips make the PINNACLE TPN Management System easy to learn and easy to operate



PINNACLE® TPN Management System

The PINNACLE TPN Management system streamlines all functional areas of operation and manages all aspects of the TPN process for improved efficiency and safety.



Saves Time


- Tracks NDC number, lot number and Expiration Date for each solution
- Eliminates need to input specific gravity, viscosity, and spike type
- Pooled solution calculation capability
- Networks multiple compounders

Easy to Operate

- Windows®-based program with drop-down menus and solution categories
- No specific gravity or viscosity input required
- Easy-to-follow tutorial
- Simple graphic interface
- Compounds TPN, Cardioplegia and Dialysis solutions
- Customized queues for neonatal, pediatric, adult and STAT orders

Enhances Safety

- Fast barcoded patient prescription option avoids manual entry errors
- Proven gravitational technology measures individual ingredients and final volume
- User identification required by barcoded ID badge or password



Efficient

- Fills one liter in less than a minute
- Preassembled sets for easy set-up in less than 15 minutes
- No need to prime lines
- Source containers will run until empty
- Six lead and nine lead vented and non-vented sets available

Components do not contain DEHP or natural rubber latex.

TPN Stand holder

Add additives needle-free via syringe to the mating line of the Pinnacle bags

