

Solcart B™

Bicarbonate Cartridge for use with
The Dialog⁺ Hemodialysis System



- Easy to Pierce
- Safe and Efficient
- Easy to Use
- Convenient to Store

Specifications:

Cartridge housing material: polypropylene
Solcart B cartridge volume: 760 grams sodium bicarbonate powder.

Piercing Point: Mesh filter for easy piercing



| Dialysate Flow Rate | Run Time with Solcart B* |
|---------------------|--------------------------|
|---------------------|--------------------------|

| | |
|--------------|----------|
| 300 (ml/min) | 12.2 hrs |
| 400 (ml/min) | 9.5 hrs |
| 500 (ml/min) | 7.3 hrs |
| 600 (ml/min) | 6.1 hrs |
| 700 (ml/min) | 5.2 hrs |
| 800 (ml/min) | 4.8 hrs |

* at Bicarbonate setting of 35 mmol/L

Reorder Information:

Item # 803 10 units/carton

Card Reader for Dialog⁺ Hemodialysis System

State-of-the-art technology for the storage
of patient information.



The Card Reader provided with Patient Therapy Cards efficiently and easily stores and reads the patient's data by simply inserting the Patient Therapy Card into the Card Reader.



The following data can be stored on the Patient Therapy Card –

- ❖ Patient's name and ID
- ❖ Patient's parameters
- ❖ Up to 50 Kt/V tables
- ❖ Up to 50 Crit-Line tables

The Patient Therapy Card can be used immediately as no formatting is required. As soon as a Patient Therapy Card is inserted in the card reader slot, all data is automatically read for their dialysis treatment.

For more information, please contact your B. Braun Renal Therapies Division sales representative or call B. Braun Customer Service at 800-848-2066.