








B. Braun Infusomat™ Space Pump Sets


Integrate safety across the continuum of care


Description	Product Code	REF Number	Units Per cs
Infusomat Space Pump IV Administration Set with Microbore Tubing for Epidural Administration			
	490038	490038	24
<p>(60 drops/mL) Universal spike, slide clamp, backcheck valve, two slide clamps, SPIN LOCK® connector. DEHP-free, Latex-free. Priming volume: 13 mL, Length: 127 in. (323 cm)</p>			
Infusomat Space Pump IV Administration Set with ULTRABLOCK™ UV Resistant Tubing			
	490040	490040	24
<p>(15 drops/mL) Universal spike, slide clamp, roller clamp, SUPOR® 0.2 micron air-eliminating filter, clamp, ULTRASITE injection site, SPIN LOCK connector. DEHP-free, Latex-free. Priming volume: 23 mL, Length: 119 in. (302 cm)</p>			
Infusomat Space Pump Filtered Primary Additive IV Sets			
	490041	490041	24
<p>(15 drops/mL) Universal spike, slide clamp, backcheck valve, two ULTRASITE injection sites, slide clamp, roller clamp, SUPOR 0.2 micron air-eliminating filter, slide clamp, SPIN LOCK connector. DEHP-free, Latex-free. Priming volume: 19 mL, Length: 120 in. (305 cm)</p>			
	490033	490033	24
<p>(15 drops/mL) Universal spike, slide clamp, backcheck valve, two ULTRASITE injection sites, slide clamp, roller clamp, SUPOR 1.2 micron air-eliminating filter, slide clamp, SPIN LOCK connector. DEHP-free, Latex-free. Priming volume: 22 mL, Length: 120 in. (305 cm)</p>			
Infusomat Space Pump Burette Set without Automatic Shutoff			
	490034	490034	20
<p>(60 drops/mL) Non-vented spike, roller clamp, 150 mL burette with 15 micron air guard filter and ULTRASITE injection site, slide clamp, two ULTRASITE injection sites, slide clamp, roller clamp, slide clamp, SPIN LOCK connector. DEHP-free, Latex-free. Priming volume: 19 mL, Length: 120 in. (305 cm)</p>			
Infusomat Space Pump Y-Type Blood Set			
	490043	490043	24
<p>Two non-vented spikes, two slide clamps, 170 micron blood filter, slide clamp, roller clamp, ULTRASITE injection site, SPIN LOCK connector. DEHP-free, Latex-free. Priming volume: 39 mL, Length: 120 in. (305 cm)</p>			


B. Braun Infusomat™ Space Pump Sets


Description	Product Code	REF Number	Units Per cs
Infusomat Space Pump Add-On Set			
			
Slide clamp, roller clamp, ULTRASITE® injection site, SPIN LOCK® connector. DEHP-free, Latex-free. Priming volume: 13 mL, Length: 88 in. (224 cm)	490030	490030	24

Infusomat Space Pump IV Administration Sets


			
(15 drops/mL) Universal spike, slide clamp, roller clamp, ULTRASITE injection site, SPIN LOCK connector. DEHP-free, Latex-free. Priming volume: 23 mL, Length: 120 in. (305 cm)	490031	490031	24

			
(15 drops/mL) Universal spike, slide clamp, backcheck valve, three ULTRASITE injection sites, slide clamp, roller clamp, SPIN LOCK connector. DEHP-free, Latex-free. Priming volume: 22 mL, Length: 120 in. (305 cm)	490032	490032	24

			
(15 drops/mL) Universal spike, slide clamp, backcheck valve, two ULTRASITE injection sites, slide clamp, roller clamp, SPIN LOCK connector. DEHP-free, Latex-free. Priming volume: 22 mL, Length: 120 in. (305 cm)	490036	490036	24

			
(60 drops/mL) Universal spike, slide clamp, backcheck valve, two ULTRASITE injection sites, slide clamp, roller clamp, SPIN LOCK connector. DEHP-free, Latex-free. Priming volume: 21 mL, Length: 115 in. (292 cm)	490042	490042	24

Infusomat Space Pump Low Adsorption Set

			
(60 drops/mL) Universal spike, slide clamp, roller clamp, SPIN LOCK connector. DEHP-free, Latex-free. Priming volume: 25 mL, Length: 136 in. (345 cm)	490037	490037	24

For more information, please contact your B. Braun Sales Representative or Customer Service at **1-800-227-2862**.

B | BRAUN
SHARING EXPERTISE

www.bbraunusa.com

Supor is a registered trademark of Pall Corporation.

Rx only. ©2008 B. Braun Medical Inc. Bethlehem PA. All rights reserved. IS47-A 9/08 JH