

# Filter Straw®

5 micron particulate ampule filter

Safety in a Snap



Enhance patient safety by reducing potential for glass particle injection when using single-dose glass ampules<sup>1</sup>

**B | BRAUN**




IN THE PHARMACY

ON THE FLOOR


IN HOMECARE SETTINGS

**FILTER STRAW®**  
particulate matter filter

- 5 micron particulate filter reduces injection of contaminants<sup>2</sup>
- Reduction of contaminants lessens potential side effects, such as phlebitis, caused by infused particulates<sup>3</sup>
- Flexible straw can be angled for easy use
- Needle-free and latex-free



Shown Actual Size



**IN YOUR PROCEDURE**

CVC and Regional Anesthesia kits include the FILTER STRAW particulate matter filter as a component to enhance safety.

Also available for use in custom kits and private labeling.

**Ordering Information**

REF 415021

FS5005 - 5 micron filter with 1¾ inch flexible straw (shown)

REF 415020

FS5000 - 5 micron filter with 4 inch flexible straw

For additional information or to request a catalog, contact your B. Braun representative, or call **1-800-BBRAUN2 (1-800-227-2862)**.



B. Braun Medical Inc.  
824 Twelfth Avenue  
Bethlehem, PA 18018  
www.bbraunusa.com

1. Furgang F: Glass particles in ampules. Anesthesiology 1973; 41:525-526.  
2. Sabon R, Cheng E, Stommel K, Hennen C: Glass Particle Contamination: Influence of Aspiration Methods and Ampule Types. Anesthesiology. 1989; 70:859-862.  
3. Lusardi P: Research Corner: Myth vs Reality: Do We Need to Filter Medication Ampules? www.aacn.org. 2003.