

GRAVITY FLOW RATE DRIP CHART

Flow Rate (mL/hr)	10 drops=1mL (drops/min)	15 drops=1mL (drops/min)	20 drops=1mL (drops/min)	60 drops=1mL (drops/min)
10	2	2	3	10
25	4	6	8	25
50	8	12	17	50
75	12	19	25	75
100	17	25	33	100
125	21	31	42	125
150	25	37	50	150
200	33	50	67	200
250	42	62	83	250
500	83	125	167	500
1000	167	250	333	1000

- Confirm tubing set drip rate on set package, i.e. 10, 15, 20 or 60 drops/mL
- Recommended that all gravity infusion bags be time taped for additional flow confirmation
- Alterations of bag height distance to patient will affect flow rate

