



B. BRAUN PAIN CONTROL PRODUCTS

Product Catalog

Effective Q2 2019

CONTENTS

SPINAL NEEDLES

2

- PENCAN® Pencil Point Spinal Needles
- SPINOCAN® Spinal Needles with Quincke Bevel
- Introducer Needle



SPINAL TRAYS

3-4

- PENCAN Spinal Trays
- PENCAN Safety Engineered SPINAL Tray
- SPINOCAN Spinal Trays
- Spinal Support Tray



EPIDURAL NEEDLES AND CATHETERS

5


- PERIFIX® Epidural Needles
- EPICAN™ Caudal Needles
- PERIFIX® FX Springwound Epidural Catheters
- PERIFIX® Epidural Catheters



SYRINGES, EPIDURAL ACCESSORIES AND SETS

6

- Glass Loss-of-Resistance Syringes
- PERIFIX™ Plastic Loss-of-Resistance Syringes
- Connectors for Epidural Catheters
- Filter for Epidural Catheters and PinPad
- PERIFIX® FX Springwound Epidural Sets with 17 Ga. Needle
- PERIFIX® Epidural Sets with 17 Ga. Needle
- PERIFIX® Epidural Sets with 18 Ga. Needle
- PERIFIX® ONE Epidural Set
- PERIFIX®ONE Pediatric Epidural Sets



EPIDURAL TRAYS

7-10

- PERIFIX® Single Dose Epidural Trays
- PERIFIX® Single Dose Safety-Engineered Epidural Tray
- PERIFIX® Continuous Epidural Trays with 17 Ga. Epidural Needle
- PERIFIX® Continuous Safety-Engineered Epidural Tray with 17 Ga. Epidural Needle
- PERIFIX® Continuous Epidural Trays with 18 Ga. Epidural Needle
- PERIFIX® Continuous Safety-Engineered Epidural Tray with 18 Ga. Epidural Needle
- PERIFIX® ONE Pediatric Epidural Trays



COMBINED SPINAL/EPIDURAL SETS AND TRAYS

11-12

- ESPOCAN® Combined Spinal/Epidural Sets
- ESPOCAN® and Non-ESPOCAN Combined Spinal/Epidural Trays



SINGLE SHOT NERVE BLOCK NEEDLES, NERVE BLOCK MONITORS AND ACCESSORIES

13

- STIMUPLEX® A Insulated Needles with 30° Bevel
- STIMUPLEX Ultra 360® Insulated Echogenic Needles with 30° Bevel
- STIMUPLEX D Insulated Needles with 15° and 30° Bevel
- ULTRAPLEX 360® Non-insulated Echogenic Needles with 30° Bevel
- STIMUPLEX Nerve Stimulator Products
- B-Smart™ Injection Pressure Monitor



CONTINUOUS NERVE BLOCK SETS

14

- CONTIPLEX® C Catheter-Over-Needle Continuous Nerve Block Sets (Insulated Needle)
- CONTIPLEX Echo Continuous Nerve Block Sets (Insulated Echogenic Contiplex Ultra 360 Tuohy Needle) with Echogenic Non-Stimulating Catheter
- CONTIPLEX Ultra 360 Continuous Nerve Block Sets (Insulated Echogenic Contiplex Ultra 360 Tuohy Needle) with Non-Stimulating Catheter
- CONTIPLEX FX Continuous Nerve Block Sets (Non-Insulated Tuohy Needle) with Non-Stimulating Springwound Catheter
- CONTIPLEX STIM Continuous Nerve Block Sets (Insulated Tuohy Needle) with Stimulating Catheter
- CONTIPLEX Tuohy Continuous Nerve Block Sets (Insulated Tuohy Needle) with Non-Stimulating Catheter
- CONTIPLEX Straight Continuous Nerve Block Sets (Insulated Needle) with Non-Stimulating Catheter
- PERICAN™ Ultra Sets



NERVE BLOCK TRAYS

15-16

- Nerve Block Support Trays
- STIMUPLEX® Ultra 360 Single Shot Nerve Block Trays
- CONTIPLEX Continuous Nerve Block Trays
- CONTIPLEX Ultra 360 Continuous Nerve Block Trays



XPERIUS® ULTRASOUND SYSTEM

17



Ordering Information:

For more information on B. Braun Pain Control Products, visit www.BBraunUSA.com or call 800-227-2862

Clinical/Technical Support 1-800-854-6851

SPINAL NEEDLES

NEEDLES



Pencil Point Needles

PENCAN® Pencil Point Spinal Needles

	REF #	Product Code	Units Per Cs
PENCAN 27 Ga. x 3½ in. (90 mm)	333873	P2735	25
PENCAN 27 Ga. x 5 in. (130 mm)	333876	P27127	25
PENCAN 25 Ga. x 1 in. (25 mm)	333878	P2510	25
PENCAN 25 Ga. x 2 in. (50 mm)	333880	P2520	25
PENCAN 25 Ga. x 3½ in. (90 mm)	333853	P2535	25
PENCAN 25 Ga. x 5 in. (130 mm)	333875	P25127	25
PENCAN 25 Ga. x 6.1 in. (156 mm)	333877	P25155	25
PENCAN 24 Ga. x 4 in. (100 mm) with 20 Ga. x 1⅜ in. (35 mm) Introducer	333884	P2440	25
PENCAN 22 Ga. x 3½ in. (90 mm)	333866	P2235	25

Quincke Bevel Needles

SPINOCAN® Spinal Needles with Quincke Bevel



SPINOCAN 27 Ga. x 3½ in. (90 mm)	333290	S2735	50
SPINOCAN 26 Ga. x 3½ in. (90 mm)	333300	S2635	50
SPINOCAN 25 Ga. x 3½ in. (90 mm)	333313	S2535	50
SPINOCAN 25 Ga. x 4¾ in. (120 mm)	333308	S25475	50
SPINOCAN 22 Ga. x 1½ in. (40 mm)	333325	S2215	50
SPINOCAN 22 Ga. x 3½ in. (90 mm)	333320	S2235	50
SPINOCAN 22 Ga. x 4¾ in. (120 mm)	333315	S22475	50
SPINOCAN 22 Ga. x 5 in. (130 mm)	333355	S22127	50
SPINOCAN 22 Ga. x 7 in. (180 mm)	333357	S22178	50
SPINOCAN 20 Ga. x 3½ in. (90 mm)	333335	S2035	50
SPINOCAN 18 Ga. x 3½ in. (90 mm)	333350	S1835	50



Introducer Needle

20 Ga. x 1⅜ in. (35 mm) - for use with 25, 26, and 27 Ga. Spinal Needles	333370	SI2000	50
--	--------	--------	----



Syringes

5 mL Plastic - Luer Slip	332145	PS5S	50
5 mL Glass - Luer Slip	332147	GS5MS	10

SPINAL TRAYS

TRAYS

PENCAN® SPINAL Trays

PENCAN Safety Engineered SPINAL Tray

Units Per Cs.	10										
Product Code	P24SK	P24BK	P24BKG	P24BKDP	P25K	P25SK	P25BK	P25BKG	P25BKDP	P27BK	P25BKAY
REF #	333881	333868	560606	333867	333858	333855	333851	560603	333862	333871	333865

Item	P24SK	P24BK	P24BKG	P24BKDP	P25K	P25SK	P25BK	P25BKG	P25BKDP	P27BK	P25BKAY
PROCEDURAL NEEDLES - SPINAL											
PENCAN 27 Ga. x 3½ in. (90 mm)											•
PENCAN 25 Ga. x 3½ in. (90 mm)						•	•	•	•		•
PENCAN 24 Ga. x 4 in. (100 mm)	•	•	•	•							
INTRODUCER NEEDLES											
20 Ga. x 1¼ in. TW Introducer						•	•	•	•	•	•
20 Ga. x 1⅜ in. Introducer	•	•	•	•							
HYPODERMIC NEEDLES											
25 Ga. x 1 in. SafetyGlide®											•
25 Ga. x 1½ in.	•	•	•	•	•	•	•	•	•	•	•
22 Ga. x 1½ in. SafetyGlide											•
18 Ga. x 1½ in. SafetyGlide											•
3mL Plastic Luer Lock Syringe with Pre-Attached 18 Ga. x 1½ in. Needle	•	•	•	•	•	•	•	•	•	•	•
PROCEDURAL SYRINGES											
5 mL Plastic - Luer Slip	•	•		•	•	•	•	•	•	•	•
5 mL Glass - Luer Slip			•					•			
SYRINGE											
3 mL Plastic Luer Lock Syringe with Pre-Attached Filter Straw®											•
PROCEDURAL DRUGS											
Bupivacaine 0.75% with Dextrose 8.25%, 2 mL			•	•	•			•	•	•	•
Lidocaine HCl (1%), 5 mL (skin wheal)	•	•	•	•		•	•	•	•	•	•
DRAPE											
Drape, Clear plastic, 4 in. Fenestrated, 30 in. x 24 in.	•	•	•	•	•	•	•	•	•	•	•
ACCESSORIES											
Packet Povidone-Iodine Solution	•	•	•	•	•	•	•	•	•	•	•
Prep Sponges (3 per tray)	•	•	•	•	•	•	•	•	•	•	•
DuraPrep® Solution				•					•		
Needle Block Foam with Holder	•	•	•	•	•	•	•	•	•	•	•
NeedleVISE®											•
Filter Straw® Particulate Matter Filter	•	•	•	•	•	•	•	•	•	•	•
Paper Towel, 3-ply	•	•	•	•	•	•	•	•	•	•	•
4 in. x 4 in. Gauze Sponges (4 per tray)	•	•	•	•	•	•	•	•	•	•	•
Drug Identification Label Set	•	•	•	•	•	•	•	•	•	•	•

SAFETY-ENGINEERED TRAY

DURAPREP is a registered trademark of 3M
 SafetyGlide is a registered trademark of Becton Dickinson and Company
 NeedleVISE is a registered trademark of Atrion Medical Products, Inc.

SPINAL TRAYS

TRAYS

SPINOCAN® Spinal Trays Spinal Support Tray

Units Per Cs.	10					
Product Code	S22K	S22SK	S22BK	S25SK	S25BK	S26BK
REF #	333235	333234	333231	333214	333211	333201
						333400

	333235	333234	333231	333214	333211	333201	333400
PROCEDURAL NEEDLES – SPINAL							
SPINOCAN 26 Ga. x 3½ in. (90 mm)							●
SPINOCAN 25 Ga. x 3½ in. (90 mm)				●	●		
SPINOCAN 22 Ga. x 3½ in. (90 mm)	●	●	●				▶
INTRODUCER NEEDLES							
20 Ga. x 1¼ in. TW Introducer							●
18 Ga. x 1¼ in. Introducer				●	●		
HYPODERMIC NEEDLES							
25 Ga. x 5/8 in.			●		●	●	
25 Ga. x 1½ in.	●	●	●	●	●	●	●
22 Ga. x 1½ in.	●	●	●	●	●	●	●
3mL Plastic Luer Lock Syringe with Pre-Attached 18 Ga. x 1½ in. Needle	●	●	●	●	●	●	●
SYRINGE							
5 mL Plastic – Luer Slip	●	●	●	●	●	●	●
PROCEDURAL DRUGS							
Bupivacaine 0.75% with Dextrose 8.25%, 2 mL			●		●	●	▶
Lidocaine HCl (1%), 5 mL (skin wheal)		●	●	●	●	●	●
DRAPE							
Drape, Clear plastic, 4 in. Fenestrated, 30 in. x 24 in.	●	●	●	●	●	●	●
ACCESSORIES							
Packet Povidone-Iodine Solution	●	●	●	●	●	●	●
Prep Sponges (3 per tray)	●	●	●	●	●	●	●
Needle Block Foam with Holder	●	●	●	●	●	●	●
Filter Straw® Particulate Matter Filter	●	●	●	●	●	●	●
Paper Towel, 3-ply	●	●	●	●	●	●	●
4 in. x 4 in. Gauze Sponges (4 per tray)	●	●	●	●	●	●	●
Drug Identification Label Set	●	●	●	●	●	●	●

EPIDURAL NEEDLES AND CATHETERS

NEEDLES



PERIFIX® Epidural Needles	REF #	Product Code	Units Per Cs
22 Ga. x 2 in. (50 mm) Tuohy Needle – non-winged and clear hub	332173	E2220TN	25
22 Ga. x 3.15 in. (80 mm) Tuohy Needle – with wings and clear hub	332171	E2230T	25
22 Ga. x 4 in. (100 mm) Tuohy Needle – with wings and clear hub	332194	E2240T	25
20 Ga. x 2 in. (50 mm) Tuohy Needle – with wings and clear hub	332178	E2020T	25
20 Ga. x 3 1/2 in. (90 mm) Tuohy Needle – non-winged and clear hub	332167	E2035TN	25
20 Ga. x 3 1/2 in. (90 mm) Tuohy Needle – with wings and clear hub	332168	E2035T	25
20 Ga. x 4 in. (100 mm) Tuohy Needle – with wings and clear hub	332189	E2040T	25
20 Ga. x 6 in. (150 mm) Tuohy Needle – with wings and clear hub	332193	E2060T	25
18 Ga. x 2 in. (50 mm) Tuohy Needle – with wings and clear hub	332172	E1820TN	25
18 Ga. x 3 1/2 in. (90 mm) Tuohy Needle – with wings and clear hub	332166	E1835T	25
18 Ga. x 3 1/2 in. (90 mm) Hustead Needle – with wings and clear hub	332165	E1835H	25
18 Ga. x 4 3/4 in. (120 mm) Tuohy Needle – with wings and clear hub	332186	E18475T	25
18 Ga. x 6 in. (150 mm) Tuohy Needle – with wings and clear hub	332160	E1860T	25
17 Ga. x 3 1/2 in. (90 mm) Tuohy Needle – with wings and clear hub	332180	E1735T	25
17 Ga. x 3 1/2 in. (90 mm) Dull Tuohy Needle – with wings and clear hub	332176	E1735DT	25
17 Ga. x 4 in. (100 mm) Dull Tuohy Needle – non-winged and metal hub	332195	E1740DTN	25
17 Ga. x 4 1/2 in. (115 mm) Tuohy Needle – with wings and clear hub	332181	E1745T	25
17 Ga. x 6 in. (150 mm) Tuohy Needle – with wings and clear hub	332192	E1760T	25

EPICAN™ Caudal Needles



EPICAN 22 Ga. x 1 3/8 in. (35 mm) Caudal Needle – non-winged and clear hub	332177	E22138CN	25
EPICAN 20 Ga. x 2 in. (50 mm) Caudal Needle – non-winged and clear hub	332175	E2020CN	25

CATHETERS

PERIFIX® FX Springwound Epidural Catheters




19 Ga. x 36 in. (91 cm) Radiopaque Polyamide Springwound PERIFIX FX Catheter (MR unsafe) CLOSED TIP with three lateral side ports, with clamp style catheter connector and thread assist guide	333512	EC19CF	25
--	--------	--------	----



19 Ga. x 36 in. (91 cm) Radiopaque Polyamide Springwound PERIFIX FX Catheter (MR unsafe) OPEN TIP with clamp style catheter connector and threading assist guide	333514	EC19OF	25
--	--------	--------	----





PERIFIX® Epidural Catheters



CLOSED TIP with six lateral side ports, with clamp style catheter connector and threading assist guide, 20 Ga. x 40 in. (101 cm) Polyamide/Polyurethane PERIFIX ONE Catheter –	333541	EC20CNG	25
CLOSED TIP with six lateral side ports, with clamp style catheter connector and threading assist guide, 20 Ga. x 40 in. (101 cm) Polyamide PERIFIX SOFT TIP Catheter – CLOSED TIP with three lateral side ports, with clamp style catheter connector and threading assist guide	333521	EC20CNG	25
20 Ga. x 40 in. (101 cm) Polyamide STYLETED PERIFIX Catheter – CLOSED TIP with three lateral side ports, with clamp style catheter connector and threading assist guide	333532	EC20CS	25
20 Ga. x 40 in. (101 cm) Polyamide PERIFIX Catheter – CLOSED TIP with three lateral side ports, with clamp style catheter connector and threading assist guide	333540	EC20C	25
20 Ga. x 40 in. (101 cm) Polyamide PERIFIX Catheter – OPEN TIP with clamp style catheter connector and threading assist guide	333520	EC200	25
19 Ga. x 40 in. (101 cm) Polyamide PERIFIX Catheter – CLOSED TIP with three lateral side ports, with clamp style catheter connector and threading assist guide	333511	EC19C	25

SYRINGES, EPIDURAL ACCESSORIES AND SETS

SYRINGES

	REF #	Product Code	Units Per Cs
	Glass Syringes		
	5 mL Glass - Luer Lock, Metal Tip	332155	GLR5ML 10
	5 mL Glass - Luer Slip, Metal Tip	332156	GLR5MS 10
	10 mL Glass - Luer Lock, Metal Tip	332157	GLR10ML 10
	10 mL Glass - Luer Slip, Metal Tip	332158	GLR10MS 10
	PERIFIX™ Plastic Loss-of-Resistance Syringes		
	8 mL Plastic - Luer Lock	332152	PLR10L 50
	8 mL Plastic - Luer Slip	332150	PLR10S 50
	Connectors for Epidural Catheters		
	PERIFIX™ Clamp Style Catheter Connector – for use with PERIFIX® 20–24 Ga. epidural catheters	332283	PCC2000 50
	PERIFIX™ Clamp Style Catheter Connector – for use with PERIFIX® 19 Ga. epidural catheters	332285	PCC1900 50
	Filter for Epidural Catheters and PinPad		
	PERIFIX® 0.2 µm Flat Epidural Filter – for use with all PERIFIX® catheters and PERIFIX™ PinPad	415000	FEF100 25
	PERIFIX PinPad – securement device for PERIFIX 0.2 µm Flat Epidural Filter	415004	FPP100 25

PERIFIX® EPIDURAL SETS

	REF #	CE17TF	CE17TFC	CE17T	CE18T	CE18H	CEP18T0	CEP20T0	
		332234	332238	332202	332200	332210	332291	332293	
		Units Per Cs. 12							
		Product Code							
EPIDURAL NEEDLES - WINGED									
20 Ga. x 2 in. (50 mm) Tuohy Needle								●	
18 Ga. x 2 in. (50 mm) Tuohy Needle							●		
18 Ga. x 3½ in. (90 mm) Tuohy Needle					●				
18 Ga. x 3½ in. (90 mm) Hustead Needle						●			
17 Ga. x 3½ in. (90 mm) Tuohy Needle		●	●	●					
EPIDURAL CATHETERS									
PERIFIX ONE 24 Ga. x 28 in. (72 cm) Polyamide/Polyurethane Pediatric closed tip catheter								●	
PERIFIX ONE 20 Ga. x 40 in. (101 cm) Polyamide/Polyurethane Pediatric closed tip catheter							●		
PERIFIX 20 Ga. x 40 in. (101 cm) closed tip catheter				●	●	●			
PERIFIX FX 19 Ga. x 36 in. (91 cm) Springwound closed tip catheter (MR unsafe)			●						
PERIFIX FX 19 Ga. x 36 in. (91 cm) Springwound open tip catheter (MR unsafe)		●							
ACCESSORIES									
Clamp Style Catheter Connector		●	●	●	●	●	●	●	
Threading Assist Guide		●	●	●	●	●	●	●	

PEDIATRIC SETS

SINGLE DOSE EPIDURAL TRAYS

PERIFIX® Single Dose Epidural Trays
PERIFIX® Single Dose Safety-Engineered Epidural Tray

Units Per Cs. 10

Product Code

REF #

Product Code	REF #
SD17TD	332248
SD18TB	332250
SD18TD	332251
SD18HD	332249
SD18TK	332256
SD18DPS	332247
SD20DPL	332258
SD20DPS	332259
SD20TD	332254
SD22TD	332255
SD18TDY	332252

PROCEDURAL NEEDLES – EPIDURAL

PERIFIX 22 Ga. x 3 in. (80 mm) Tuohy, Metal Stylet, Winged
 PERIFIX 20 Ga. x 3½ in. (90 mm) Tuohy, Metal Stylet, Winged
 PERIFIX 18 Ga. x 3½ in. (90 mm) Tuohy, Winged
 PERIFIX 18 Ga. x 3½ in. (90 mm) Hustead, Winged
 PERIFIX 17 Ga. x 3½ in. (90 mm) Tuohy, Winged

HYPODERMIC NEEDLES

27 Ga. x 1¼ in.
 25 Ga. x 1½ in.
 25 Ga. x 5/8 in.
 25 Ga. x 1 in. SafetyGlide®
 22 Ga. x 1½ in.
 22 Ga. x 1½ in. SafetyGlide
 18 Ga. x 1½ in. SafetyGlide
 3mL Plastic Luer Lock Syringe with Pre-Attached 18 Ga. x 1½ in. Needle

LOSS-OF-RESISTANCE SYRINGES

5 mL Glass – Luer Lock, Metal Tip
 PERIFIX™ 8 mL Plastic – Luer Slip
 PERIFIX™ 8 mL Plastic – Luer Lock

SYRINGES

1 mL Plastic – Luer Slip
 3 mL Plastic – Luer Lock
 10 mL Plastic – Luer Lock
 20 mL Plastic – Luer Lock

PROCEDURAL DRUGS

Lidocaine HCl (1%), 5 mL (skin wheal)
 Lidocaine HCl (1.5%) with Epi, 5 mL (test dose)

DILUENT/LOSS-OF-RESISTANCE

Sodium Chloride (0.9%), 10 mL

DRAPES

Drape, Clear Plastic, 4 in. Fenestrated, 30 in. x 24 in.
 Drape, Blue, 4 in. Fenestrated, 30 in. x 24 in.

ACCESSORIES

Packet Povidone – Iodine Solution
 Prep Sponges (3 per tray)
 Needle Block Foam with Holder
 NeedleVISE®
 MP1000 MICRO PIN® Dispensing Pin
 Filter Straw® Particulate Matter Filter
 4 in. x 4 in. Gauze Sponges (4 per tray)
 Paper Towel, 3 Ply
 Drug Identification Label Set

SAFETY-ENGINEERED TRAY

SafetyGlide is a registered trademark of Becton Dickinson and Company
 NeedleVISE is a registered trademark of Atrion Medical Products, Inc.

CONTINUOUS EPIDURAL TRAYS (17 Ga.)

PERIFIX® Continuous Epidural Trays with 17 Ga. Epidural Needle

PERIFIX Continuous Safety-Engineered Epidural Tray with 17 Ga. Epidural Needle

Units Per Cs. 10

Product Code

REF #

Product Code	CE17TKF	CE17TKFPS	CE17TBFC	CE17TKFCPS	CE17TKFC	CE17TKFSA	CE17TKFCS	CE17TKFS	CE17TKFSDT	CE17TKFSDTD	CE17TKFMHDT	CE17TB	CE17TOD	CE17TK	CE17TKPS	CE17TKCDLS	CE17TKST	CE17TB19C	CE17TK19CDP	CE17TKFCY
332086																				
332087																				
332077																				
332078																				
332079																				
332080																				
332097																				
332098																				
332093																				
332095																				
332096																				
332224																				
332235																				
332221																				
332233																				
332212																				
332082																				
332242																				
332262																				
332075																				

PROCEDURAL NEEDLES – EPIDURAL

PERIFIX® 17 Ga. x 3½ in. (90 mm) Tuohy, Winged																				
17 Ga. x 3½ in. (90 mm) Dull Tuohy with Plastic Hub, Winged																				
PERIFIX 17 Ga. x 4 in. (100 mm) Tuohy, Winged																				
17 Ga. x 4 in. (100 mm) Dull Tuohy with Metal Hub and Stylet, Non-Winged																				

HYPODERMIC NEEDLES

27 Ga. x 1¼ in.																				
25 Ga. x 1½ in.																				
25 Ga. x 5/8 in.																				
25 Ga. x 1 in. SafetyGlide®																				
22 Ga. x 1½ in.																				
22 Ga. x 1½ in. SafetyGlide																				
18 Ga. x 1½ in.																				
18 Ga. x 1½ in. SafetyGlide																				
3mL Plastic Luer Lock Syringe with Pre-Attached 25 Ga. x 5/8 in. Needle																				
3mL Plastic Luer Lock Syringe with Pre-Attached 25 Ga. x 1½ in. Needle																				
3mL Plastic Luer Lock Syringe with Pre-Attached 18 Ga. x 1½ in. Needle																				

LOSS-OF-RESISTANCE SYRINGES

5 mL Glass – Luer Slip, Metal Tip																				
5 mL Glass – Luer Lock, Metal Tip																				
PERIFIX™ 8 mL Plastic – Luer Slip																				

SYRINGES

1 mL Plastic – Luer Slip																				
3 mL Plastic – Luer Lock																				
10 mL Plastic – Luer Slip																				
20 mL Plastic – Luer Lock																				

EPIDURAL CATHETERS

PERIFIX® 20 Ga. Closed Tip																				
PERIFIX 20 Ga. Open Tip																				
PERIFIX 20 Ga. Soft Tip, Closed Tip																				
PERIFIX 19 Ga. Closed Tip																				
PERIFIX FX 19 Ga. Springwound, Open Tip (MR unsafe)																				
PERIFIX FX 19 Ga. Springwound, Open Tip (MR unsafe) with dispenser																				
PERIFIX FX 19 Ga. Springwound, Closed Tip (MR unsafe)																				

PROCEDURAL DRUGS

Lidocaine HCl (1%), 5 mL (skin wheal)																				
Lidocaine HCl (1.5%), with Epi, 5 mL (test dose)																				

DILUENT/LOSS-OF-RESISTANCE

Sodium Chloride (0.9%), 10 mL																				
-------------------------------	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--

DRAPES

Drape, Blue, Non-Fenestrated																				
Drape, Blue, 4 in. Fenestrated, 30 in. x 24 in.																				
Drape, Clear plastic, 4 in. Fenestrated, 30 in. x 24 in.																				
Drape, Clear plastic, 4 in. Fenestrated, 30 in. x 24 in. with peel patch																				

ACCESSORIES

Packet Povidone-Iodine Solution																				
Prep Sponges (3 per tray)																				
DuraPrep® Solution																				
PERIFIX® 0.2 µm Flat Epidural Filter																				
PERIFIX™ PinPad (securement device for PERIFIX 0.2 µm Flat Epidural Filter)																				
Catheter Support Pad																				
Tegaderm® 4 in. x 4¾ in. Dressing																				
Needle Block Foam with Holder																				
NeedleVISE®																				
Female Luer Lock Connector (converts Luer Lock to Luer Slip)																				
MP1000 MICRO PIN® Dispensing Pin																				
Filter Straw® Particulate Matter Filter																				
Paper Towel, 3-ply																				
4 in. x 4 in. Gauze Sponges (4 per tray)																				
Drug Identification Label Set																				
Threading Assist Guide																				
Clamp Style Catheter Connector																				

SafetyGlide is a registered trademark of Becton Dickinson and Company. DuraPrep and Tegaderm are registered trademarks of 3M. NeedleVISE is a registered trademark of Atrion Medical Products, Inc.

CONTINUOUS EPIDURAL TRAYS (18 Ga.)

PERIFIX® Continuous Epidural Trays with 18 Ga. Epidural Needle

Units Per Cs.

10

PERIFIX Continuous Safety-Engineered Epidural Tray with 18 Ga. Epidural Needle

Product Code

REF #

	CE18TB	CE18HB	CE18HK	CCE18HK10L	CE18HKST	CE18HKPS	CE18TBS	CE18TOD	CE18TK	CE18TKPS	CE18TKDLS	CE18TKD10L	CE18TKST	CE18TKDP	CE18TKO	CE18TKY
PROCEDURAL NEEDLES – EPIDURAL																
PERIFIX® 18 Ga. x 3½ in. (90 mm) Tuohy, Winged	•							•	•	•	•	•	•	•	•	•
PERIFIX 18 Ga. x 3½ in. (90 mm) Hustead, Winged		•	•	•	•	•										
HYPODERMIC NEEDLES																
27 Ga. x 1¼ in.					•	•				•	•	•	•			
25 Ga. x 1½ in.	•	•	•	•	•	•		•	•	•	•	•	•	•	•	•
25 Ga. x 5⁄8 in.						•	•			•						
25 Ga. x 1 in. SafetyGlide®																•
22 Ga. x 1½ in.					•	•	•			•			•			
22 Ga. x 1½ in. SafetyGlide																•
18 Ga. x 1½ in. SafetyGlide																•
3 mL Plastic Luer Lock Syringe with Pre-Attached 18 Ga. x 1½ in. Needle	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
LOSS-OF-RESISTANCE SYRINGES																
5 mL Glass – Luer Slip, Metal Tip								•			•					
5 mL Glass – Luer Lock, Metal Tip	•	•	•		•				•	•			•	•	•	•
10 mL Glass – Luer Lock, Metal Tip				•								•				
PERIFIX™ 8 mL Plastic – Luer Slip						•				•						
SYRINGES																
3 mL Plastic – Luer Lock																
20 mL Plastic – Luer Lock	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
EPIDURAL CATHETERS																
PERIFIX® 20 Ga. Closed Tip	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
PERIFIX 20 Ga. Open Tip								•								
PERIFIX 20 Ga. Soft Tip, Closed Tip					•								•			
PERIFIX ONE 20 Ga. Closed Tip (with six lateral sideports)															•	
PROCEDURAL DRUGS																
Lidocaine HCl (1%), 5 mL (skin wheal)			•	•	•	•	•	•	•	•	•	•	•	•	•	•
Lidocaine HCl (1.5%), with Epi, 5 mL (test) dose			•	•	•	•	•	•	•	•	•	•	•	•	•	•
DILUENT/LOSS-OF-RESISTANCE																
Sodium Chloride (0.9%), 10 mL			•	•	•	•	•	•	•	•	•	•	•	•	•	•
DRAPES																
Drape, Blue, Non-Fenestrated			•						•		•					
Drape, Blue, 4 in. Fenestrated, 30 in. x 24 in.				•	•				•	•	•	•	•	•	•	•
Drape, Clear plastic, 4 in. Fenestrated, 30 in. x 24 in.	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
ACCESSORIES																
Packet Povidone-Iodine Solution	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Prep Sponges (3 per tray)	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
DuraPrep® Solution														•		
PERIFIX® 0.2 µm Flat Epidural Filter			•	•	•	•			•	•	•	•	•	•	•	•
PERIFIX™ PinPad (securement device for PERIFIX 0.2 µm Flat Epidural Filter)																•
Catheter Support Pad	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Tegaderm® 4 in. x 4¾ in. Dressing										•						
Needle Block Foam with Holder	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
NeedleVISE®										•						
MP1000 MICRO PIN® Dispensing Pin																•
Filter Straw® Particulate Matter Filter	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Paper Towel, 3-ply	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
4 in. x 4 in. Gauze Sponges (4 per tray)	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Drug Identification Label Set	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Threading Assist Guide	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
Clamp Style Catheter Connector	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•

SAFETY-ENGINEERED TRAY

SafetyGlide is a registered trademark of Becton Dickinson and Company. DuraPrep and Tegaderm are registered trademarks of 3M. NeedleVISE is a registered trademark of Atrion Medical Products, Inc.

PEDIATRIC EPIDURAL TRAYS

PERIFIX® ONE Pediatric Epidural Trays	Units Per Cs.	10	
	Product Code	CEP18TKO	CEP20TKO
	REF #	332290	332292

PROCEDURAL NEEDLES – EPIDURAL		
PERIFIX® 20 Ga. x 2 in. (50 mm) Tuohy, Winged		●
PERIFIX 18 Ga. x 2 in. (50 mm) Tuohy, Winged	●	
HYPODERMIC NEEDLES		
25 Ga. x 1½ in.	●	●
3mL Plastic Luer Lock Syringe with Pre-Attached 18 Ga. x 1½ in. Needle	●	●
LOSS-OF-RESISTANCE SYRINGES		
5 mL Glass – Luer Slip, Metal Tip	●	●
SYRINGE		
20 mL Plastic – Luer Lock	●	●
EPIDURAL CATHETERS		
PERIFIX® ONE 24 Ga. Closed Tip (with six lateral sideports)		●
PERIFIX ONE 20 Ga. Closed Tip (with six lateral sideports)	●	
PROCEDURAL DRUG		
Lidocaine HCl (1.5%), with Epi, 5 mL (test dose)	●	●
DILUENT/LOSS-OF-RESISTANCE		
Sodium Chloride (0.9%), 10 mL	●	●
DRAPE		
Drape, Clear plastic, 4 in. Fenestrated, 15 in. x 15 in.	●	●
ACCESSORIES		
Packet Povidone-Iodine Solution	●	●
Prep Sponges (3 per tray)	●	●
PERIFIX® 0.2 µm Flat Epidural Filter	●	●
Catheter Support Pad	●	●
Benzoin Swab Stick	●	●
Needle Block Foam with Holder	●	●
Tegaderm® 4 in. x 4¾ in. Dressing	●	●
Filter Straw® Particulate Matter Filter	●	●
Paper Towel, 3-ply	●	●
4 in. x 4 in. Gauze Sponges (4 per tray)	●	●
Drug Identification Label Set	●	●
Threading Assist Guide	●	●
Clamp Style Catheter Connector	●	●

COMBINED SPINAL/EPIDURAL SETS AND TRAYS

ESPOCAN® Combined Spinal/Epidural Anesthesia Sets

	Units Per Cs.	12						
	Product Code	ES1725	ES1827DS	ES1827	ES1725L	ES1727L	ES1827L	ES1827Q
REF #		333161	333174	333160	333150	333151	333152	333180
EPIDURAL NEEDLES								
ESPOCAN 18 Ga. x 3½ in. (90 mm) Tuohy Needle with backeye lumen			●	●				●
ESPOCAN 18 Ga. x 4½ in. (110 mm) Tuohy Needle with backeye lumen							●	
ESPOCAN 17 Ga. x 3½ in. (90 mm) Tuohy Needle with backeye lumen		●						▶
ESPOCAN 17 Ga. x 4½ in. (110 mm) Tuohy Needle with backeye lumen					●	●		
SPINAL NEEDLES								
SPINOCAN® 27 Ga. x 5 in. (130 mm) Needle with centering sleeve								●
PENCAN® 27 Ga. x 5 in. (130 mm) Pencil Point Needle with centering sleeve				●				
PENCAN 27 Ga. x 5½ in. (140 mm) Pencil Point Needle with centering sleeve			●					
PENCAN 27 Ga. x 6 in. (150 mm) Pencil Point Needle with centering sleeve						●	●	
PENCAN 25 Ga. x 5 in. (130 mm) Pencil Point Needle with centering sleeve		●						▶
PENCAN 25 Ga. x 6 in. (150 mm) Pencil Point Needle with centering sleeve					●			
SYRINGES								
5 mL Clear Plastic Luer Slip		●	●	●	●	●	●	●
ACCESSORIES								
ESPOCAN Docking System			●					

COMBINED SPINAL/EPIDURAL SETS AND TRAYS

ESPOCAN® and Non-ESPOCAN Combined Spinal/Epidural Trays

Units Per Cs.	10						
Product Code	ES1725K	ES1725KFX	NES1727KFX	ES1827K	ES1827QK	ES1827KDS	ES1827KO
REF #	333165	333196	333197	333192	333182	333194	333193

	333165	333196	333197	333192	333182	333194	333193
PROCEDURAL NEEDLES – EPIDURAL							
ESPOCAN 18 Ga. x 3½ in. (90 mm) Tuohy, Winged				•	•	•	•
ESPOCAN 17 Ga. x 3½ in. (90 mm) Tuohy, Winged	•	•					
PERIFIX® 17 Ga. x 3½ in. (90 mm) Tuohy, Winged			•				
PROCEDURAL NEEDLES – SPINAL							
SPINOCAN® 27 Ga. x 5 in. (130 mm) with Centering Sleeve					•		
PENCAN® 27 Ga. x 5 in. (130 mm) with Centering Sleeve				•			•
PENCAN 27 Ga. x 5 in. (130 mm) without Centering Sleeve			•				
PENCAN 27 Ga. x 5½ in. (140 mm) with Centering Sleeve						•	
PENCAN 25 Ga. x 5 in. (130 mm) with Centering Sleeve	•	•					
HYPODERMIC NEEDLES							
27 Ga. x 1¼ in.		•					
25 Ga. x 5⁄8 in.					•		
25 Ga. x 1½ in.	•	•	•	•		•	•
22 Ga. x 1½ in.		•	•		•		
3mL Plastic Luer Lock Syringe with Pre-Attached 18 Ga. x 1½ in. Needle	•	•	•	•	•	•	•
LOSS-OF-RESISTANCE SYRINGES							
5 mL Glass – Luer Slip, Metal Tip		•	•				
5 mL Glass – Luer Lock, Metal Tip	•			•	•	•	•
SYRINGES							
1 mL Plastic – Luer Slip	•	•		•	•	•	•
3 mL Plastic – Luer Slip			•				
5 mL Clear Plastic – Luer Slip	•	•	•	•	•	•	•
20 mL Plastic – Luer Lock	•	•	•	•	•	•	•
EPIDURAL CATHETERS							
PERIFIX® ONE 20 Ga. Closed Tip							•
PERIFIX 20 Ga. Soft Tip, Closed Tip						•	
PERIFIX 20 Ga. Closed Tip	•			•	•		
PERIFIX FX 19 Ga. Springwound, Closed Tip (MR unsafe)		•					
PERIFIX FX 19 Ga. Springwound, Open Tip (MR unsafe)			•				
PROCEDURAL DRUGS							
Lidocaine HCl (1%), 5 mL (Skin Wheal)	•	•	•	•	•	•	•
Lidocaine HCl (1.5%), with Epi, 5 mL (Test Dose)	•	•	•	•	•	•	•
DILUENT/LOSS-OF-RESISTANCE							
Sodium Chloride (0.9%), 10 mL	•	•	•	•	•	•	•
DRAPE							
Drape, Clear plastic, 4 in. Fenestrated, 30 in. x 24 in.	•	•	•	•	•	•	•
ACCESSORIES							
Packet Povidone-Iodine Solution	•	•	•	•	•	•	•
Prep Sponges (3 per tray)	•	•	•	•	•	•	•
PERIFIX® 0.2 µm Flat Epidural Filter	•	•	•	•	•	•	•
PERIFIX™ PinPad (securement device for PERIFIX 0.2 µm Flat Epidural Filter)		•	•			•	
Catheter Support Pad	•	•	•	•	•	•	•
Needle Block Foam with Holder	•	•	•	•	•	•	•
ESPOCAN® Docking System						•	
Filter Straw® Particulate Matter Filter	•	•	•	•	•	•	•
Paper Towel, 3-ply	•	•	•	•	•	•	•
4 in. x 4 in. Gauze Sponges (4 per tray)	•	•	•	•	•	•	•
Drug Identification Label Set	•	•	•	•	•	•	•
Threading Assist Guide	•	•	•	•	•	•	•
Clamp Style Catheter Connector	•	•	•	•	•	•	•

SINGLE SHOT NERVE BLOCK NEEDLES, NERVE BLOCK MONITORS AND ACCESSORIES

NEEDLES



STIMUPLEX® A Insulated Needles with 30° Bevel (for use with STIMUPLEX nerve stimulators)

	REF #	Product Code	Units Per Cs
24 Ga. x 1 in. (25 mm) Insulated Needle with Extension Set	4894251		25
22 Ga. x 2 in. (50 mm) Insulated Needle with Extension Set	4894502		25
21 Ga. x 4 in. (100 mm) Insulated Needle with Extension Set	4894260		25
20 Ga. x 6 in. (150 mm) Insulated Needle with Extension Set	4894278		25



STIMUPLEX Ultra 360® Insulated Echogenic Needles with 30° Bevel (for use with Ultrasound and STIMUPLEX nerve stimulators)

22 Ga. x 1 3/8 in. (35 mm) Insulated Echogenic Needle with Extension Set	333640	STIMU32235	25
22 Ga. x 2 in. (50 mm) Insulated Echogenic Needle with Extension Set	333642	STIMU32250	25
22 Ga. x 3 1/8 in. (80 mm) Insulated Echogenic Needle with Extension Set	333644	STIMU32280	25
20 Ga. x 4 in. (100 mm) Insulated Echogenic Needle with Extension Set	333648	STIMU320100	25
20 Ga. x 6 in. (150 mm) Insulated Echogenic Needle with Extension Set	333646	STIMU320150	25



STIMUPLEX D Insulated Needles with 30° and 15° Bevel (for use with STIMUPLEX nerve stimulators)

25 Ga. x 1 3/8 in. (35 mm) Insulated Needle with Extension Set, 15° Bevel	333670	STIMD2535/15	25
23 Ga. x 2 3/4 in. (70 mm) Insulated Needle with Extension Set, 15° Bevel	333672	STIMD2370/15	25
22 Ga. x 2 in. (50 mm) Insulated Needle with Extension Set, 30° Bevel	333689	STIMD2250/30	25
22 Ga. x 3 1/8 in. (80 mm) Insulated Needle with Extension Set, 15° Bevel	333674	STIMD2280/15	25
22 Ga. x 3 1/8 in. (80 mm) Insulated Needle with Extension Set, 30° Bevel	333684	STIMD2280/30	25
22 Ga. x 4 3/4 in. (120 mm) Insulated Needle with Extension Set, 15° Bevel	333676	STIMD22120/15	25



ULTRAPLEX 360® Non-Insulated Echogenic Needles with 30° Bevel (for use with Ultrasound)

22 Ga. x 1 3/8 in. (35 mm) Non-Insulated Echogenic Needle with Extension Set	333638	UPLX32235	25
22 Ga. x 2 in. (50 mm) Non-Insulated Echogenic Needle with Extension Set	333630	UPLX32250	25
22 Ga. x 3 1/8 in. (80 mm) Non-Insulated Echogenic Needle with Extension Set	333632	UPLX32280	25
20 Ga. x 4 in. (100 mm) Non-Insulated Echogenic Needle with Extension Set	333634	UPLX320100	25
20 Ga. x 6 in. (150 mm) Non-Insulated Echogenic Needle with Extension Set	333636	UPLX320150	25

NERVE BLOCK MONITORS & ACCESSORIES



STIMUPLEX Nerve Stimulator Products

STIMUPLEX HNS12 Nerve Stimulator with SENSE*, with four-pin electrode connecting cable	4892098	HNS12	1
STIMUPLEX PEN Percutaneous nerve mapping device (for use with HNS11 or HNS12)	333544	STIMPEN	1
STIMUPLEX HNS11 replacement electrode connecting cable (for use with HNS11 and HNS12 only)	333576	HEC	1



B-Smart™ Injection Pressure Monitor

B-Smart Disposable Manometer for Measuring PNB Injection Pressure	BS015-30		25
---	----------	--	----

*Sequential Electrical Nerve Stimulation

NERVE BLOCK SETS

CONTIPLEX® Sets PERICAN™ Ultra Sets

	Units Per Cs.	5	10	12	10	12	10	12	10	12	10	10								
REF #	Product Code	333697	333696	331766	331670	331672	331757	331758	331695	331691	331693	331694	331698	331767	331768	331640	331641	331642	331771	331773
PROCEDURAL NEEDLES - PERIPHERAL NERVE BLOCK																				
CONTIPLEX Tuohy 18 Ga. x 1½ in. (40 mm), Stimulating																				
CONTIPLEX Tuohy 18 Ga. x 2 in. (50 mm), Stimulating																				
CONTIPLEX Tuohy 18 Ga. x 4 in. (100 mm), Stimulating																				
CONTIPLEX Tuohy 18 Ga. x 6 in. (150 mm), Stimulating																				
CONTIPLEX Ultra 360 18 Ga. x 2 in. (50 mm), Stimulating, Echogenic																				
CONTIPLEX Ultra 360 18 Ga. x 4 in. (100 mm), Stimulating, Echogenic																				
CONTIPLEX Ultra 360 18 Ga. x 6 in. (150 mm), Stimulating, Echogenic																				
CONTIPLEX Straight 18 Ga. x 4 in. (100 mm), Stimulating																				
17 Ga. x 3½ in. (90 cm) Non-Stimulating Tuohy Needle																				
PERICAN™ Ultra 18 Ga. x 3 ½ in. (8.9 cm)																				
PERIPHERAL NERVE BLOCK CATHETERS																				
20 Ga. x 39 in. (100 cm) Polyamide Nylon Catheter with Closed Tip																				
CONTIPLEX Echo 20 Ga. x 22 in. (55 cm) Echogenic Catheter with Stylet, Closed Tip (MR Unsafe)																				
CONTIPLEX Echo 20 Ga. x 22 in. (55 cm) Echogenic Catheter with Stylet, Open Tip (MR Unsafe)																				
CONTIPLEX STIM 20 Ga. x 22 in. (55 cm) Stimulating Catheter with Integrated Wire (MR unsafe)																				
CONTIPLEX FX 19 Ga. x 36 in. (91 cm) Polyamide Springwound FX Catheter, Closed Tip (MR unsafe)																				
CONTIPLEX FX 19 Ga. x 36 in. (91 cm) Polyamide Springwound FX Catheter, Open Tip (MR unsafe)																				
CATHETER-OVER-NEEDLE SYSTEM																				
CONTIPLEX C 19 Ga. x 7⅜ in. (18.7 cm) catheter over 25 Ga. x 7½ in. (190 mm) needle with 30° bevel																				
CONTIPLEX C 19 Ga. x 7⅜ in. (18.7 cm) catheter over 25 Ga. x 7½ in. (190 mm) needle with 15° bevel																				
ACCESSORIES																				
0.2 µm Procedural Filter																				
PinPad™ Securement Device																				
30 cm Connection Tubing																				
Sideport Valve Assembly with 30 in. (76 cm) Extension Tubing																				
Threading Assist Guide																				
Clamp Style Catheter Connector																				
Microbore Extension Set																				

PERICAN NERVE BLOCK SETS



NERVE BLOCK SUPPORT TRAYS

Nerve Block Support Trays

	Units Per Cs.	10			
	Product Code	SNB000U	NB100	NB300	NB400
	REF #	332138	332100	332105	332103
HYPODERMIC NEEDLES					
27 Ga. x 1¼ in.			•	•	•
25 Ga. x 5⁄8 in.				•	•
25 Ga. x 1½ in.		•	•	•	•
18 Ga. x 1½ in.		•	•		
3mL Plastic Luer Lock Syringe with Pre-Attached 18 Ga. x 11/2 in. Needle				•	•
SYRINGES					
1 mL Plastic – Luer Slip			•	•	•
3 mL Plastic – Luer Lock			•		
5 mL Plastic – Luer Lock		•	•		
10 mL Plastic – Luer Lock				•	•
20 mL Plastic – Luer Lock		•	•		
20 mL Plastic – Luer Lock (2 per tray)				•	•
PROCEDURAL DRUGS					
Lidocaine HCl (1%), 5 mL (skin wheal)				•	•
DRAPES					
Clear Drape with Circular Fenestration				•	
Clear Drape with Rectangular Fenestration			•		•
ACCESSORIES					
Prep Sponges (3 per tray)			•		
ChloroPrep® with Orange Tint-3 mL Applicator		•		•	•
Prep Solution Well		•			
Skin Marker			•	•	•
Paper Ruler				•	•
Filter Straw® Particulate Matter Filter				•	•
CIV-Flex™ Transducer Cover and Ultrasound Gel – 30 cm		•			
CIV-Flex Transducer Cover and Ultrasound Gel – 147 cm					•
0.2 µm Procedural Filter				•	•
PinPad™ (Securement Device for 0.2 µm Procedural Filter)				•	•
Benzoin Swabstick				•	•
Steri-Strip® Closure (card of 3)				•	•
Tegaderm® 4 in. x 4¾ in. Dressing				•	•
Needle Block Foam with Holder		•	•	•	•
4 in. x 4 in. Gauze Sponges (4 per tray)		•	•	•	•
Drug Identification Label Set				•	•
Paper Towel, 3-ply				•	•
Paper Towel, 3-ply (2 per tray)		•			

PROCEDURAL NERVE BLOCK TRAYS

STIMUPLEX® Ultra 360® Single Shot Nerve Block Trays
CONTIPLEX® Continuous Nerve Block Trays

Units Per Cs. 10

Product Code

REF #

PROCEDURAL NEEDLES – PERIPHERAL NERVE BLOCK	
STIMUPLEX Ultra 360 22 Ga. x 2 in. (50 mm)	
STIMUPLEX Ultra 360 20 Ga. x 4 in. (100 mm)	
CONTIPLEX Tuohy 18 Ga. x 2 in. (50 mm)	
CONTIPLEX Tuohy 18 Ga. x 4 in. (100 mm)	
CONTIPLEX Ultra 360 18 Ga. x 2 in. (50 mm)	
CONTIPLEX Ultra 360 18 Ga. x 4 in. (100 mm)	
HYPODERMIC NEEDLES	
27 Ga. x 1¼ in.	
25 Ga. x 1½ in.	
18 Ga. x 1½ in.	
3mL Plastic Luer Lock Syringe with Pre-Attached 18 Ga. x 1½ in. Needle	
SYRINGES	
3 mL Plastic – Luer Lock	
5 mL Plastic – Luer Lock	
10 mL Plastic – Luer Lock	
10 mL Plastic – Luer Slip	
20 mL Plastic – Luer Lock	
PERIPHERAL NERVE BLOCK CATHETERS	
CONTIPLEX STIM Stimulating Catheter–20 Ga., Closed Tip with Integrated Wire (MR unsafe)	
CONTIPLEX Polyamide Catheter–20 Ga., Closed Tip	
CONTIPLEX Echo 20 Ga. x 22 in. (55 cm) Echogenic Catheter with Stylet, Closed Tip (MR Unsafe)	
PROCEDURAL DRUGS	
Lidocaine HCl (1%), 5 mL (skin wheal)	
Lidocaine HCl (1.5%), with Epi, 5 mL (test dose)	
DRAPES	
Clear Drape with Circular Fenestration	
Clear Drape with Rectangular Fenestration	
ACCESSORIES	
ChloraPrep® with Orange Tint–3mL Applicator	
Prep Solution Well	
Paper Ruler	
Filter Straw® Particulate Matter Filter	
CIV-Flex™ Transducer Cover and Ultrasound Gel – 30 cm	
CIV-Flex Transducer Cover and Ultrasound Gel – 147 cm	
Sideport Valve Assembly with Extension Tubing	
0.2 µm Procedural Filter	
PinPad (Securement Device for 0.2 µm Procedural Filter)	
Benzoin Swabstick	
Steri-Strip® Closure (card of 3)	
Steri-Strip Closure (card of 3; 2 cards per tray)	
Tegaderm® 4 in. x 4¾ in. Dressing	
Tegaderm 4 in. x 4¾ in. Dressing (2 per tray)	
NeedleVISE®	
Needle Block Foam with Holder	
4 in. x 4 in. Gauze Sponges (4 per tray)	
Paper Towel, 3 Ply	
Drug Identification Label Set	
Threading Assist Guide	
Clamp Style Catheter Connector	
30 cm Extension Set	

REF #	Product Code	Units Per Cs.
332135	SNB200U3	10
332137	SNB400U3	10
331701	CNB200TK	10
331706	CNB400TK	10
332127	CNB200TKU	10
332128	CNB400TKU	10
332142	CNB200TU3KU	10
332143	CNB400TU3KU	10
331755	SCNB2TK	10
331756	SCNB4TK	10
332170	CNB400TU3KUE	10






STIMUPLEX SINGLE SHOT NERVE BLOCK TRAYS

CONTIPLEX CONTINUOUS NERVE BLOCK TRAYS

NeedleVISE is a registered trademark of Atrion Medical Products, Inc. ChloraPrep is a registered trademark of CareFusion Corporation. Tegaderm and Steri-Strip are registered trademarks of 3M. CIV-Flex is a trademark of CIVCO Medical Instruments Co., Inc.

Xperius® Ultrasound System

Designed around you to seamlessly fit into your care area

Description		REF #	Units Per Cs
Xperius Cart System with 1 transducer L12-4	<ul style="list-style-type: none"> • Type of unit: Compact Cart based Point of Care Ultrasound System with fully articulating arm and Linear L12-4 transducer 	9896 0547 0201-02	1
Xperius Cart System with printer and 1 transducer L12-4	 <ul style="list-style-type: none"> • Common to both Xperius Cart System and Xperius Mobile Workspot • Imaging mode keys: 2D, Color Doppler, Color Power Angio (CPA), M-Mode • Image Processing: AutoSCAN, XRES, Duplex Imaging, 5x Pan/Zoom, Dynamic range up to 170 db (full time input) • Needle Visualization Software: NeedleViz • Image enhancement controls: SonoCT, XRES • Patient-specific optimization keys: AutoSCAN, iSCAN • Connectivity: HDMI, 5 USB ports on cart, 1 USB port & 1 ethernet port on tablet, Wi-Fi, DICOM networking and Philips Remote Services • Display: 15.6 inch (39.6 cm) 1920 x 1080 LED-backlit IPS display • Dimensions: 52.2 in. (132.5 cm) x 22 in. (56 cm) x 26 in. (66 cm) (H/W/D) • Battery operating time: over 3 hrs. of scanning on full charge 	9896 0547 0211-02	1
Xperius Tablet System with 1 transducer L12-4	 <ul style="list-style-type: none"> • Type of unit: Point of Care Tablet Ultrasound System with Kickstand and Linear 12-4 transducer • Dimensions: 11.0 in. (28.0 cm) x 16.5 in. (42.0 cm) x 1.9 in. (5.0 cm) (H/W/D) 	9896 0547 6671-02	1
Xperius Accessory Items			
Single Transducer Linear	 <ul style="list-style-type: none"> • L12-4 broadband linear array transducer • 12 to 4 MHz extended operating frequency range • Aperture size: 1.3 in. (34 mm) • High resolution imaging for shallow applications • Center line marker • USB transducer with replaceable cable 	9896 0545 1171-02	1
Single Transducer Curved	 <ul style="list-style-type: none"> • C5-2 broadband curved array transducer • 5 to 2 MHz extended operating frequency range • 2 in. (50 mm) radius of curvature • High resolution imaging for deeper applications • Center line marker • USB transducer with replaceable cable 	9896 0545 1181-02	1
Peripherals	 <ul style="list-style-type: none"> • Printer Sony UP-D711MD • High-quality digital graphic printer • Monochrome, black/white, A7 paper size 	9896 0545 2822-02	1

Take Away the Guesswork.



Objective injection pressure monitoring with the **B-Smart™** pressure monitor provides a visual indication of high injection pressure.

"High injection pressures at the onset of injection may indicate an intraneural needle placement and lead to severe fascicular injury and persistent neurologic deficits."¹

- *The first disposable objective injection pressure monitor for peripheral nerve blocks*
- *Helps standardize nerve block injections and documentation*

Find out more: **PNBmonitoring.com** and **BBraunUSA.com/BSmart**

B. Braun Medical Inc. | Bethlehem, PA
800-277-2862 | BBraunUSA.com

1. Hadzic A, Dilberovic F, Shah S, Kulenovic A, Kapur E, Zaciragic A, Cosovic E, Vuckovic I, Divanovic K-A, Mornjakovic Z, Thys DM, Santos AC: Combination of intraneural injection and high injection pressure leads to fascicular injury and neurologic deficits in dogs. *Reg Anesth Pain Med* 2004; 29:417-23.