

IHE Integration Statement			
Vendor	Product Name	Version	Date
B. Braun Medical	DoseLink	3.0	9/9/19

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Alert Communication Management	Alert Reporter	None
Consistent Time	Time Client	None
Device Enterprise Communication	Device Observation Filter	None
Device Enterprise Communication	Device Observation Reporter	None
Infusion Pump Event Communication	Device Observation Reporter	None
Medical Equipment Management Device Management Communication	Device Management Information Reporter	None
Medical Equipment Management Location Services	Location Observation Consumer	None
Patient Administration Management	Patient Demographics Consumer	Merge
Patient Administration Management	Patient Encounter Consumer	Inpatient / Outpatient Encounter Management
Point-of-Care Infusion Verification	Infusion Order Consumer	None

Internet address for vendor's IHE information :

[https://www.bbraunusa.com/content/dam/b-braun/us/website/products-and-therapies/infusion\\_pumps/pumps\\_docs/B%20Braun%20DoseLink%20-%20IntegrationStatement%20ver%203.0.pdf.bb-.88146454/B%20Braun%20DoseLink%20-%20IntegrationStatement%20ver%203.0.pdf](https://www.bbraunusa.com/content/dam/b-braun/us/website/products-and-therapies/infusion_pumps/pumps_docs/B%20Braun%20DoseLink%20-%20IntegrationStatement%20ver%203.0.pdf.bb-.88146454/B%20Braun%20DoseLink%20-%20IntegrationStatement%20ver%203.0.pdf)

Link to general information on IHE : [www.ihe.net](http://www.ihe.net)