

Improving Medication Infusion Safety Through Technology and Practice

Rosanne Raso, RN, MS, SVP Nursing
Joan Velletri, RN, MSN, AVP Nursing Education
Steven DiCrescento, RPh, MS Lutheran Medical Center

Presented at the American Organization of Nurse Executives 42nd Annual Meeting, April 2009



Objective

This study was conducted to determine the incidence of averted infusion-related medication events after implementing smart pump technology with practice changes and to identify potential areas of risk for continuous process improvement.


Methods

- 275 smart pumps were implemented in 2006.
- P&Ts were established to support smart pump use, to require independent double-checks for all dose overrides, and to provide random compliance checks.
- Retrospective analyses of pump data logs were conducted at three and 27 months post implementation. Data logs were analyzed for incidence of averted events, dose corrections and overrides.
- Based on analysis findings, process improvements were implemented.

Results

One year after implementation of smart pump technology and evidence-based protocols, we experienced a 57% reduction in infusion-related medication events. By year three, our infusion-related medication events were zero.

Infusion-related Medication Events 2005–2008



Analysis of the data logs demonstrated there were a total of 8 averted medication errors out of 1,536 programmed doses, representing an error rate of only 0.5%. Averted errors were associated with midazolam, propofol, vasopressin and amiodarone. No errors were associated with insulin or heparin.

Additionally, we found:

- 100% compliance with dose mode
- 100% wt-based ordering practice
- Low incidence of programming error (0.5%)

Impact of Data Analysis

- Identified dosing trends and issues
- Provided practice recommendations
- Provided drug library recommendations
- Documented 0 alerts with insulin or heparin
- Documented low incidence of error
- Documented positive results = positive staff

Discussion

Our concrete reduction in infusion-related medication events was due to:

- Easy-to-use smart pump technology
- Multi-disciplinary development of a drug library that reflected true clinical practice
- Adoption of protocols for dose mode use, override verification and compliance audits
- Ongoing education of nursing and resident staff
- Vendor support for data analysis to identify improvement opportunities and library changes
- Total engagement of management and staff
- Sharing of positive results

References

Raso, R., Velletri, J., DiCrescento, S. (2007). Making the most of data for patient safety. *Patient Safety & Quality Healthcare*, May/June 2007.

Acknowledgements

Lutheran Medical Center utilizes B. Braun Medical Inc's Outlook® Safety Infusion System and DoseTrac™ Analysis Services.

