



HRS Control Panel Spare Parts

Hot Rinse Smart: 10-20



	Article Number	Description
1	34200050	Control circuit board power until LT1plus
	34200051	Micro-fuse 5x20mm 2.0A, 10 pcs
2	34200177	Emergency module 24 VDC BG 5924.48
3	34200015	Relay for plug 24V DC 8 pin/1 pc
4	34200016	Assembly for relay base, 1 pc
5	34200196	Input/Output module 14pol
6	34200197	Motor protection switch PKZM0-2.5
7	34200056	Fault current circuit breaker 4-pole UL
8	34200060	Main contactor DILM 12-10
9	34200061	Relay for plug (230V AC) 8 pin/1 pc
10	34200059	Relay base derive resistor, 1 pc
11	34200049	Relais semiconductor 20A 600VAC



Article Number	Description
12	34200049 Control board CPU II with program (please specify for HRS or RO)
13	34200048 Lithium battery coin cell 3V
14	34200104 LED element with adapter, 2 pcs
15	34200065 Connecting cable control board CPU
16	34200052 Power unit 10A Pri 85 264VAC Sec 24V
17	34200195 Mains filter FAI DB 2150 ZC
18	34200175 Micro-fuse 5x20mm 6.3A delayed, 10 pcs
	34200180 Flashlight element red 24 VDC