




Solcart B®
Biocarbonate Cartridge for use with the
Dialog+® Hemodialysis System

Solcart B®

A better dialysis experience doesn't just come from a box. It's built from better products, better support, and better partnerships.



DESIGNED FOR EASE OF USE



CONVENIENT TO STORE



LIGHTWEIGHT



POLYPROYLENE CARTRIDGE HOUSING



	Solcart B 650g	Solcart B 760g	Solcart B 1100g
Item #	802	803	805
Units/Carton	10	10	8
Dialysate Flow Rate	Run Time		
500 ml/min	6.2 hours	7.3 hours	10.6 hours
600 ml/min	5.2 hours	6.1 hours	8.8 hours
700 ml/min	4.4 hours	5.2 hours	7.5 hours
800 ml/min	4.1 hours	4.8 hours	6.9 hours

Bicarbonate setting of 35 mmol/L. Run time values based on testing conducted using 760g Solcart B Cartridge. Data extrapolated for other sizes and represent approximate run time values.

B. Braun Medical Inc. | Bethlehem PA | braunusa.com/solcart

Rx Only. ©2019 B. Braun Medical Inc., Bethlehem, PA. All Rights Reserved. AVT.033_12/19