

BULLETIN

B|BRAUN
SHARING EXPERTISE



B. Braun Medical Renal Therapies Quarterly Newsletter

The Benefits of Upgrading Dialog+® Hemodialysis System to 9.1B Comfort Software

Automatic Blood Pressure Module: Designed for Comfort

- Ability to reset without turning machine on/off
- Measures during inflation (only with M2015 ABPM)

Increased Regulation: additional sensor added to monitor dialysis fluid temperature

Optimized Flow Control: added balance chamber calibration

Optimized Dialysate Flow: air is purged seamlessly from DF filter every 150 seconds

- Alarms if air separator sensor opens more than 2 times in one hour



Ready to upgrade? Call customer service at **1-800-848-2066**.

Technical Website Developments



Service Portal and LMS are Live!

The Service Portal and Learning Management System are now available for you to access troubleshooting documents and equipment manuals.

[Learn More](#)

Non-Emergency Question?

If you have a non-emergency question, you no longer have to call us! You can contact our knowledgeable technicians by sending a message online.

[Contact Us](#)

Need to Register for Classes?

You can now register for all our available, upcoming Dialog+® hemodialysis system and Aquaboss® water treatment system technical classes online! Find itineraries for the classes, lodging and transportation details, and our yearly class schedule.

[Learn More](#)

Missed last quarter's bulletin? View all our archived new bulletins at bbraunusa.com/renalbulletin.

Our mailing address is:

824 12th Avenue
Bethlehem, PA 18018

Want to change how you receive these emails?

You can [update your preferences](#) or [unsubscribe from this list](#).