

# BULLETIN

**B|BRAUN**  
SHARING EXPERTISE



## B. Braun Medical Renal Therapies Quarterly Newsletter

---

### Clinical Support Update






During the COVID-19 crisis, our clinical education team continued to provide the necessary Dialog+ hemodialysis system training for clinics through online webinars. Currently, our clinical support team is providing both online education as well as clinical in-service training.

---

### Pulse Therapy on Dialog+ Hemodialysis Machines

The need for CRRT continues to increase due to COVID-19. Pulse therapy on the Dialog+ hemodialysis machines can provide an alternative option to using a CRRT machine.

Three profiles are designed for extended treatment times. These profiles allow for a gentler removal of fluid with periods of pure ultrafiltration that allow for brief periods of solute re-saturation. The profiles can also be edited, individualizing the treatment to best meet the needs of the patient.



## Clinical Library Update

We now have FMC and Dialog Streamline bloodline troubleshooting videos live on our website!

[View Now](#)

## Upcoming Webinars

- 7/14/20: Introducing Dialog+
- 7/15/20: Adimea, Accurate Dialysis Measurement
- 7/16/20: Basics of ESRD and Hemodialysis
- 8/18/20: Troubleshooting Techniques for Dialog+
- 8/19/20: Profiling Capabilities of Dialog+
- 8/20/20: Dialysis Quality Water and Ultrapure Dialysate
- 9/15/20: Introducing Dialog+
- 9/16/20: Adimea, Accurate Dialysis Measurement
- 9/17/20: Basics of ESRD and Hemodialysis

[Sign Up](#)

## Customer Service Holidays

The B. Braun Renal Therapies customer service team will be out of the office and unavailable to take orders through the phone on July 2-3, 2020 and September 7, 2020.

*Missed last quarter's bulletin? View all our archived bulletins  
at [bbraunusa.com/renalbulletin](http://bbraunusa.com/renalbulletin).*

---

*Copyright © B. Braun Medical Inc., All rights reserved.*

**Our mailing address is:**

824 12th Avenue  
Bethlehem, PA 18018

Want to change how you receive these emails?  
You can [update your preferences](#) or [unsubscribe from this list](#).