

## Concentrate Supply

Part Number	Part Name	Technical Description	Qty.
<b>Wands (Central &amp; Non-Central)</b>			
34562260	Male Connector for Suction Rod	Coupling w/ Red-Blue Caps [Non-Central]	1
38245442	Seal Nipple	Wall Box Connector for Concentrate [Central]	1
34560920	Suction Coupling	Coupling w/ Red-Blue Caps [Central]	5
34570438	O-ring [8.0 x 1.25]	O-Ring for Suction Coupling	5
34560939	Suction Nipple	Attaches to PE Tube [Central]	5
34570004	O-ring [9.3 x 2.4]	O-Rings for Suction Nipple [Central]	10
3451728A	PE Tube [4x1]	Hard Plastic Tubing That Goes into Jugs	3
34565280	Silicone tubing [3.2 x 1.6]	Attaches PE tube to Suction Nipple	5
S20110	Delivery Coil – BICARB	Blue Coil Tubing [Central]	1
S20100	Delivery Coil – ACID	Red Coil Tubing [Central]	1
34565213	Silicone tubing [2.9 x 3]	Silicone Tubing for Wands [Non-Central]	20
34561498	Disposable Filter [2.6 mm]	Filters for Wands	1
3451726A	Feed-Through Connector	Connector Under Front Door for Wands	1
<b>Retainer Tubes</b>			
3456077A	Retainer Tubes (green) for Suction Rods	Green Tubes that the PE Tubing Rests In	1 pair
34565345	Silicone Tubing [12 x 2.5]	Holds Retainer Tubes on Bracket	1m
34514120	Retainer for Suction Rods	White Housing Wand Couplings Attach To	1
34517138	Reed Sensor [with contacts & O-ring 3.68 x 1.78]	Sensor To Make Sure Wands are Attached	1
<b>Misc.</b>			
604010	Concentrate Container 10 L NATURE 365G	Empty Jug no Label	1
34514562	Suction Tubing "YELLOW" Complete	Yellow Lid, Tubing, and Metal Filter	1
604033	Canister Cap, Yellow	Disinfection Jug Lid	1
39971260	Disinfection Tubing [3.6 x 5.8]	Tubing only	1
34512675	Suction Pipe	Metal Filter for Citric Jug	1
S20066	Label Disinfection Jug	Label for Disinfection Jug	1
34517774	Canister Lid, Blue	Blue Lid w/ Hole for Bicarb	1
604032	Stopper For Canister Blue 45 OV	Solid Blue Lid for Bicarb	1
S20064	Label BICARB Jug	Blue for Bicarb	1
34517766	Canister Lid, Red	Red Lid w/ Hole for Acid	1
604029	Stopper for Canister Red 45OV AX-D	Solid Red Lid for Acid	1
S20065	Label Acid Jug	RED for Acid	1

Part Numbers Current as of 11/29/2018