

Top & Low Level

Part Number	Part Name	Technical Description	Qty.
Boards			
34560610	LX 800 Motherboard	TLC	1
34519327	Battery 2032 / 3 Volt	Battery for Mother Board	1
3456088A	Compact Flash Card, 1 GB	CF Card	1
34562168	Basic Board V3	LLBB	1
34562150	Power Board (Valves)	PBV	1
34560750	Power Board (Motors)	PBM	1
3456074B	Digital Board V4	DB for ≤ SW 9.12	1
34562141	Digital Board V5	DB for ≥ SW 9.18	1
34560734	Analog Board	AB	1
3451487E	ABPM Module M2015	ABPM	1
3456127A	Temperature Switch	Turns on Fan When Hot	1
Ribbon Cables			
34562225	Ribbon Cable 26 Conductors (FPB-MB)	Motherboard → Front Panel Board	1
34560912	Ribbon Cable 12 Conductors (FPB-OSD)	Front Panel Board → OSD	1
34562221	Ribbon Cable 10 Conductors (TCB-MB)	Motherboard → Touch Controller Board	1
34560971	Ribbon Cable (FPB-BIB)	Front Panel Board → Backlight Inverter Board	1
34562222	Ribbon Cable 16 Conductors (FPB-BB)	Front Panel Board → Basic Board (P14)	1
34562224	Ribbon Cable Level Regulation	Level Regulation → Basic Board (P11)	1
34562228	Ribbon Cable for SAD	Basic Board → Safety Air Detector	1
Communications			
7102231	WAN - BSL sw. 9.XX	Entire BSL kit	1
34514554	Bedside Link Board / BSL	BSL Board only	1
7107232	External MFC-Ian Connection Cable	External cable for BSL	1
7107218	DCI Cable - All	DCI	1
710725A	ReBell Site License	For DCI	1
7107251	Server SW ReBell Tool	For DCI	1
7105230	Card Reader - All	Full Card Reader Kit	1
34560548	Card Reader Board	Replacement Board for Card Reader	1
7105232	Card Pack	Pack of 5 Patient ID Cards	5
Misc.			
34514929	Loudspeaker+	Loudspeaker by Motherboard	1

Part Numbers Current as of 11/29/2018