

# Sharps Safety Engineered Epidural & Spinal Products



Safety engineered components for you and your patients

# Why choose B. Braun for your sharps safety needs?

As the worldwide leader in Regional Anesthesia products, you can feel confident and secure in using quality sharps safety products designed to provide a consistent tactile feel and comply with OSHA and JCAHO guidelines.

*"Risks to anesthesia personnel could be greatly reduced by eliminating the use of unnecessary sharps devices, reducing the use of necessary sharp devices to a minimum, and increasing the use of safety devices, including needleless or shielded needle devices and protective barriers."<sup>1</sup>*

Contaminated percutaneous injuries (CPI) were analyzed during a 2-year period at 11 hospitals. There were 138 CPIs, 22 occurred during regional anesthesia procedures. After data was corrected for CPI underreporting, the table below shows the incidence per 100,000 units used.<sup>1</sup>

Procedure	Incidence
Epidural	7.91
Spinal	10.09
Peripheral Nerve Block	187.48

1. Greene, ES, et. al. "Multicenter study of contaminated percutaneous injuries in anesthesia personnel". Anesthesiology 1998; 89: 1362-1372



Filter Straw® particulate matter filter for needle-free aspiration of fluids from ampules

MICRO PIN® for needle-free aspiration of fluids from vial

NeedleVise® designed to protect up to five needles that cannot be removed once placed

25G x 1 in. and 22G x 1.5 in. SafetyGlide™ hypodermic needles designed to fit on conventional syringes. The shielding mechanism requires a single finger stroke to activate.



Perifix® 0.2µm Flat Filter filters out particulates and PinPad secures the catheter connector to the patient

Packet Povidone-Iodine Solution and Prep Sponges

# Ordering Information

Reorder Number	Product Code	Description	Case Qty.
<b>Sharps Safety Engineered Epidural Products</b>			
332252	SD18TDY	18G x 3.5 in. Perifix® Tuohy Needle with 25G x 1 in. SafetyGlide™ hypodermic needle, 22G x 1.5 in. SafetyGlide hypodermic needle, 5cc glass luer lock LOR syringe, 1cc plastic luer slip syringe, 3cc plastic luer lock syringe, 10cc plastic luer lock syringe, 20cc plastic luer lock syringe, Lidocaine 1% (Skin Wheal), Sodium Chloride 0.9%, Povidone-Iodine Solution, prep sponges, clear drape, NeedleVise®, MICRO PIN® dispensing pin, Filter Straw®, gauze, paper towel, syringe labels	10
332075	CE17TKFCY	17G x 3.5 in. Perifix Tuohy Needle with 25G x 1 in. SafetyGlide hypodermic needle, 22G x 1.5 in. SafetyGlide hypodermic needle, 5cc glass luer slip LOR syringe, Perifix FX 19G closed tip catheter, Lidocaine 1% (Skin Wheal), Lidocaine 1.5% with Epinephrine, Sodium Chloride 0.9%, blue drape, clear drape, Povidone-Iodine Solution, prep sponges, Perifix 0.2µm flat epidural filter, Perifix PinPad, catheter pad, Tegaderm® 4 in. x 4¾ in. dressing, NeedleVise, 3cc plastic luer lock syringe, 20cc plastic luer lock syringe, MICRO PIN dispensing pin, Filter Straw, gauze, paper towel, syringe labels	10
332266	CE18TKY	18G x 3.5 in. Perifix Tuohy Needle with 25G x 1 in. SafetyGlide hypodermic needle, 22G x 1.5 in. SafetyGlide hypodermic needle, 5cc glass luer lock LOR syringe, Perifix 20G closed tip nylon catheter, Lidocaine 1% (Skin Wheal), Lidocaine 1.5% with Epinephrine, Sodium Chloride 0.9%, blue drape, clear drape, Povidone-Iodine Solution, prep sponges, Perifix 0.2µm flat epidural filter, Perifix PinPad, catheter pad, Tegaderm 4 in. x 4¾ in. dressing, NeedleVise, 3cc plastic luer lock syringe, 20cc plastic luer lock syringe, MICRO PIN dispensing pin, Filter Straw, gauze, paper towel, syringe labels	10
<b>Sharps Safety Engineered Spinal Products</b>			
333236	S22LKAY	22G x 3.5 in. Spinocan® Needle with 25G x 1 in. SafetyGlide hypodermic needle, 22G x 1.5 in. SafetyGlide hypodermic needle, 3cc plastic luer lock syringe with pre-attached 1.75 in. Filter Straw, Epinephrine 0.1%, Lidocaine 5% with Dextrose 7.5%, Lidocaine 1% (Skin Wheal), Povidone-Iodine Solution, prep sponges, clear drape, NeedleVise, 5cc plastic spinal syringe, paper towel, gauze, syringe labels	10
333239	S25BKAY	25G x 3.5 in. Spinocan Needle with 18G x 1.25 in. introducer needle, 25G x 1 in. SafetyGlide hypodermic needle, 22G x 1.5 in. SafetyGlide hypodermic needle, 3cc plastic luer lock syringe with pre-attached 1.75 in. Filter Straw, Epinephrine 0.1%, Bupivacaine 0.75% with Dextrose 8.25%, Lidocaine 1% (Skin Wheal), Povidone-Iodine Solution, prep sponges, clear drape, NeedleVise, 5cc plastic spinal syringe, paper towel, gauze, syringe labels	10
333864	P25LKAY	25G x 3.5 in. Pencan® Needle with 20G x 1.25 in. introducer needle, 25G x 1 in. SafetyGlide hypodermic needle, 22G x 1.5 in. SafetyGlide hypodermic needle, 3cc plastic luer lock syringe with pre-attached 1.75 in. Filter Straw, Epinephrine 0.1%, Lidocaine 5% with Dextrose 7.5%, Lidocaine 1% (Skin Wheal), Povidone-Iodine Solution, prep sponges, clear drape, NeedleVise, 5cc plastic spinal syringe, paper towel, gauze, syringe labels	10
333865	P25BKAY	25G x 3.5 in. Pencan Needle with 20G x 1.25 in. introducer needle, 25G x 1 in. SafetyGlide hypodermic needle, 22G x 1.5 in. SafetyGlide hypodermic needle, 3cc plastic luer lock syringe with pre-attached 1.75 in. Filter Straw, Epinephrine 0.1%, Bupivacaine 0.75% with Dextrose 8.25%, Lidocaine 1% (Skin Wheal), Povidone-Iodine Solution, prep sponges, clear drape, NeedleVise, 5cc plastic spinal syringe, paper towel, gauze, syringe labels	10

Call **1-800-227-2862** for your B. Braun representative and for more information on our full line of regional anesthesia products.

