



Clik-FIX® Epidural/PNB Catheter Securement Device

Designed to be low-profile,
soft and secure

The Clik-FIX Epidural/PNB Catheter Securement Device is built to be low-profile, soft and secure to help ensure patient comfort for epidural and peripheral nerve block procedures.



ECS-1000B*

* Follow your facility's protocols for catheter dressing

B | BRAUN
SHARING EXPERTISE

Securement

- All-in-one version is engineered to protect the insertion site and secure the catheter using only soft, low-profile materials
- Universal securement designed to secure 18-24 gauge epidural and PNB catheters

Comfort

- Made with patient comfort in mind
- Adhesive designed to be skin friendly, yet strong
- Foam base pad cushions under the catheter tubing to help reduce the risk of pressure injuries

Design

- All-in-one version is built to be easy-to-use by replacing the need to use a separate securement device and transparent dressing
- Built-in tension loop designed to reduce the risk of dislodgement and movement of the catheter
- Helps to reduce the risk of epidural and PNB catheter related complications

Seven reasons to switch

1. **All-in-one:** bordered bio-occlusive dressing, securement and cushion. Designed to be less bulky and easy to place.
2. **Bordered dressing:** provides a clear view of the insertion site and surrounding skin.
3. **Securement channels:** feature strategically placed contact points to guide and hold the catheter in place.
4. **Foam base pad:** designed to cushion the skin and help reduce the risk of pressure injuries.
5. **Solvent-free adhesive:** designed to be skin friendly, comfortable, stable and easy to remove.
6. **Cut-out:** allows for clear visualization of the insertion site.
7. **Soft, low profile design:** no hard plastic parts to catch on objects, bruise skin or cause discomfort.



ECS-2000B

Ordering Information



ECS-1000B

Clik-FIX® Epidural/PNB Catheter Securement Device. Standalone.
400 units per case (each case contains 8 boxes with 50 units)



ECS-2000B

All-in-one Clik-FIX Epidural/PNB Catheter Securement Device.
400 units per case (each case contains 8 boxes with 50 units)

B. Braun Medical Inc. | Bethlehem, PA | 800-227-2862 | bbraunusa.com