

# IV Sets and Accessories

CATALOG | APRIL 2026



# TABLE OF

# CONTENTS

|                                                                                       |                              |
|---------------------------------------------------------------------------------------|------------------------------|
| <a href="#">Stopcocks and Manifolds</a>                                               | <a href="#">1–2</a>          |
| <a href="#">Anesthesia Extension Sets</a>                                             | <a href="#">3–5</a>          |
| <a href="#">Anesthesia IV Administration Sets</a>                                     | <a href="#">6–7</a>          |
| <a href="#">Blood Administration Sets</a>                                             | <a href="#">8</a>            |
| <a href="#">Basic IV Administration Sets</a>                                          | <a href="#">9</a>            |
| <a href="#">Extension Sets</a>                                                        | <a href="#">10–12</a>        |
| <a href="#">CARESAFE™ IV Administration Sets and Extension Sets</a>                   | <a href="#">13–14</a>        |
| <a href="#">STEADYCare™ Extension Sets</a>                                            | <a href="#">15</a>           |
| <a href="#">Rate Flow® Regulator II Administration Sets</a>                           | <a href="#">16</a>           |
| <a href="#">Irrigation/Urology Sets</a>                                               | <a href="#">17–18</a>        |
| <a href="#">Measured Volume Solution Burette Sets</a>                                 | <a href="#">19</a>           |
| <a href="#">Special Application Accessories</a>                                       | <a href="#">19</a>           |
| <a href="#">Injection Caps and Replacement Caps</a>                                   | <a href="#">20</a>           |
| <b><a href="#">CARESITE® Needleless Connector</a></b>                                 | <b><a href="#">21–22</a></b> |
| <a href="#">CARESITE Needleless Connector</a>                                         | <a href="#">21</a>           |
| <a href="#">CARESITE Small Bore Extension Sets</a>                                    | <a href="#">21–22</a>        |
| <a href="#">CARESITE Standard Bore Extension Sets</a>                                 | <a href="#">23</a>           |
| <a href="#">Filterflow® Filtered Extension Sets with CARESITE Injection Sites</a>     | <a href="#">24</a>           |
| <a href="#">Filterflow Filtered Administration Sets with CARESITE Injection Sites</a> | <a href="#">24</a>           |
| <a href="#">Gravity IV Administration Sets with CARESITE Injection Sites</a>          | <a href="#">24–25</a>        |








## **[CARESITE Micro® Needleless Connector](#)** . . **[26–28](#)**







|                                                                     |                    |
|---------------------------------------------------------------------|--------------------|
| <a href="#">CARESITE Micro Needleless Connector</a>                 | <a href="#">26</a> |
| <a href="#">CARESITE Micro Small Bore Extension Sets</a>            | <a href="#">26</a> |
| <a href="#">CARESITE Micro Standard Bore Extension Sets</a>         | <a href="#">27</a> |
| <a href="#">CARESITE Micro Small Bore Bifurcated Extension Sets</a> | <a href="#">28</a> |
| <a href="#">CARESITE Micro STEADYCare Extension Sets</a>            | <a href="#">28</a> |

## **[SAFELINE® Needleless Connector](#)** . . . . . **[29–32](#)**









|                                                                                                                     |                    |
|---------------------------------------------------------------------------------------------------------------------|--------------------|
| <a href="#">Extension Sets with SAFELINE Split Septum Injection Sites</a>                                           | <a href="#">29</a> |
| <a href="#">Filterflow Filtered Add-on Extension Set with SAFELINE Split Septum Injection Site</a>                  | <a href="#">30</a> |
| <a href="#">Basic Administration Sets with SAFELINE Split Septum Injection Site</a>                                 | <a href="#">30</a> |
| <a href="#">Secondary Set for use with SAFELINE Split Septum Injection Site</a>                                     | <a href="#">30</a> |
| <a href="#">ADDitIV® Primary Administration Sets with Backcheck Valve and SAFELINE Split Septum Injection Sites</a> | <a href="#">31</a> |
| <a href="#">Filterflow Filtered Gravity Administration Sets with SAFELINE Split Septum Injection Sites</a>          | <a href="#">31</a> |
| <a href="#">SAFELINE Needleless Split Septum Injection Site</a>                                                     | <a href="#">32</a> |
| <a href="#">Accessories for SAFELINE Split Septum Injection Site</a>                                                | <a href="#">32</a> |

# Stopcocks and Manifolds

| Description                                                                         |                                                                                                                                                                                                   | REF Number | Case Quantity |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| DISCOFIX® Stopcocks                                                                 |                                                                                                                                                                                                   |            |               |
|    | One-way stopcock with female Luer lock port and SPIN-LOCK® Connector. Not made with DEHP or natural rubber latex.<br>Stopcock volume: 0.20 mL. Length: 1 in. (2.54 cm.).                          | 455980     | 100           |
|    | Three-way stopcock with two female Luer lock ports and male Luer slip connector. Not made with DEHP or natural rubber latex.<br>Stopcock volume: 0.26 mL. Length: 1 in. (2.54 cm.).               | 455991     | 100           |
|    | Three-way stopcock with two female Luer lock ports and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex.<br>Stopcock volume: 0.26 mL. Length: 1 in. (2.54 cm.).                    | 456003     | 100           |
|    | Four-way stopcock with two female Luer lock ports and male Luer slip connector.<br>Not made with DEHP or natural rubber latex. Stopcock volume: 0.26 mL. Length: 1 in. (2.54 cm.).                | 456006     | 100           |
|    | Four-way stopcock with two female Luer lock ports and SPIN-LOCK Connector.<br>Not made with DEHP or natural rubber latex. Stopcock volume: 0.26 mL. Length: 1 in. (2.54 cm.).                     | 456020     | 100           |
|   | Two-gang, three-way stopcocks with three female Luer lock ports and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Stopcock volume: 0.86 mL. Length: 3 in. (7.62 cm.).          | 455992     | 100           |
|  | High-flow, four-way stopcock, two female Luer lock ports and SPIN-LOCK Connector. Port covers.<br>Not made with DEHP or natural rubber latex. Stopcock volume: 0.32 mL. Length: 1 in. (2.54 cm.). | 456061     | 100           |







| Description                                                                                                                                                                                                                                                                                                                                                          | REF Number | Case Quantity |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| <b>ULTRAPORT™ ZerO and ULTRAPORT™ Luer-Activated Stopcocks</b>                                                                                                                                                                                                                                                                                                       |            |               |
|  <p>High-flow, four-way stopcock with Luer-activated ULTRAPORT ZerO Port, female Luer lock port and SPIN-LOCK® Connector. Not made with DEHP or natural rubber latex.<br/>Max pressure: 29 psi (2 bar). Valve port volume: 0.1 mL. Stopcock volume: 0.4 mL.</p>                     | 456064     | 100           |
|  <p>High-flow, four-way stopcock with two Luer-activated ULTRAPORT ZerO Ports and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Max pressure: 29 psi (2 bar). Valve port volume: 0.1 mL. Stopcock volume: 0.4 mL.</p>                                            | 456065     | 100           |
|  <p>High-flow, four-way stopcock with Luer-activated ULTRAPORT Port, female Luer lock port and SPIN-LOCK Connector. Max pressure: 29 psi (2 bar). Not made with DEHP or natural rubber latex. Max pressure: 29 psi (2 bar). Valve port volume: 0.1 mL. Stopcock volume: 0.5 mL.</p> | 456066     | 100           |
|  <p>High-flow, four-way stopcock with two Luer-activated ULTRAPORT Ports and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Max pressure: 29 psi (2 bar). Valve port volume: 0.1 mL. Stopcock volume: 0.4 mL.</p>                                                 | 456063     | 100           |
| <b>Manifolds</b>                                                                                                                                                                                                                                                                                                                                                     |            |               |
|  <p>Tri-connector with female Luer lock connector, two bonded CARESITE® Luer Access Devices and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex.<br/>Priming volume: 1.3 mL. Length: 2.5 in. (6.35 cm.).</p>                                                      | 470152     | 50            |
|  <p>Tri-connector with two bonded normally closed backcheck valves and SPIN-LOCK connector. Not made with DEHP or natural rubber latex. Priming volume: 1.14 mL. Length: 2.5 in. (6.35 cm.).</p>                                                                                  | 472026     | 50            |







# Anesthesia Extension Sets

| Description                                                                                                                                                                                                                                             | REF Number | Case Quantity |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| <b>Open Lumen Stopcocks</b>                                                                                                                                                                                                                             |            |               |
| Standard bore extension set with four-way open lumen stopcock and SPIN-LOCK® Connector. Not made with DEHP or natural rubber latex. Priming volume: 1.4 mL Length: 14 in. ( 35.5 cm.).                                                                  |            |               |
|                                                                                                                                                                       | 454320     | 100           |
| Standard bore extension set with four-way open lumen stopcock and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 1.8 mL Length: 23 in. ( 58 cm.).                                                                     |            |               |
|                                                                                                                                                                       | 456025     | 50            |
| Standard bore extension set with four-way open lumen stopcock and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 2.6 mL Length: 35 in. ( 89 cm.).                                                                     |            |               |
|                                                                                                                                                                       | 456030     | 50            |
| Standard bore extension set with four-way, high-flow open lumen stopcock and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 3.5 mL Length: 22 in. ( 56 cm.).                                                          |            |               |
|                                                                                                                                                                      | 473541     | 50            |
| Standard bore extension set with four-way, high-flow open lumen stopcock and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 5.9 mL Length: 37 in. ( 94 cm.).                                                          |            |               |
|                                                                                                                                                                     | 473543     | 50            |
| Standard bore extension set with three, four-way, high-flow open lumen stopcocks and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 6 mL Length: 47 in. ( 119 cm.).                                                   |            |               |
|                                                                                                                                                                     | 473054     | 100           |
| Standard bore extension set with four-way, high-flow open lumen stopcock, pinch clamp, CARESITE® Luer Access Device and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 6.1 mL Length: 40 in. (102 cm.).               |            |               |
|                                                                                                                                                                     | 490318     | 50            |
| Standard bore extension set with two removable, four-way, high-flow open lumen stopcocks, slide clamp, CARESITE Luer Access Device and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 5.3 mL Length: 33 in. (84 cm.). |            |               |
|                                                                                                                                                                     | 490418     | 50            |









 ULTRAPORT Zero Stopcock

 Needleless Valve (ULTRAPORT / ULTRAPORT Zero)







| Description                                                                                                                                                                                                                                                                           | REF Number | Case Quantity |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| <b>Needleless Stopcocks</b>                                                                                                                                                                                                                                                           |            |               |
| Standard bore extension set with four-way, high-flow ULTRAPORT™ Stopcock with Luer activated needleless port and SPIN-LOCK® Connector. Not made with DEHP or natural rubber latex. Priming volume: 1.5 mL Length: 10 in. (25 cm.).                                                    |            |               |
|  8"                                                                                                                                                                                                 | 456507     | 100           |
| Standard bore extension set with four-way, high-flow ULTRAPORT Stopcock with Luer activated needleless port, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 4.6 mL Length: 29 in. (74 cm.).                                         |            |               |
|  20" 6" 6"                                                                                                                                                                                          | 456503     | 50            |
| Standard bore extension set with four-way, high-flow ULTRAPORT Stopcock with Luer activated needleless port, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 6.9 mL Length: 43 in. (109 cm.).                                        |            |               |
|  35" 6" 6"                                                                                                                                                                                          | 456508     | 50            |
| Standard bore extension set with four-way, high-flow ULTRAPORT™ ZerO Stopcock with Luer activated needleless port, CARESITE® Luer Access Device, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 6.7 mL Length: 42 in. (107 cm.).    |            |               |
|  26" 8" 6"                                                                                                                                                                                         | 457505     | 100           |
| Standard bore extension set with four-way, high-flow ULTRAPORT Stopcock with Luer activated needleless port, CARESITE needleless connector, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 6.3 mL Length: 42 in. (107 cm.).         |            |               |
|  26" 8" 6"                                                                                                                                                                                        | 354222     | 100           |
| Standard bore extension set with four-way, high-flow ULTRAPORT ZerO Stopcock with Luer activated needleless port, two CARESITE Luer Access Devices, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 6.7 mL Length: 42 in. (107 cm.). |            |               |
|  12" 20" 6"                                                                                                                                                                                       | 490460     | 50            |
| Standard bore extension set with three four-way, high-flow ULTRAPORT Stopcocks with Luer activated needleless ports and SPIN-LOCK Connector. Not made with natural rubber latex. Priming volume: 8 mL Length: 50 in. (127 cm.).                                                       |            |               |
|  44"                                                                                                                                                                                              | 490139     | 50            |
| Standard bore extension set with five four-way, high-flow ULTRAPORT Stopcocks with Luer activated needleless ports, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 2.4 mL Length: 15 in. (38 cm.).                                  |            |               |
|                                                                                                                                                                                                   | 490295     | 24            |

| Description                                                                                                                                                                                                                                                                                        | REF Number | Case Quantity |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| <b>Extension Sets with SAFEPORT™</b>                                                                                                                                                                                                                                                               |            |               |
| Standard bore extension set with two port SAFEPORT Open Lumen Manifold, backcheck valve and SPIN-LOCK® Connector. Not made with natural rubber latex. Priming volume: 1.4 mL. Length: 7 in. (18 cm.).                                                                                              |            |               |
|                                                                                                                                                                                                                  | 456520     | 50            |
| Standard bore extension set with three port SAFEPORT Open Lumen Manifold with in-line backcheck valve and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 1.9mL. Length: 8 in. (20.3cm).                                                                          |            |               |
|                                                                                                                                                                                                                  | 456522     | 50            |
| Standard bore extension set with three port SAFEPORT Manifold with Luer activated needleless ports, backcheck valve and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 1.9 mL. Length: 8 in. (20.3 cm.).                                                         |            |               |
|                                                                                                                                                                                                                  | 456521     | 50            |
| Standard bore extension set with two port SAFEPORT Open Lumen Manifold, four-way, high-flow open lumen Stopcock and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 6.3 mL. Length: 39 in. (99 cm.).                                                              |            |               |
|                                                                                                                                                                                                                | 472036     | 50            |
| Standard bore extension set with two port SAFEPORT Open Lumen Manifold, backcheck valve, four-way, high-flow open lumen Stopcock and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 7 mL. Length: 41 in. (104 cm.).                                              |            |               |
|                                                                                                                                                                                                                | 472038     | 50            |
| Standard bore extension set with four-way, high-flow ULTRAPORT™ Stopcock with Luer activated needleless port, two port SAFEPORT Manifold with Luer activated needleless ports and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 8 mL. Length: 56 in. (142 cm.). |            |               |
|                                                                                                                                                                                                                | 456524     | 50            |
| Standard bore extension set with three port SAFEPORT Manifold with Luer activated needleless ports, CARESITE® Luer Access Device, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 9.2 mL. Length: 54 in. (137 cm.).                               |            |               |
|                                                                                                                                                                                                                | 490316     | 50            |






# Anesthesia IV Administration Sets

| Description                                                                                                                                                                                                                                                                                                                                                                                                        | REF Number | Case Quantity |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| <b>60 Drop Anesthesia IV Administration Sets</b>                                                                                                                                                                                                                                                                                                                                                                   |            |               |
| (60 drop/mL) Universal spike set with backcheck valve, roller clamp, two slide clamps, two CARESITE® Luer Access Devices, four-way, high-flow open lumen stopcock, two port SAFEPORT™ Open Lumen Manifold and SPIN-LOCK® Connector. Not made with DEHP or natural rubber latex. Priming volume: 21 mL. Length: 122 in. (310 cm.).                                                                                  | 352394     | 50            |
|                                                                                                                                                                                                                                                                                                                                  |            |               |
| (60 drop/mL) Universal spike set with backcheck valve, roller clamp, two CARESITE Luer Access Devices, four-way, high-flow ULTRAPORT™ ZerO Stopcock with Luer activated needleless port, slide clamp and SPIN-LOCK Connector. Not made with natural rubber latex. Priming volume: 26 mL. Length: 141 in. (358 cm.).                                                                                                | 490349     | 24            |
|                                                                                                                                                                                                                                                                                                                                  |            |               |
| (60 drop/mL) Universal spike set with backcheck valve, roller clamp, two CARESITE Luer Access Devices, two four-way, high-flow ULTRAPORT ZerO stopcocks with Luer activated needleless ports and SPIN-LOCK Connector. Not made with natural rubber latex. Priming volume: 20 mL. Length: 106 in. (269 cm.).                                                                                                        | 490392     | 24            |
|                                                                                                                                                                                                                                                                                                                                  |            |               |
| <b>15 Drop Anesthesia IV Administration Sets with Stopcocks</b>                                                                                                                                                                                                                                                                                                                                                    |            |               |
| (15 drop/mL) Universal spike set with backcheck valve, roller clamp, two slide clamps, three CARESITE Luer Access Devices, two four-way, high-flow open lumen Stopcocks and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 25 mL. Length: 139 in. (353 cm.).                                                                                                                     | 490535     | 24            |
|                                                                                                                                                                                                                                                                                                                                |            |               |
| (15 drop/mL) Universal spike set with backcheck valve, roller clamp, two CARESITE Luer Access Devices, four-way, high-flow ULTRAPORT™ Stopcock with Luer activated needleless port, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 23 mL. Length: 131 in. (333 cm.).                                                                                             | 354209     | 50            |
|                                                                                                                                                                                                                                                                                                                                |            |               |
| (15 drop/mL) Universal spike set with two backcheck valves, roller clamp, four CARESITE Luer Access Devices, four-way, high-flow ULTRAPORT Stopcock with Luer activated needleless port, SPIN-LOCK Connector and standard bore extension set with CARESITE Luer Access Device, slide clamp and SPIN-LOCK® Connector. Not made with DEHP or natural rubber latex. Priming volume: 24 mL. Length: 132 in. (335 cm.). | 490067     | 50            |
|                                                                                                                                                                                                                                                                                                                                |            |               |
| (15 drop/mL) Universal spike set with backcheck valve, roller clamp, slide clamp, two CARESITE Luer Access Devices, two four-way, high-flow ULTRAPORT Stopcocks with Luer activated needleless port and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 24 mL. Length: 130 in. (330 cm.).                                                                                         | 354210     | 50            |
|                                                                                                                                                                                                                                                                                                                                |            |               |
| (15 drop/mL) Universal spike set with backcheck valve, roller clamp, two CARESITE Luer Access Devices, four-way, high-flow ULTRAPORT ZerO Stopcock with Luer activated needleless port, slide clamp and SPIN-LOCK Connector. Not made with natural rubber latex. Priming volume: 23 mL. Length: 131 in. (333 cm.).                                                                                                 | 457506     | 50            |
|                                                                                                                                                                                                                                                                                                                                |            |               |

 ULTRAPORT ZerO Stopcock
  Needleless Valve (ULTRAPORT / ULTRAPORT ZerO)

| Description                                                                                                                                                                                                                                                                                                                                                                                                                | REF Number | Case Quantity |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| <b>15 Drop Anesthesia IV Administration Sets with Stopcocks (continued)</b>                                                                                                                                                                                                                                                                                                                                                |            |               |
| (15 drop/mL) Universal spike set with backcheck valve, roller clamp, four CARESITE® Luer Access Devices, four-way, high-flow ULTRAPORT™ ZerO Stopcock with Luer activated needleless port, two slide clamps, and SPIN-LOCK® Connector. Not made with DEHP or natural rubber latex. Priming volume: 23 mL. Length: 120 in. (304 cm.).                                                                                       | 490427     | 24            |
|                                                                                                                                                                                                                                                                                                                                          |            |               |
| (15 drop/mL) Universal spike set with backcheck valve, roller clamp, two CARESITE Luer Access Devices, two four-way, high-flow ULTRAPORT ZerO Stopcocks with Luer activated needleless ports, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 23 mL. Length: 130 in. (330 cm.).                                                                                           | 457504     | 50            |
|                                                                                                                                                                                                                                                                                                                                          |            |               |
| (15 drop/mL) Universal spike set with backcheck valve, roller clamp, two CARESITE Luer Access Devices, three four-way, high-flow ULTRAPORT ZerO Stopcocks with Luer activated needleless ports, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 22 mL. Length: 117 in. (297 cm.).                                                                                         | 490468     | 24            |
|                                                                                                                                                                                                                                                                                                                                          |            |               |
| <b>15 Drop Anesthesia IV Administration Sets with SAFEPORT™ Manifolds</b>                                                                                                                                                                                                                                                                                                                                                  |            |               |
| (15 drop/mL) Universal spike set with backcheck valve, roller clamp, three SAFELINE® Needleless Connectors, three slide clamps, four-way, high-flow open lumen stopcock, two port SAFEPORT Open Lumen Manifold and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 24 mL. Length: 134 in. (340 cm.).                                                                                      | NF1270A    | 50            |
|                                                                                                                                                                                                                                                                                                                                        |            |               |
| (15 drop/mL) Universal spike set with backcheck valve, roller clamp, two CARESITE Luer Access Devices, two port SAFEPORT Open Lumen Manifold, three slide clamps, removable 6 in. extension set and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 23 mL. Length: 131 in. (333 cm.).                                                                                                     | 354208     | 50            |
|                                                                                                                                                                                                                                                                                                                                        |            |               |
| (15 drop/mL) Universal spike set with backcheck valve, roller clamp, three CARESITE Luer Access Devices, three slide clamps, four-way, high-flow open lumen stopcock, two port SAFEPORT Open Lumen Manifold and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 24 mL. Length: 137 in. (348 cm.).                                                                                         | 354211     | 50            |
|                                                                                                                                                                                                                                                                                                                                        |            |               |
| (15 drop/mL) Universal spike set with backcheck valve, roller clamp, three CARESITE Luer Access Devices, two slide clamps, high-flow ULTRAPORT ZerO stopcock with Luer activated needleless port, two port SAFEPORT Manifold with Luer activated needleless ports, removable 6 in. extension tubing and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 24 mL. Length: 136 in. (345 cm.). | 457515     | 50            |
|                                                                                                                                                                                                                                                                                                                                        |            |               |

# Blood Administration Sets







| Description                                                                                                                                                                                                                                                                                              | REF Number | Case Quantity |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| Blood Administration Sets                                                                                                                                                                                                                                                                                |            |               |
| (10 drop/mL) Drip chamber with 170 micron blood filter, roller clamp, slide clamp and SPIN-LOCK® Connector, and bag hanger. Not made with natural rubber latex. Priming volume: 23 mL. Length: 40 in. (101 cm.).                                                                                         | 470032     | 50            |
|                                                                                                                                                                                                                        |            |               |
| (10 drop/mL) Y-type blood set with two non-vented spikes, 170 micron filter, slide clamp, three roller clamps, CARESITE® Luer Access Device and SPIN-LOCK Connector. Not made with natural rubber latex. Priming volume: 31 mL. Length: 110 in. (279 cm.).                                               | 354217     | 50            |
|                                                                                                                                                                                                                        |            |               |
| (10 drop/mL) Y-type blood set with two non-vented spikes, 170 micron filter, hand pump, slide clamp, three roller clamps, CARESITE Luer Access Device and SPIN-LOCK Connector. Not made with natural rubber latex. Priming volume: 75 mL. Length: 89 in. (226 cm.).                                      | 354218     | 50            |
|                                                                                                                                                                                                                      |            |               |
| (10 drop/mL) Y-type blood set with two non-vented spikes, 170 micron filter, hand pump, three roller clamps, three CARESITE Luer Access Devices, four-way high-flow open lumen stopcock and SPIN-LOCK Connector. Not made with natural rubber latex. Priming volume: 74 mL. Length: 105 in. (267 cm.).   | 490425     | 24            |
|                                                                                                                                                                                                                      |            |               |
| (10 drop/mL) Y-type blood set with two non-vented spikes, 170 micron filter, slide clamp, three roller clamps, three CARESITE Luer Access Devices, four-way high-flow open lumen stopcock and SPIN-LOCK Connector. Not made with natural rubber latex. Priming volume: 38 mL. Length: 135 in. (343 cm.). | 490293     | 24            |
|                                                                                                                                                                                                                      |            |               |

## Basic IV Administration Sets

| Description                                                                                                                                                                                                                             | REF Number | Case Quantity |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| <b>Basic IV Administration Sets</b>                                                                                                                                                                                                     |            |               |
| (15 drops/mL) Universal spike, roller clamp, slide clamp and SPIN-LOCK® Connector. Not made with DEHP or natural rubber latex. Priming Volume: 13 mL. Length: 85 in. (215.9 cm.).                                                       | 352601     | 50            |
| (60 drops/mL) Universal spike, roller clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming Volume: 13 mL. Length: 67 in. (170.2 cm.).                                                                     | V1425-15   | 50            |
| Secondary Administration Set (15 drops/mL) Universal spike, roller clamp and SPIN-LOCK Connector with 10.5" plastic bag hanger. Not made with DEHP or natural rubber latex. Priming Volume: 9.9 mL. Length: 40 in. (101.6 cm.).         | V1921      | 50            |
| <b>Standard Bore Extension Sets</b>                                                                                                                                                                                                     |            |               |
| Female Luer connector, distal SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 3 mL. Length 21 in. (52.5 cm.).                                                                                          | V5406      | 50            |
| Female Luer connector, large bore tubing, roller clamp and distal SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 4.3 mL. Length: 31 in. (78.5 cm.). Priming volume: 74 mL. Length: 105 in. (267 cm.). | V5484      | 50            |
| Female Luer connector, slide clamp and distal SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 1 mL. Length: 7 in. (18 cm.).                                                                            | 471975     | 100           |
| Female Luer connector, slide clamp and distal SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 4.3 mL. Length: 30 in. (76 cm.).                                                                         | 473012     | 100           |


# Extension Sets

| Description                                                                                                                                                                                                                                                                            | REF Number | Case Quantity |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| <b>Small Bore Extension Sets</b>                                                                                                                                                                                                                                                       |            |               |
| Female Luer connector, slide clamp and SPIN-LOCK® Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.2 mL. Length 6 in. (15.2 cm.).                                                                                                                              | 471960     | 100           |
| <p><i>For use with power injectors - maximum 300 psi and maximum flow rate 10 mL/sec</i></p> <p>Female Luer connector, slide clamp and SPIN-LOCK connector. Not made with DEHP or natural rubber latex. Priming volume: 0.3 mL. Length: 7 in. (18 cm.). <b>High Pressure Rated</b></p> | 470064     | 50            |
| Female Luer connector, slide clamp and distal SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.44 mL. Length 13 in. (31.0 cm.).                                                                                                                      | 473105     | 100           |
| Female Luer connector, slide clamp and distal SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.2 mL. Length: 6 in. (15.2cm.).                                                                                                                        | 471991     | 100           |
| Female Luer connector, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 2 mL. Length: 61 in. (155 cm.).                                                                                                                                | 471973     | 100           |
| Y-type with two female Luer lock connectors, two slide clamps and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.46 mL. Length: 6 in. (15 cm.).                                                                                                    | 471957     | 100           |
| Y-type with two female Luer lock connectors, two pinch clamps and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.5 mL. Length: 8 in. (20 cm.).                                                                                                     | 471959     | 100           |
| Triple-leg set with three female Luer connectors, three slide clamps and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.5 mL. Length: 9 in. (23 cm.).                                                                                              | 473041     | 100           |

| Description                                                                                                                                                                                                                                                                     | REF Number | Case Quantity |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| <b>Syringe Pump PE-Lined Fluid Path Microbore Extension Sets</b>                                                                                                                                                                                                                |            |               |
| Female Luer and SPIN-LOCK® Connector. Microbore extension set, 0.02 in. ID. tubing. Not made with DEHP or natural rubber latex. Priming volume: 0.3 mL. Length: 36 in. (91 cm.).                                                                                                | V6220      | 50            |
|                                                                                                                                                                                               |            |               |
| Female Luer, slide clamp and SPIN-LOCK Connector. Microbore extension set, 0.03 in. ID. tubing. Use with higher flow rates and viscous fluids. Not made with DEHP or natural rubber latex. Priming volume: 0.6 mL. Length: 36 in. (91 cm.).                                     | V6203      | 50            |
|                                                                                                                                                                                               |            |               |
| Female Luer, slide clamp and SPIN-LOCK Connector. Microbore extension set, 0.02 in. ID. tubing with fluid path that does not contain PVC. Not made with DEHP or natural rubber latex. Priming volume: 0.5 mL. Length: 60 in. (152 cm.).                                         | V6222      | 50            |
|                                                                                                                                                                                               |            |               |
| Female Luer, slide clamp and SPIN-LOCK Connector. Microbore extension set, 0.03 in. ID. with fluid path that does not contain PVC. Use with higher flow rates and viscous fluids. Not made with DEHP or natural rubber latex. Priming volume: 0.8 mL. Length: 60 in. (152 cm.). | V6223      | 50            |
|                                                                                                                                                                                              |            |               |
| Female Luer, slide clamp and SPIN-LOCK Connector. Microbore extension set, 0.2 micron filter, 0.03 in. ID. tubing with fluid path that does not contain PVC. Not made with DEHP or natural rubber latex. Priming volume: 1.3 mL. Length: 60 in. (152 cm.).                      | V6215      | 50            |
|                                                                                                                                                                                             |            |               |
| <b>Syringe Pump Microbore Extension Sets with Proximal and Distal Male Luer Locks</b>                                                                                                                                                                                           |            |               |
| Proximal and SPIN-LOCK Connector. Microbore extension set, 0.04 in. ID., slide clamp. Not made with natural rubber latex. Priming volume: 0.6 mL. Length: 31 in. (79 cm.).                                                                                                      | V5450      | 50            |
|                                                                                                                                                                                             |            |               |
| 1.2 micron air-eliminating filter with female Luer lock connector, male SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming Volume: 4.3 mL. Length: 10 in. (25 cm.).                                                                                       | 473994     | 50            |
|                                                                                                                                                                                             |            |               |

| Description             | REF<br>Number | Case<br>Quantity |
|-------------------------|---------------|------------------|
| Filtered Extension Sets |               |                  |


5.0 micron air-eliminating in-line filter with female Luer connector, SPIN-LOCK® connector. Not made with DEHP or natural rubber latex. Priming volume: 3.6 mL. Length: 10 in. (25.4 cm.).



474002

50


0.2 micron air eliminating filter with pinch clamp, female Luer connector, male SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 4 mL. Length: 9 in. (23 cm.).



470073

50


0.2 micron low volume SUPOR® Filter, slide clamp, female Luer connector, male SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 1.2 mL. Length 8 in. (20 cm.).



352205

50





1.2 micron low volume SUPOR Filter, slide clamp, female Luer connector, male SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.9 mL. Length 8 in. (20 cm.).









352212

50

# CARESAFE™ IV Administration Sets and Extension Sets





| Description                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | REF Number | Case Quantity |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| CARESAFE IV Administration Sets with AirStop                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |            |               |
| <p>(20 drops/mL.) Drip chamber with AirStop, backcheck valve, CARESITE Injection Sites 87 in. and 6 in. above distal end, roller clamp, 0.2 micron air eliminating filter, slide clamp and SPIN-LOCK Connector. Components are not made with PVC, DEHP or natural rubber latex. Priming volume: 28 mL. Length: 110 in. (279 cm.).</p>                                                                                                                                                   | 354306     | 50            |
| <p>(20 drops/mL.) Drip chamber with AirStop, backcheck valve, two CARESITE Injection Sites, roller clamp, two removable high-flow, four-way ULTRAPORT™ Stopcocks with Luer activated needless ports, female Luer lock port and SPIN-LOCK Connector, removable 6 in. (15.2 cm.). extension tubing with Luer lock connection, slide clamp and SPIN-LOCK Connector. Components are not made with PVC, DEHP or natural rubber latex. Priming volume: 30 mL. Length: 134 in. (340 cm.).</p>  | 354308     | 50            |
| CARESAFE IV Administration Sets                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |            |               |
| <p>(20 drops/mL.) Drip chamber with 15µm filter, slide clamp, backcheck valve, CARESITE Luer Access Device Injection sites 87 in. and 12 in. above distal end, roller clamp and SPIN-LOCK Connector. Components are not made with PVC, DEHP or natural rubber latex. Priming volume: 24 mL. Length: 104 in. (264 cm.).</p>                                                                                                                                                             | 354303     | 50            |
| <p>(20 drops/mL.) Drip chamber with 15 µm filter, backcheck valve, three CARESITE Luer Access Devices, roller clamp, slide clamp and SPIN-LOCK Connector. Components are not made with PVC, DEHP or natural rubber latex. Priming volume: 26 mL. Length: 114 in. (290 cm.).</p>                                                                                                                                                                                                       | 354305     | 50            |

| Description                                                                                                                                                                                                                                                                                                                                                                                                                                                  | REF Number | Case Quantity |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| <b>CARESAFE IV Administration Sets (continued)</b>                                                                                                                                                                                                                                                                                                                                                                                                           |            |               |
| (20 drops/mL.) Drip chamber with 15µm filter, backcheck valve, CARESITE® Injection Sites 87 in. and 6 in. above distal end, roller clamp, 0.2 micron filter, slide clamp and SPIN-LOCK® Connector. Components are not made with PVC, DEHP or natural rubber latex. Priming volume: 28 mL. Length: 110 in. (279 cm.).                                                                                                                                         | 354307     | 50            |
|                                                                                                                                                                                                                                                                                                                                                                            |            |               |
| (20 drops/mL.) Drip chamber with 15µm filter, backcheck valve, two CARESITE Injection Sites, roller clamp, two removable high-flow, four-way ULTRAPORT™ Stopcocks with Luer activated needless ports and SPIN-LOCK Connector, removable 6 in. (15.2 cm.). extension tubing with Luer lock connection, slide clamp and SPIN-LOCK Connector. Components are not made with PVC, DEHP or natural rubber latex. Priming volume: 30 mL. Length: 134 in. (340 cm.). | 354310     | 50            |
|                                                                                                                                                                                                                                                                                                                                                                            |            |               |
| (20 drops/mL.) Drip chamber with 15µm filter, roller clamp and SPIN-LOCK Connector with plastic bag hanger. Components are not made with PVC, DEHP or natural rubber latex. Priming Volume: 12 mL. Length: 41 in. (104 cm.).                                                                                                                                                                                                                                 | V1931      | 50            |
|                                                                                                                                                                                                                                                                                                                                                                           |            |               |
| <b>CARESAFE Extension Sets</b>                                                                                                                                                                                                                                                                                                                                                                                                                               |            |               |
| Standard bore extension set with two CARESITE Luer Access Device injection sites 28 in. and 6 in. above distal end, slide clamp and SPIN-LOCK Connector. Components are not made with PVC, DEHP or natural rubber latex. Priming volume: 6.3 mL. Length: 35 in. (89 cm.).                                                                                                                                                                                    | 354320     | 50            |
|                                                                                                                                                                                                                                                                                                                                                                          |            |               |
| Standard bore extension set with bonded CARESITE Luer Access Device, slide clamp and SPIN-LOCK Connector. Components are not made with PVC, DEHP or natural rubber latex. Priming volume: 1.3 mL. Length: 8 in. (20 cm.).                                                                                                                                                                                                                                    | 471858     | 50            |
|                                                                                                                                                                                                                                                                                                                                                                          |            |               |
| Standard bore extension set with removable CARESITE Valve, slide clamp and SPIN-LOCK Connector. Components are not made with PVC, DEHP or natural rubber latex. Priming volume: 1.3 mL. Length: 8 in. (20 cm.).                                                                                                                                                                                                                                              | 471859     | 50            |
|                                                                                                                                                                                                                                                                                                                                                                          |            |               |





# STEADYCare™ Extension Sets



| Description                                                                                                                                                                                                                                        | REF Number | Case Quantity |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| STEADYCare Extension Sets                                                                                                                                                                                                                          |            |               |
| Small bore tubing, slide clamp, bonded CARESITE® Luer Access Device, SPIN-LOCK® Connector, STEADYCare Extension Set with Wedge Catheter Stabilizer and slide clamp, 8 in. Not made with DEHP or natural rubber latex.                              | 470185     | 100           |
| Standard bore tubing, slide clamp, bonded CARESITE Luer Access Device, STEADYCare Extension Set with Wedge Catheter Stabilizer and SPIN-LOCK Connector, 7 in. Not made with DEHP or natural rubber latex.                                          | 470186     | 100           |
| Standard bore tubing, slide clamp, bonded CARESITE Luer Access Device, STEADYCare Extension Set with Wedge Catheter Stabilizer and SPIN-LOCK Connector, 8 in. Not made with DEHP or natural rubber latex.                                          | 470188     | 100           |
| Small bore tubing, slide clamp, female Luer lock connector, STEADYCare Extension Set with Wedge Catheter Stabilizer and SPIN-LOCK Connector, 7 in. Not made with DEHP or natural rubber latex.                                                     | 470193     | 50            |
| Small bore tubing. STEADYCare Extension Set with Wedge™ Catheter Stabilizer, removable CARESITE Luer Access Device, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.6 mL. Length 9 in (23 cm.). | 477033     | 50            |

# Rate Flow® Regulator II Administration Sets









| Description                                                                                                                                                                                                                                      | REF Number | Case Quantity |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| Rate Flow Regulator II Administration Sets                                                                                                                                                                                                       |            |               |
| Rate Flow Extension Set with CARESITE® Luer Access Device and SPIN-LOCK® Connector. Not made with DEHP or natural rubber latex. Priming Volume: 2.8 mL. Length: 17 in.                                                                           |            |               |
|                                                                                                                                                                | 375256     | 50            |
| Rate Flow IV Set with 15 drop/mL drip chamber with 15 micron filter, CARESITE Luer Access Device, pinch clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming Volume: 15 mL. Length: 82 in.                         |            |               |
|                                                                                                                                                                | 375257     | 50            |
| Rate Flow IV Set with 15 drop/mL drip chamber with 15 micron filter, backcheck valve, two CARESITE® Luer Access Devices, pinch clamp and SPIN-LOCK® Connector. Not made with DEHP or natural rubber latex. Priming Volume: 17 mL. Length: 92 in. |            |               |
|                                                                                                                                                              | 375258     | 50            |
| Rate Flow IV Set with 60 drop/mL drip chamber with 15 micron filter, injection site, pinch clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming Volume: 15 mL. Length: 82 in.                                      |            |               |
|                                                                                                                                                              | 375252     | 50            |
| Rate Flow IV Set with 10 drop/mL drip chamber, backcheck valve, two CARESITE Luer Access Devices, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming Volume: 21 mL. Length: 108 in.                        |            |               |
|                                                                                                                                                              | 375260     | 50            |

# Irrigation/Urology Sets





| Description                                                                                                                                                                                                                                                      | REF Number | Case Quantity |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| Irrigation/Urology Sets                                                                                                                                                                                                                                          |            |               |
| Two-bag irrigation set, two non-vented spikes, 0.281 in. ID large bore tubing, two on/off clamps, drip chamber, roller clamp and 4½ in. distal flexible connector. Not made with natural rubber latex. Length: 94 in. (238.8 cm.). For use with irrigation bags. |            |               |
|                                                                                                                                                                                | 313003     | 10            |
| Continuous overhead irrigation set, non-vented spike, 0.188 in. ID tubing, drip chamber, roller clamp and 4½ in. distal flexible connector. Not made with natural rubber latex. Length: 82 in. (208.3 cm.). For use with irrigation bags.                        |            |               |
|                                                                                                                                                                                | V4608-20   | 20            |
| Four-bag irrigation set, four non-vented spikes, four on/off clamps, drip chamber, roller clamp, 0.281 in. ID large bore tubing, 4½ in. distal flexible connector. Not made with natural rubber latex. Length: 96 in. (243.8 cm.). For use with irrigation bags. |            |               |
|                                                                                                                                                                               | 313004     | 10            |
| UniSpike Connector, 0.188 in. ID tubing, drip chamber, roller clamp and distal flexible connector with a 4½ in. distal flexible connector. Not made with natural rubber latex. Length: 79 in. (200.7 cm.). For use with plastic irrigation containers.           |            |               |
|                                                                                                                                                                              | 313005     | 10            |

| Description                                                                                                                                                                                                                                                                                           | REF Number | Case Quantity |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| <b>Irrigation/Urology Sets</b>                                                                                                                                                                                                                                                                        |            |               |
| <p>SuperSet™ Two-Bottle Irrigation Set, Two UniSpike™ Connectors, 0.281 in. ID large bore tubing, two on/off clamps, drip chamber, roller clamp and 4½ in. distal flexible connector. Not made with natural rubber latex. Length: 94 in. (238.8 cm.). For use with plastic irrigation containers.</p> |            |               |
|                                                                                                                                                                                                                     | 313000     | 10            |
| <p>T.U.R./Cysto set Non-Vented Spike, 0.281 in. ID large bore tubing, on/off clamp and 4½ in. distal flexible connector. Not made with natural rubber latex. Length: 80 in. (203.20cm.). For use with irrigation bags.</p>                                                                            |            |               |
|                                                                                                                                                                                                                     | 313002     | 10            |
| <p>Cysto Post-Operative Irrigation Set, UniSpike Connector, 0.188 in. ID tubing, on/off clamp and 4½ in. distal flexible connector. Not made with natural rubber latex. Length: 80 in. (203.2 cm.). For use with plastic irrigation containers.</p>                                                   |            |               |
|                                                                                                                                                                                                                   | V4500      | 20            |
| <p>PIC™ Irrigation Container Spike Adapter — UniSpike Adapter for B. Braun. Permits the connection of a non-vented irrigation bag spike to a plastic irrigation container.</p>                                                                                                                        |            |               |
|                                                                                                                                                                                                                   | N2150      | 50            |

# Measured Volume Solution Burette Sets and Special Application Accessories

| Description                                                                                                                                                                                                                                                                                         | REF Number | Case Quantity |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| <b>Measured Volume Solution Burette Sets</b>                                                                                                                                                                                                                                                        |            |               |
| Non-vented spike, roller clamp, 150 mL. burette chamber without automatic shutoff with CARESITE® Luer Access Device, slide clamp and spike adapter. Not made with natural rubber latex. Priming volume: 3.0 mL. Length: 20 in. (51.0 cm.).                                                          | 470118     | 20            |
|                                                                                                                                                                                                                   |            |               |
| Non-vented spike, roller clamp, 150 mL. burette chamber, slide clamp and spike adapter. Components are not made with natural rubber latex. Not made with natural rubber latex. Priming Volume: 3 mL. Length: 20 in. (51.0 cm.).                                                                     | 375112     | 20            |
|                                                                                                                                                                                                                   |            |               |
| (60 drops/mL.) Non-vented spike, roller clamp, 150 mL. burette chamber with CARESITE Luer Access Device, roller clamp, slide clamp, CARESITE Injection Site 6 in. above distal end and SPIN-LOCK® Connector. Not made with natural rubber latex. Priming volume: 16.0 mL. Length: 89 in. (226 cm.). | 470116     | 20            |
|                                                                                                                                                                                                                   |            |               |
| <b>Special Application Accessories</b>                                                                                                                                                                                                                                                              |            |               |
| 0.2 micron SUPOR Low Volume Air-Eliminating IV Filter with female Luer lock and SPIN-LOCK connector. Not made with natural rubber latex. Priming Volume: 0.8 mL.                                                                                                                                    | 473036     | 50            |
|                                                                                                                                                                                                                  |            |               |
| Spike adapter. Not made with natural rubber latex.                                                                                                                                                                                                                                                  | N2050      | 24            |
|                                                                                                                                                                                                                 |            |               |
| Vented spike adapter. Not made with DEHP or natural rubber latex.                                                                                                                                                                                                                                   | N2061      | 50            |
|                                                                                                                                                                                                                 |            |               |
| Vented spike adapter—permits connection of non-vented IV spike to IV bottle requiring a vent. Not made with natural rubber latex.                                                                                                                                                                   | SA2000     | 20            |
|                                                                                                                                                                                                                 |            |               |
| Dual male adapter—male Luer lock connector and distal SPIN-LOCK Connector. Connects two female Luer connections. Not made with DEHP or natural rubber latex. Priming volume: 0.1 mL.                                                                                                                | 456080     | 100           |
|                                                                                                                                                                                                                  |            |               |
| Normally closed check valve/antisiphon valve—requires syringe pressure to open, Luer lock. Not made with DEHP or natural rubber latex. Priming volume: 0.12 mL.                                                                                                                                     | 415062     | 100           |
|                                                                                                                                                                                                                  |            |               |


# Injection Caps and Replacement Caps

| Description                                                                                                                                                                                                                                                                            | REF Number                 | Case Quantity      |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------|--------------------|
| <b>Injection Caps</b>                                                                                                                                                                                                                                                                  |                            |                    |
| Intermittent injection cap-opaque yellow body, proximal self-sealing injection port, distal male Luer lock connector with recessed Luer taper. Not made with DEHP or natural rubber latex. Priming volume: 0.16 mL. Length: 0.875 in. ( 2.2 cm.).                                      | 654237                     | 1000               |
|                                                                                                                                                                                                       |                            |                    |
| Intermittent injection cap - proximal injection port, clear body, distal male Luer lock connector. Not made with DEHP or natural rubber latex. Priming volume: 0.2 mL. Length: 1.125 in. ( 2.9 cm.).                                                                                   | 418020                     | 100                |
|                                                                                                                                                                                                       |                            |                    |
| Intermittent injection cap-proximal injection port, distal Luer slip connector. Not made with DEHP or natural rubber latex. Priming volume: 0.2 mL. Length: 1.125 in. ( 2.9 cm.).                                                                                                      | 418025                     | 100                |
|                                                                                                                                                                                                       |                            |                    |
| Intermittent injection cap-proximal injection port, distal male Luer lock connector. Not made with DEHP or natural rubber latex. Priming volume: 0.1 mL. Length: 0.75 in. ( 1.9 cm.).                                                                                                  | 418030                     | 100                |
|                                                                                                                                                                                                       |                            |                    |
| <b>Replacement Caps</b>                                                                                                                                                                                                                                                                |                            |                    |
| Blue replacement cap-male Luer lock. Not made with DEHP or natural rubber latex. Blue replacement cap-bulk pack of 1,000 in dispenser boxes of 100 each. Not made with DEHP or natural rubber latex. White replacement cap-male Luer lock. Not made with DEHP or natural rubber latex. | 418017<br>418027<br>418015 | 100<br>1000<br>100 |
|                                                                                                                                                                                                     |                            |                    |
| Replacement caps - blue male Luer lock cap and white female Luer lock cap. Not made with DEHP or natural rubber latex.                                                                                                                                                                 | BW1000                     | 100                |
|                                                                                                                                                                                                     |                            |                    |
| White replacement cap - female Luer lock. Not made with DEHP or natural rubber latex.                                                                                                                                                                                                  | 474900                     | 100                |
|                                                                                                                                                                                                     |                            |                    |
| RED CAP™ Luer Cap-Dual function red cap with male and female end. Bulk pack of 1,000 in dispenser. Not made with DEHP or natural rubber latex.                                                                                                                                         | R2000B                     | 1000               |
|                                                                                                                                                                                                     |                            |                    |
| Blue cap Luer cap-Dual function blue cap with male and female end. Bulk pack of 1,000 in dispenser boxes of 100 each. Not made with DEHP or natural rubber latex.                                                                                                                      | B2000B                     | 1000               |
|                                                                                                                                                                                                     |                            |                    |

# CARESITE® Needleless Connectors

| Description                   | REF Number | Case Quantity |
|-------------------------------|------------|---------------|
| CARESITE Needleless Connector |            |               |

CARESITE Needleless Connector for aspiration, injection or gravity flow of fluid upon insertion of a male Luer fitting. Priming volume: 0.22 mL. Power injection rated to 400 psi and 15 mL/sec. Not made with DEHP or natural rubber latex.




415122

200

## CARESITE Small Bore Extension Sets

For use with power injectors - maximum 300 psi and maximum flow rate 10 mL/sec


Small bore extension set with bonded CARESITE Needleless Connector, slide clamp and SPIN-LOCK® Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.5 mL. Length: 8 in. (20 cm.). Power injection rated to 300 psi and 10 mL/sec. **High Pressure Rated**



470100

100


Small bore extension set with removable CARESITE Needleless Connector, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.5 mL. Length: 8 in. (20 cm.).



470101

100


Small bore extension set with bonded CARESITE Needleless Connector, slide clamp and male Luer slip connector. Not made with DEHP or natural rubber latex. Priming volume: 0.5 mL. Length: 8 in. (20 cm.).



470102

100


Small bore extension set with removable CARESITE Needleless Connector, slide clamp and male Luer slip connector. Not made with DEHP or natural rubber latex. Priming volume: 0.5 mL. Length: 8 in. (20 cm.).



470103





100

Small bore extension set with bonded CARESITE Needleless Connector, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.7 mL. Length: 14 in. (36 cm.).

















470104



100

| Description                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | REF Number | Case Quantity |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| <b>CARESITE Small Bore Extension Sets (continued)</b>                                                                                                                                                                                                                                                                                                                                                                                                                                                     |            |               |
| <p>Small bore extension set with removable CARESITE Needleless Connector, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.7 mL. Length: 14 in. (36 cm.).</p>                                                                                                                                                                                                          | 470105     | 100           |
| <p>Small bore T-port extension set with CARESITE Needleless Connector, SAFELINE® Injection Site on T-port, slide clamp and PE fluid path tubing. Not made with DEHP or natural rubber latex. Priming volume: 0.27 mL. Length: 7 in. (18 cm.).</p>                                                                                                                                                                        | 470128     | 50            |
| <p>Small bore Y-extension set with one removable CARESITE Needleless Connector, two slide clamps, SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.9 mL. Length: 7 in. (18 cm.).</p>                                                                                                                                                                                                 | 470107     | 50            |
| <p><i>For use with power injectors - maximum 300 psi and maximum flow rate 10 mL/sec</i></p> <p>Small bore high pressure rated Y-extension set with two bonded CARESITE Needleless Connectors, two slide clamps and SPIN-LOCK Connector. Set rated to 300 psi and max flow rate 10 mL/sec. Not made with DEHP or natural rubber latex. Priming volume: 1.0 mL. Length: 8 in. (20 cm.). <b>High Pressure Rated</b></p>  | 354250     | 50            |






# CARESITE® Needleless Connectors

| Description                                                                                                                                                                                                                                                                                                                                                                                      | REF Number | Case Quantity |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| <b>CARESITE Standard Bore Extension Sets</b>                                                                                                                                                                                                                                                                                                                                                     |            |               |
| <p>Standard bore extension set with bonded CARESITE Needleless Connector, slide clamp and SPIN-LOCK® Connector. Not made with DEHP or natural rubber latex.</p> <p>Priming volume: 0.9 mL. Length: 6 in. (15 cm.).</p>                                                                                                                                                                           | 470108     | 100           |
|                                                                                                                                                                                                                                                                                                                |            |               |
| <p>Standard bore extension set with removable CARESITE Needleless Connector, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex.</p> <p>Priming volume: 0.9 mL. Length: 6 in. (15 cm.).</p>                                                                                                                                                                         | 470109     | 100           |
|                                                                                                                                                                                                                                                                                                                |            |               |
| <p><i>For use with power injectors - maximum 300 psi and maximum flow rate 10 mL/sec</i></p> <p>Standard bore extension set with bonded CARESITE Needleless Connector, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.5 mL. Length: 8 in. (20 cm.). Power injection rated to 300 psi and 10 mL/sec. <b>High Pressure Rated</b></p>           | 470124     | 100           |
|                                                                                                                                                                                                                                                                                                                |            |               |
| <p>Standard bore extension set with slide clamp, CARESITE Needleless Connector injection site and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 1.5 mL. Length: 8 in. (20 cm.).</p>                                                                                                                                                                           | 354219     | 50            |
|                                                                                                                                                                                                                                                                                                               |            |               |
| <p>Standard bore extension set with two CARESITE Needleless Connector injection sites 28 in. and 6 in. above distal end, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 5.5 mL. Length: 35 in. (89 cm.).</p>                                                                                                                                   | 354220     | 50            |
|                                                                                                                                                                                                                                                                                                              |            |               |
| <p><i>For use with power injectors - maximum 300 psi and maximum flow rate 10 mL/sec</i></p> <p>Standard bore extension set. Female Luer Connector, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.4 mL. Length: 8 in. (20 cm.). Power injection rated to 300 psi and max. flow rate 10 mL/sec. <b>High Pressure Rated</b></p>               | 480223     | 50            |
|                                                                                                                                                                                                                                                                                                              |            |               |
| <p><i>For use with power injectors - maximum 300 psi and maximum flow rate 10 mL/sec</i></p> <p>Standard bore extension set with removable CARESITE Needleless Connector, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.5 mL. Length: 8 in. (20 cm.). Set rated to 300 psi and max. flow rate 10 mL/sec. <b>High Pressure Rated</b></p>     | 470159     | 100           |
|                                                                                                                                                                                                                                                                                                              |            |               |
| <p><i>For use with power injectors - maximum 300 psi and maximum flow rate 10 mL/sec</i></p> <p>Standard bore extension set with non-removable CARESITE Needleless Connector, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.5 mL. Length: 8 in. (20 cm.). Set rated to 300 psi and max. flow rate 10 mL/sec. <b>High Pressure Rated</b></p> | 470170     | 50            |
|                                                                                                                                                                                                                                                                                                              |            |               |

| Description                                                                                                                                                                                                                                                                                                  | REF Number | Case Quantity |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| <b>Filterflow® Filtered SUPOR® Extension Sets with CARESITE Needleless Connector Injection Sites</b>                                                                                                                                                                                                         |            |               |
| 0.2 micron filter, CARESITE Needleless Connector injection site 6 in. above distal end, on/off clamp and SPIN-LOCK® Connector. Not made with DEHP or natural rubber latex. Priming volume: 5 mL. Length: 16 in. (40.6 cm.).                                                                                  | 354221     | 50            |
|                                                                                                                                                                                                                            |            |               |
| 1.2 micron filter, CARESITE Needleless Connector injection site 6 in. above distal end, on/off clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 5 mL. Length: 16 in. (40.6 cm.).                                                                                   | 470117     | 50            |
|                                                                                                                                                                                                                            |            |               |
| <b>Filterflow Filtered Administration Sets with CARESITE Needleless Connector Injection Sites</b>                                                                                                                                                                                                            |            |               |
| (15 drops/mL.) Universal spike, backcheck valve, CARESITE Needleless Connector injection sites 87 in. and 6 in. above distal end, roller clamp, 0.2 micron filter, two slide clamps and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 19 mL. Length: 104 in. (264.2 cm.). | 354207     | 50            |
|                                                                                                                                                                                                                            |            |               |
| (15 drops/mL.) Universal spike, roller clamp, 1.2 micron filter, slide clamp, CARESITE Needleless Connector injection site 6 in. above distal end, SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 23 mL. Length: 98 in. (248.9 cm.).                                       | 470115     | 50            |
|                                                                                                                                                                                                                          |            |               |
| <b>Gravity IV Administration Sets with CARESITE Needleless Injection Sites</b>                                                                                                                                                                                                                               |            |               |
| (15 drops/mL.) Universal spike, slide clamp, roller clamp, CARESITE Needleless Connector injection site 6 in. above distal end and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 20 mL. Length: 111 in. (282 cm.).                                                        | 354201     | 50            |
|                                                                                                                                                                                                                          |            |               |
| (15 drops/mL.) Universal spike, backcheck valve, CARESITE Needleless Connector injection sites 68 in. and 6 in. above distal end, roller clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 17 mL. Length: 84 in. (210 cm.).                                         | 354202     | 50            |
|                                                                                                                                                                                                                          |            |               |
| (15 drops/mL.) Universal spike, slide clamp, backcheck valve, CARESITE Needleless Connector injection sites 87 in. and 12 in. above distal end, roller clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 20 mL. Length: 104 in. (264 cm.).                          | 354203     | 50            |
|                                                                                                                                                                                                                          |            |               |


| Description                                                                                                                                                                                                                                                                                                                                                                           | REF<br>Number | Case<br>Quantity |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|------------------|
| Gravity IV Administration Sets with CARESITE Needleless Connector Injection Sites (continued)                                                                                                                                                                                                                                                                                         |               |                  |
| <p>(60 drops/mL.) Universal spike, slide clamp, backcheck valve, CARESITE Needleless Connector injection sites 87 in. and 12 in. above distal end, roller clamp and SPIN-LOCK® Connector. Not made with DEHP or natural rubber latex. Priming volume: 19 mL. Length: 103 in. (262 cm.).</p>         | 354204        | 50               |
| <p>(15 drops/mL.) Universal spike, backcheck valve, CARESITE Needleless Connector injection sites 85 in., 41 in. and 6 in. above distal end, roller clamp, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 20 mL. Length: 112 in. (285 cm.).</p>   | 354205        | 50               |
| <p>(60 drops/mL.) Universal spike, backcheck valve, CARESITE Needleless Connector injection sites 61 in., 26 in. and 6 in. above distal end, roller clamp, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 17 mL. Length: 88 in. (223.5 cm.).</p>  | 354206        | 50               |

# CARESITE Micro<sup>®</sup> Needleless Connectors

| Description                                                                                                                                                                                                                                                                                                    | REF Number | Case Quantity |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| CARESITE Micro Needleless Connectors                                                                                                                                                                                                                                                                           |            |               |
| <p>CARESITE Micro Needleless Connector for aspiration, injection or gravity flow of fluid upon insertion of a male Luer fitting. Not made with DEHP PVC, BPA, or natural rubber latex. Priming Volume: 0.04 mL. Power injection rated to 400 psi and 10 mL/sec. <b>High Pressure Rated</b></p>                 |            |               |
|                                                                                                                                                                                                                              | 415400     | 200           |
| CARESITE Micro Small Bore Extension Sets                                                                                                                                                                                                                                                                       |            |               |
| <p>Small bore extension set with non-removable CARESITE Micro Needleless Connector, slide clamp and SPIN-LOCK<sup>®</sup> Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.3 mL. Length: 8 in (20 cm.). Power injection rated to 325 psi and 10 mL/sec. <b>High Pressure Rated</b></p> |            |               |
|                                                                                                                                                                                                                              | 354265     | 100           |
| <p>Small bore extension set with removable CARESITE Micro Needleless Connector, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.3 mL. Length: 8 in (20 cm.). Power injection rated to 325 psi and 10 mL/sec. <b>High Pressure Rated</b></p>                 |            |               |
|                                                                                                                                                                                                                            | 354269     | 100           |
| <p>Small bore extension set with non-removable CARESITE Micro Needleless Connector, slide clamp and SPIN-LOCK<sup>®</sup> Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.8 mL. Length: 17 in (43 cm.).</p>                                                                           |            |               |
|                                                                                                                                                                                                                            | 354270     | 100           |
| <p>Small bore extension set with removable CARESITE Micro Needleless Connector, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.8 mL. Length: 17 in (43 cm.).</p>                                                                                           |            |               |
|                                                                                                                                                                                                                            | 354271     | 100           |

| Description                                 | REF Number | Case Quantity |
|---------------------------------------------|------------|---------------|
| CARESITE Micro Standard Bore Extension Sets |            |               |


Standard bore extension set with non-removable CARESITE Micro Needleless Connector, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.4 mL. Length: 8 in (20 cm.). Power injection rated to 325 psi and 10 mL/sec. **High Pressure Rated**



354267

100


Standard bore extension set with removable CARESITE Micro Needleless Connector, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.4 mL. Length: 8 in (20 cm.). Power injection rated to 325 psi and 10 mL/sec. **High Pressure Rated**



354268

100


Standard bore extension set with non-removable CARESITE Micro Needleless Connector, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 1.0 mL. Length: 8 in (20 cm.).



354274







100

Standard bore extension set with removable CARESITE Micro Needleless Connector, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 1.0 mL. Length: 8 in (20 cm.).












354266

100

| Description                                                                                                                                                                                                                                                                                                              | REF Number | Case Quantity |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| <b>CARESITE Micro Bifurcated Extension Sets</b>                                                                                                                                                                                                                                                                          |            |               |
| Smallbore high pressure rated Y-extension set with two non-removable CARESITE Micro Needleless Connectors, two slide clamps and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.5 mL. Length 8 in. (20cm.) Power injection rated to 300 psi and 10 mL/sec. <b>High Pressure Rated</b> | 354272     | 50            |
|                                                                                                                                                                                                                                        |            |               |
| Smallbore Y-extension set with one removable CARESITE Micro Needleless Connector, one female Luer access port, two slide clamps and SPIN-LOCK Connector, 8 in., Priming volume: 0.5 mL. <b>High Pressure Rated</b>                                                                                                       | 354273     | 50            |
|                                                                                                                                                                                                                                        |            |               |
| <b>CARESITE Micro STEADYCare™ Extension Sets</b>                                                                                                                                                                                                                                                                         |            |               |
| Small bore extension set with non-removable CARESITE Micro Needleless Connector, slide clamp, SPIN-LOCK Connector and STEADYCare Wedge. Not made with DEHP or natural rubber latex. Priming volume: 0.3 mL. Length: 8 in. (20cm.). Power injection rated to 325 psi and 10 mL/sec. <b>High Pressure Rated</b>            | 354276     | 100           |
|                                                                                                                                                                                                                                      |            |               |
| Standard bore extension set with non-removable CARESITE Micro Needleless Connector, slide clamp, SPIN-LOCK Connector, and STEADYCare Wedge. Not made with DEHP or natural rubber latex. Priming volume: 0.4 mL. Length: 8 in. (20cm.). Power injection rated to 325 psi and 10 mL/sec. <b>High Pressure Rated</b>        | 354278     | 100           |
|                                                                                                                                                                                                                                      |            |               |
| Standard bore extension set with non-removable CARESITE Micro Needleless Connector, slide clamp, SPIN-LOCK Connector and STEADYCare Wedge. Not made with DEHP or natural rubber latex. Priming volume: 1.0 mL. Length: 8 in. (20cm.).                                                                                    | 354275     | 100           |
|                                                                                                                                                                                                                                      |            |               |
| Standard bore extension set with removable CARESITE Micro Needleless Connector, slide clamp, SPIN-LOCK Connector and STEADYCare Wedge. Not made with DEHP or natural rubber latex. Priming volume: 1.0 mL. Length: 8 in. (20cm.).                                                                                        | 354277     | 100           |
|                                                                                                                                                                                                                                      |            |               |

# SAFELINE® Needleless Connector

| Description                                                                                                                                                                                                                                                                     | REF Number | Case Quantity |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| Extension Sets with SAFELINE Split Septum Injection Sites                                                                                                                                                                                                                       |            |               |
| SAFELINE Needleless Connector split septum injection site, slide clamp, microbore tubing and male Luer lock connector. Not made with DEHP or natural rubber latex. Priming volume: 0.42 mL. Length: 7 in. (17.8 cm.).                                                           |            |               |
|                                                                                                                                                                                               | NF1310     | 50            |
| SAFELINE Needleless Connector split septum injection site, slide clamp, microbore tubing and male Luer slip connector. Not made with DEHP or natural rubber latex. Priming volume: 0.5mL. Length: 8 in. (20.32 cm.).                                                            |            |               |
|                                                                                                                                                                                               | NF1312     | 50            |
| Y-type small bore extension set with attached SAFELINE Needleless Connector split septum injection site, female adapter, two slide clamps and male Luer lock connector. Not made with DEHP or natural rubber latex. Priming volume: 0.5 mL. Length: 7 in. (17.8 cm.).           |            |               |
|                                                                                                                                                                                              | NF1318     | 100           |
| Small bore extension set with SAFELINE Needleless Connector split septum injection site, slide clamp and SPIN-LOCK® Connector. Not made with DEHP or natural rubber latex. Priming volume: 0.63 mL. Length: 7 in. (17.8 cm.).                                                   |            |               |
|                                                                                                                                                                                             | NF1320     | 50            |
| Female adapter, microbore tubing, removable slide clamp, T-connector with SAFELINE Needleless Connector split septum injection site and SPIN-LOCK Connector. Not made with natural rubber latex. Priming volume: 0.27 mL. Length: 6 in. (15.2 cm.).                             |            |               |
|                                                                                                                                                                                             | NF1330     | 50            |
| Female adapter, large bore tubing, slide clamp, SAFELINE Needleless Connector split septum injection sites 28 in. and 6 in. above distal end, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 6.0 mL. Length: 35 in. (89 cm.). |            |               |
|                                                                                                                                                                                             | NF1370     | 50            |


| Description                                                                                                                                                                                                                                                                            | REF Number | Case Quantity |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------|
| <b>Filterflow® Filtered Add-on Extension Set with SAFELINE Split Septum Injection Site</b>                                                                                                                                                                                             |            |               |
| <p>Female Luer lock, 0.2 micron filter, SAFELINE Needleless Connector split septum injection site 6 in. above distal end, pinch clamp and SPIN-LOCK® Connector. Not made with DEHP or natural rubber latex. Priming volume: 5.3 mL. Length 16 in. (41 cm.).</p>                        |            |               |
|                                                                                                                                                                                                      | 473993     | 50            |
| <p>Female Luer lock, 1.2 micron filter, SAFELINE Needleless Connector split septum injection site 6 in. above distal end, pinch clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 5.1 mL. Length 16 in. (40.0 cm.).</p>                       |            |               |
|                                                                                                                                                                                                      | 473998     | 50            |
| <b>Basic Administration Sets with SAFELINE Split Septum Injection Site</b>                                                                                                                                                                                                             |            |               |
| <p>(15 drops/mL.) Universal spike, roller clamp, SAFELINE Needleless Connector split septum injection sites 40 in. and 6 in. above distal end, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 19.0 mL. Length: 90 in. (228 cm.).</p> |            |               |
|                                                                                                                                                                                                      | NF1150     | 50            |
| <p>(60 drops/mL.) Universal spike, roller clamp, SAFELINE Needleless Connector split septum injection sites 38 in. and 6 in. above distal end, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 18 mL. Length: 94 in. (239 cm.).</p>   |            |               |
|                                                                                                                                                                                                    | NF1190     | 50            |
| <b>Secondary Set for use with SAFELINE Split Septum Injection Site</b>                                                                                                                                                                                                                 |            |               |
| <p>(15 drops/mL.) Universal spike, roller clamp and SPIN-LOCK Connector, SAFELINE Clip Lock® Cannula, 10.5" plastic bag hanger. Not made with DEHP or natural rubber latex. Priming volume: 10.2 mL. Length: 49 in. (125 cm.).</p>                                                     |            |               |
|                                                                                                                                                                                                    | NF1430     | 50            |

| Description                                                                                                | REF Number | Case Quantity |
|------------------------------------------------------------------------------------------------------------|------------|---------------|
| <b>ADDitIV® Primary Administration Sets with Backcheck Valve and SAFELINE Split Septum Injection Sites</b> |            |               |

(15 drops/mL.) Universal spike, backcheck valve, SAFELINE Needleless Connector split septum injection sites 80 in. and 6 in. above distal end, roller clamp and SPIN-LOCK® Connector. Not made with DEHP or natural rubber latex. Priming volume: 19 mL. Length: 104 in. (264 cm.).

352604


50



(15 drops/mL.) Universal spike, backcheck valve, SAFELINE Needleless Connector split septum injection sites 88 in., 30 in. and 6 in. above distal end, roller clamp, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 19.68 mL. Length: 111 in. (282 cm.).

NF1251


50



(60 drops/mL.) Universal spike, backcheck valve, SAFELINE Needleless Connector split septum injection sites 62 in., 27 in. and 6 in. above distal end, roller clamp, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 17 mL. Length: 88 in. (223.5 cm.).

NF1290

50




### Filterflow® Filtered Gravity Administration Sets with SAFELINE Split Septum Injection Sites

(15 drops/mL.) Universal spike, 0.2 micron filter, roller clamp, slide clamp, SAFELINE Needleless Connector split septum injection site 6 in. above distal end and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 21.0 mL. Length: 89 in. (226 cm.).

352434


50



(60 drops/mL.) Universal spike, backcheck valve, 0.2 micron filter, roller clamp, SAFELINE Needleless Connector split septum injection sites 38 in. and 6 in. above distal end, slide clamp and SPIN-LOCK Connector. Not made with DEHP or natural rubber latex. Priming volume: 18 mL. Length: 94 in. (239 cm.).


NF1210F

50



| Description                                     | REF<br>Number | Case<br>Quantity |
|-------------------------------------------------|---------------|------------------|
| SAFELINE Needleless Split Septum Injection Site |               |                  |

SAFELINE Needleless Connector split septum injection site, male Luer lock connector. Not made with DEHP or natural rubber latex. Priming volume: 0.25 mL.




NF9100

400

### Accessories for SAFELINE Split Septum Injection Site


SAFELINE Clip Lock® Cannula; locking cannula with female Luer. Not made with DEHP or natural rubber latex.



NF9200

400


SAFELINE Blunt Cannula; blunt cannula with female Luer. Not made with DEHP or natural rubber latex.



NF9210

400

Medic® Plastic Anti-Stick Needle/Connector. Not made with DEHP or natural rubber latex.



M8-5007

1,000

# General Ordering Information and Terms of Sale

Thank you for your interest in B. Braun Medical Inc. ("B. Braun") products. The following provisions set forth the terms and conditions on which B. Braun sells its products.

**Terms of Sale:** Purchase of any products sold by B. Braun shall be subject to and expressly limited by the terms and conditions contained herein. No changes to, waiver of, or addition to any of these terms and conditions shall be effective unless agreed to in writing and signed by B. Braun. Buyer acknowledges and agrees that these terms and conditions supersede the terms and conditions of any purchase order or other documentation used by Buyer and, except for delivery and billing addresses and quantities prices and items ordered, any conflicting or additional terms are void and have no effect, but that Buyer may place orders by use of purchase orders and other documentation for its convenience purposes only. Notwithstanding the foregoing, B. Braun reserves the right at any time to amend these terms and conditions and Buyer shall be deemed to accept such amended terms and conditions by ordering products herein offered after the date of such amendment. Additional special terms and conditions of B. Braun may be applicable with respect to certain products.

**Minimum Order Quantity:** All products must be ordered in full case quantities. Orders received with a value of less than \$250.00 will be subject to a handling charge of \$40.00. Wholesalers and Distributors are subject to a handling charge of \$75.00 for orders less than \$1,000.00.

**Prices:** All pricing quotes must be documented in writing and signed by B. Braun to be valid. Prices quoted, unless otherwise stated, refer to full case quantities. All prices are firm for 30 days from the date quoted. B. Braun reserves the right to change the prices and specifications of its products at any time without notice.

Products purchased from distributors will be at the prices set by those distributors and subject to the terms and conditions of sale of the distributor.

**Tax Information:** Any tax, assessment, duty, custom or other fee of any nature imposed upon the products, their sale, transportation, delivery, use or consumption shall be paid by Buyer in addition to the price quoted or invoiced. If B. Braun is required to prepay any such tax or fee, Buyer will reimburse B. Braun. Buyer must provide B. Braun with a resale/exemption certificate in order to avoid the withholding of applicable taxes. No refund or adjustment to previously withheld taxes will be made by B. Braun sixty (60) days after the invoice date. Proof of certification should be mailed to: B. Braun Medical Inc., 824 Twelfth Avenue, Bethlehem, PA 18018-0027.

**Payment:** Payment terms are net 30 days from date of invoice. Credit card payments will not be accepted without prior approval from B. Braun. Buyer's obligation to pay outstanding invoices and all other amounts is absolute and unconditional and is not subject to any abatement, reduction, set-off, defense, counterclaim, interruption, deferment or recoupment for any reason whatsoever. Balances remaining unpaid at due date are subject to a interest charge of 1.5% per month or the highest rate permitted by law, whichever is lower, until paid.

Any discounts, rebates, administrative fees, credits, or other fees due or owed to Buyer will be applied against delinquent balances before payment or reimbursement is made.

Any disputed amounts should be reported immediately and remitted with the undisputed amount by the payment due date. If B. Braun agrees with the billing dispute, B. Braun will credit Buyer the amount of the agreed-upon billing dispute. All billing disputes must be made within six (6) months of the applicable invoice date, or will be deemed to be waived.

B. Braun reserves the right in its sole discretion to require prepayment from any Buyer at any time and may refuse to sell and/or withhold further shipment until all overdue balances are made current. Buyer shall be liable for and shall reimburse B. Braun for all costs and expenses it may incur in connection with collection of any amounts owed to B. Braun or enforcement of its rights, including without limitation, reasonable attorneys' fees and expenses, court costs and cost of collection agencies.

**Returned Goods Policy:** Product returns are subject to the following conditions:

1. All returns are subject to the prior authorization of B. Braun, in its discretion. Buyer must notify Customer Support at (800) 227-2862 and complete a Return Request Form. The Return Request Form requires lot numbers, quantities and catalog numbers along with a specific reason for return. Customer Support will either authorize or deny the request for return. Only items appearing on an approved Return document are acceptable for return. Product returns will only be accepted from the original purchaser. Product returns will not be accepted from third party return companies. Unauthorized returns will be destroyed and no credit issued. All authorized returned goods must be shipped freight prepaid to the B. Braun location indicated on the Return document, except B. Braun will pay freight costs for Product shipped-in-error.
2. All Products returned within 30 days of delivery are subject to a 25% restocking charge, except for Products shipped-in-error. Products returned after 30 days are considered excess stock and will be eligible for 50% credit. No Products will be accepted for return after 90 days from the date of delivery.

3. Subject to paragraph 2 above, credit will be issued at the net purchase price for all authorized returns provided:
  - a. Products have at least one year shelf-life remaining, or products with original expiration dating of eighteen months or less have at least six months shelf-life remaining.
  - b. Proper authorization has been obtained prior to return of products.
  - c. Products are in their original packaging.
  - d. Products are current inventory items.
  - e. The Products have been shipped and billed to Buyer by B. Braun and Buyer has paid for said Products.
4. Certain Products are not eligible for return. These are:
  - a. Products that have deteriorated because of improper storage, handling, abuse or other factors.
  - b. Products that have been opened, partly used or which the labels or seals have been removed or tampered.
  - c. Products that have been involved in a special promotion sale.
  - d. Broken, damaged or opened cases. Resealed cartons are not eligible for return.
  - e. Special products made to Buyer's specification.
  - f. Sets over two (2) years old, flush syringe products, infusion systems devices, related accessories and IV poles.
5. Products that are to be returned for repair (e.g., medical device equipment) must be accompanied by an approved Repair Notification.

**Expiration Date:** On dated products, the expiration date is shown as a month and year, e.g., August 2012, 8/12. The date of expiration is the last day of the given month, e.g., August 31, 2012.

Credits for returned goods are conditioned upon B. Braun's inspection and approval of such goods upon their return. If B. Braun determines, in its discretion, that any returned goods are not eligible for return due to any of the reasons provided in paragraph 4 above, Buyer will not receive a credit, even if an Return document was issued. No advance credits will be accepted.

**Warranty:** With respect to disposable products, B. Braun warrants to the original purchaser that, at time of delivery, each standard product manufactured by B. Braun shall be free of defects in material and workmanship and, when used for the purposes and indications described on the labeling, is fit for the purposes and indications described on the labeling. All warranties for a product shall expire as of product expiration date, or if none, after one (1) year from the date of shipment from B. Braun. B. Braun's warranty hereunder shall not apply if: (i) a product is not used in accordance with its instructions or if it is used for a purpose not indicated on the labeling, (ii) any repairs, alterations or other work has been performed by Buyer or others on such item, other than work performed with B. Braun's authorization and according to its approved procedures; or (iii) the alleged defect is a result of abuse, misuse, improper maintenance, accident or the negligence of any party other than B. Braun. The warranty set forth herein is conditioned upon proper storage, installation, use and maintenance in accordance with applicable written recommendations of B. Braun. The warranty furnished hereunder does not extend to damage to items purchased hereunder resulting in whole or in part from the use of components, accessories, parts or supplies not furnished by B. Braun.

B. Braun's sole obligation shall be to repair or replace, at B. Braun's option, any defective component or item and pay transportation expenses for such replacement. Buyer shall provide labor for the removal of the defective component or item and installation of its replacement at no charge to B. Braun. Buyer shall bear all risk of loss or damage to returned goods while in transit. In the event no defect or breach of warranty is discovered by B. Braun upon receipt of any returned item, the item will be returned to Buyer at Buyer's expense and Buyer will reimburse B. Braun for the transportation charges, labor and associated charges incurred in testing the allegedly defective item.

Except as expressly provided herein, B. Braun makes no representation or warranty of any kind, expressed or implied with respect to any products, parts or services provided by B. Braun including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose. Products distributed, but not manufactured, by B. Braun are not warranted by B. Braun and Buyer must instead rely on the representations and warranties, if any, provided directly to Buyer by the manufacturer of such product. The sole and exclusive remedy for breach of any warranty is limited to the remedies provided in the paragraph above.

Warranties for infusion system devices are provided in a separate warranty document.

#### Miscellaneous:

1. Any required notices shall be given in writing, in the case of B. Braun, at the address set forth below and in the case of Buyer, at the address designated on Buyer's purchase order or to such other address as either party may substitute by written notice to the other and shall be deemed given upon personal delivery, overnight delivery or three days following deposit in the mail.
2. Except as expressly provided herein, no changes or modifications to, or waiver of, any of these terms and conditions shall be valid or binding on either party unless in writing and signed by an authorized representative of each party. B. Braun's failure or delay to exercise or enforce any of its rights hereunder shall not constitute or be deemed to be a waiver of such rights or forfeiture of such rights and B. Braun may, at its option, from time to time, exercise any of its rights or remedies.
3. These Terms bind Buyer and its successors and permitted assigns.
4. B. Braun will use its reasonable efforts to fill orders, but B. Braun shall not be liable for nonperformance or delays caused by a shortage of raw materials, manufacturing problems, delivery or labor problems, priorities, acts of regulatory agencies or judicial bodies, discontinuation of a product line, acts of God or third parties, infringement claims, or other causes beyond its reasonable control. Buyer agrees that in such events B. Braun may allocate products among all purchasers as it deems reasonable, without liability.  
B. Braun reserves the right from time to time to substitute a product with a product that has the same function as such product, or to delete a product.
5. The products are sold subject to Pennsylvania law. These terms and any dispute or claim relating to these terms or the sale of products ("Claim") shall be governed by and construed under Pennsylvania law, notwithstanding its law of conflicts of law. If any Claim cannot be settled amicably between the parties, such Claim shall be tried by a court and not a jury. Buyer expressly and unconditionally waives its rights to a jury trial in any such Claim.
6. B. Braun shall promptly deliver written notice or verbal, followed by written, notice of any recall of Product. Should the recall require Products be removed from the market, B. Braun shall, to the extent reasonably possible and at B. Braun's option, replace any such recalled Products as soon as practicable with comparable Products not subject to such recall or repair any such recalled Products and return them to Buyer.
7. Buyer may not change, adulterate, obscure, remove or deface trademarks, tradenames or labels appearing on any Product of B. Braun.
8. If the pricing offered by B. Braun to a Buyer constitutes a discount or other reduction in price under Section 1128(b)(3)(a) of the Social Security Act 42 U.S.C. 1320a-7b(b) (3)(a) and C.F.R. § 1001.952(h), Buyer shall disclose the discount or reduction in price to the fullest extent required under any state or federal program that provides cost or charge-based reimbursement to Buyers for products. This act requires, among other things, that Buyer fully and accurately report on any claim or request for payment it submits to Medicare and Medicaid the actual purchase price paid by Buyer for products, net of any discounts, rebates or allowances provided hereunder. Buyer may also be required, upon request, to provide documentation of the discount or other reduction in price to the Secretary of Health and Human Services.
9. In all cases, all intellectual property rights in and to and all technology relating to, the Products supplied to Buyer, their design and all improvements thereto or thereof, whether or not such Product, design or improvement is made pursuant to Buyer's specifications or at Buyer's expense, shall be and remain the exclusive property of B. Braun.
10. Any knowledge or information that Buyer may disclose to B. Braun shall not be deemed to be confidential or proprietary information and shall be acquired by B. Braun free from any restriction.
11. B. Braun shall not in any event be liable to buyer for any indirect, incidental, special, punitive or consequential damages (including any damage for lost profits), or otherwise arising out of or in connection with furnishing of products, parts or service hereunder, or the performance, use of, or inability to use any products, parts or service, or otherwise, whether based in contract, warranty, tort, including without limitation, negligence, or any other legal or equitable theory. B. Braun's total liability for any claim or action shall not exceed the purchase price of the products out of which such claim or action arose.
12. If any of the products are medical devices, Buyer acknowledges that it is familiar with the Safe Medical Devices Act of 1990 (the "Devices Act") and the reporting obligations imposed on device users thereunder. In this regard, Buyer agrees to notify B. Braun within ten (10) days of the occurrence of any event identified in the Devices Act imposing a reporting obligation on Buyer and/or B. Braun (except for events representing an imminent hazard that require notification to the United States Food and Drug Administration (the "FDA") within seventy-two hours, in which case, such notice will be delivered to the other party immediately).

Buyer shall maintain adequate tracking for the products to enable B. Braun to meet the FDA requirements applicable to the tracking of medical devices.

13. Buyer shall purchase the products for its own use only and shall not resell the products to any other party. Buyer represents it has examined the products and that they are acceptable and clinically suitable for its intended purposes.

#### Shipping:

1. Freight charges will be prepaid when orders are shipped via a transportation mode and carrier selected by B. Braun. Unless otherwise stated in an agreement, shipment of infusion systems devices will be made F.O.B. Origin.
2. When expedited delivery, specialized service, or alternate transportation modes are requested, or if requests are inconsistent with efficient distribution practices, an additional charge to cover the premium expense will be added to the invoice. Inside deliveries are subject to additional charges based on current common carrier rates or hourly rates. Palletized deliveries are encouraged.
3. Products will be shipped on a scheduled order and delivery basis. Orders shipped outside Buyer's order and delivery schedule may be subject to a handling charge of \$75.00. B. Braun will provide proof of delivery upon request.

**Damage or Loss in Transit:** Identity of items and extent of damage or loss must be noted on Buyer's copy of delivery document by the agent of the transportation company. If damage is discovered after receipt of shipment, notify the transportation company immediately and request that inspection be made and an inspection report rendered. Buyer must report concealed shortages or damages within palletized shipments to B. Braun Customer Support within 3 business days of delivery or credit will not be allowed. In addition, Buyer must provide B. Braun with a copy of Buyer's claim request accompanied by a delivery receipt or an inspection report upon which the transportation company has properly noted damage or loss and B. Braun will issue a credit for the loss or damage and file a claim with the carrier. If such information is not received within ten (10) days of delivery, no credit will be issued.

**Send copy of carrier freight bill to Customer Support indicating item and quantity damaged or not received.**

**Count and inspect your freight before carrier departs. Damaged merchandise should not be accepted.**

**Please forward all information to Attn: Customer Support, B. Braun Medical Inc., 901 Marcon Blvd., Allentown, PA 18109. (800) 227-2862.**

---

**For Order Placement:** B. Braun welcomes orders either by telephone, fax, mail or electronically. For telephone orders, please call (800) 227-2862. Faxed orders should be sent to: (610) 266-6122. For infusion devices, call (800) 627-7867 or fax (610) 266-2429. In Canada, please call (800) 624-2920 or fax (800) 624-2939. Mail orders should be directed to Attn: Customer Support, B. Braun Medical Inc., 901 Marcon Blvd, Allentown, PA 18109.

---

# Trademarks

Trademarks of B. Braun Medical Inc. and its affiliates:

|                       |            |              |              |                 |
|-----------------------|------------|--------------|--------------|-----------------|
| ADD-PRO™              | CONTIPLEX™ | Hyperformer™ | RED CAP™     | UniSpike™       |
| Adimea™               | Cytoguard™ | Injekt™      | SPACE™       | ULTRABLOCK™     |
| Apex™                 | DoseCom™   | Physiolyte™  | SAFEPORT™    | ULTRAPORT™ Zero |
| AutoCOMPLETE™         | DoseGuard™ | PIC™         | Ster-ASSIST™ |                 |
| B-Smart™              | DoseScan™  | PinPad™      | Superset™    |                 |
| CARESAFE™             | EPICAN™    | QC Tester™   | Titan XL™    |                 |
| Chemo Dispensing Pin™ | Flexima™   | RateGuard™   | Transofix™   |                 |

Registered trademarks of B. Braun Medical Inc. and its affiliates:

|                       |             |               |              |              |             |
|-----------------------|-------------|---------------|--------------|--------------|-------------|
| addEASE®              | DiaLines®   | FILTER STRAW® | MINI-SPIKE®  | PERIFIX®     | SPIN-LOCK®  |
| ADDitIV®              | Dialog+®    | FreAmine®     | MULTI-AD®    | Physiolyte®  | SPINOCAN®   |
| BIG TAB®              | Diapact®    | FreAmine HBC® | NephrAmine®  | PINNACLE®    | STIMUPLEX®  |
| CAPS®                 | DISCOFIX®   | H.E.L.P.®     | Nutriliplid® | Plasmat®     | TrophAmine® |
| CARESITE®             | DoseTrac®   | HepatAmine®   | Omnican®     | PLEXUFIX®    | ULTRAPORT®  |
| CARESITE Micro®       | Duosol®     | Horizon NXT®  | OMNIFIX®     | ProcalAmine® | ULTRASITE®  |
| ClearChoice-DCB®      | DUPLEX®     | HESpan®       | ONGUARD®     | Prontosan®   | VISTA®      |
| Clip Lock®            | E3®         | Hyperlyte®    | Outlook®     | Rate Flow®   | WHIN®       |
| Diacap®               | ESPOCAN®    | Infusomat®    | PAB®         | SAFELINE®    |             |
| Diacap a Polysulfone® | EXCEL®      | Isolyte®      | PENCAN®      | SAFSITE®     |             |
|                       | Filterflow® | MICRO PIN®    | Perfusor®    | Solcart B®   |             |

BodyGuard, CMExpress and BodyComm and ColorVision are registered trademarks of CME America, LLC.

ChlorPrep and SEPP are registered trademarks of CareFusion 2200, Inc.

DOCit and HubScrub are trademarks of Saxa Medical Solutions.

T-Pak is a registered trademark of Allegiance.

CIV-Flex is a trademark of CIVCO Medical Instruments Company, Inc.

DuraPrep, Steri-Strip and Tegaderm are registered trademarks of 3M.

Medic is a registered trademark of Medisystems Corporation.

Microsoft Access and Windows are registered trademarks of Microsoft Corporation.

Monoject is a registered trademark of Covidien, AG.

NeedleVISE is a registered trademark of Atrion Medical Products, Inc.

OmniPal is a trademark and Patient Pal, Trans-Pal, Versa Pole and Walk 'N' Roller are registered trademarks of Pryor Products, Inc.

OpSite is a registered trademark of T.J. Smith and Nephew, Ltd.

PeraChek, PERACIDIN and PeraSidual are registered trademarks of Angelini Pharma Inc.

Plastipak, SafetyGlide and Vacutainer are registered trademarks of Becton Dickinson and Co.

SAF-T-POLE is a trademark of American Medical Manufacturing, Inc.

SPROTTE is a registered trademark of Pajunk GmbH.

Statlock is a registered trademark of C.R. Bard, Inc.

Steripick, Locksite, Buttonhole, MasterGuard, FingerShield, SecureClip and StreamLine are registered trademarks of NxStage Medical, Inc.

SUPOR is a registered trademark of the Pall Corporation.

Tevadaptor is a registered trademark of Teva Medical Ltd.

Trissel's is a trademark of Lawrence A. Trissel.

Velcro is a registered trademark of Velcro Industries B.V.

**For more information on IV Sets and Accessories,  
call **800-227-2862** or visit **[bbraunusa.com](https://www.bbraunusa.com)**.**