





CARESAFE™ IV Administration Set Available with AirStop

Designed for Convenient and Safe IV Therapy

CARESAFE™ IV Administration Sets



Clinicians are frequently under high pressure and looking to streamline workflow

Designed to Streamline Workflow


- Nurses spend only 35% of their time in the patient's room.¹
- According to a blinded B. Braun market research study, 77% of respondents specified experiencing a high level of stress associated with the IV set running dry.⁶
- Re-priming of IV lines diverts clinician attention from more value-added patient care.

Promotes Patient Safety

- Air-in-line prevention is critical in IV therapy.^{2,3}
- The innovative AirStop filter is designed to only allow fluid to pass and thereby helps to prevent the line from running dry.⁴
- CARESAFE IV administration sets are not made with DEHP or PVC



AirStop filter won't let the line run dry.⁴





Long drop former

Transparent for easy viewing


Unique ring designed to assist easy spiking

AirStop filter helps prevent air infusion and prevents the line from running dry.⁴

Drives Environmental Sustainability

The only supplier to offer a robust portfolio of IV solutions and sets not made with DEHP and PVC .

- EPA is concerned about phthalates because of their toxicity and the evidence of pervasive human and environment exposure to them.⁷
- DEHP is one of the 20 chemicals on the EPA's High-Priority Chemical Substances proposal.⁸
- When deposited into landfills, products made with PVC can leach toxic chemicals into the groundwater.⁹
- CARESAFE IV administration set is not made with DEHP, which has been shown to be harmful to certain patient populations including pregnant women, male neonates and young children.⁵




CARESAFE™ IV Administration Sets

Ordering Information



WITH AIRSTOP



WITHOUT AIRSTOP



Ref Number

With Airstop

Administration Set, 20 drop Drip Chamber with AirStop, with 2 CARESITE® Y-Site NC	354302
Administration Set, 20 drop Drip Chamber with AirStop, with 3 CARESITE Y-Site NC	354304
Administration Set, 20 drop Drip Chamber with AirStop, with 2 CARESITE Y-Site NC and a 0.2 µm filter	354306
Administration Set, 20 drop Drip Chamber with AirStop, with 2 CARESITE Y-Site NC and two removable 4-way ULTRAPORT® Stopcocks	354308

Without Airstop

Administration Set, 20 drop Drip Chamber with 15 µm filter, with 2 CARESITE Y-Site NC	354303
Administration Set, 20 drop Drip Chamber with 15 µm filter, with 3 CARESITE Y-Site NC	354305
Administration Set, 20 drop Drip Chamber with 15 µm filter, with 2 CARESITE Y-Site NC and a 0.2 µm filter	354307
Administration Set, 20 drop Drip Chamber with 15 µm filter, with 2 CARESITE Y-Site NC and two removable 4-way ULTRAPORT Stopcocks	354310
Secondary Set, 20 drop Drip Chamber with 15 µm filter	V1931
Extension Set, standard bore with 2 CARESITE Y-Sites	354320
Extension Set, standard bore with 1 bonded CARESITE NC	471858
Extension Set, standard bore with 1 removable CARESITE NC	471859

- Yen, P. Y., Kellye, M., Lopetegui, M., Saha, A., Loversidge, J., Chipps, E. M., Gallagher-Ford, L., & Buck, J. (2018). Nurses' Time Allocation and Multitasking of Nursing Activities: A Time Motion Study. AMIA ... Annual Symposium proceedings. AMIA Symposium, 2018, 1137–1146.
- Cook, Lynda S. Infusion-Related Air Embolism. Journal of Infusion Nursing. February 2013. Volume 36 Number 1: 25-34.
- Infusion Therapy Standards of Practice 2021. January/February 2021. Volume 44 Number 1S.
- B. Braun Data on File
- Schettler, T. Polyvinyl chloride in health care. Health Care Without Harm. January 2020.
- Key Group Consultative Market Research: New PVC/DEHP-Free IV Sets: New Product Concept & Conjoint Study. January 14 2021.
- U.S Environmental Protection Agency: Phthalates Action Plan. March 2012.
- EPA [2019] EPA Finalizes List of Next 20 Chemicals to Undergo Risk Evaluation under TSCA. [Press Release]. 20/12.
- Belliveau, M., Lester, S. PVC Bad News Comes in 3's: The Poison Plastic, Health Hazards and the Looming Waste Crisis. Center for Health, Environment and Justice Environment Health Strategy Center. December 2004.

B. Braun Medical Inc. | Bethlehem, PA | 800-227-2862 | bbraunusa.com